Shanika Craig
Designated Federal Officer (DFO)
Medical Devices Advisory Committee
Center for Devices and Radiological Health
10903 New Hampshire Ave
White Oak Bldg 66 Rm. 1613
Silver Spring, MD 20993

RE: PMA P020014, BAY 1454032
Bayer's Final Executive Summary for FDA Panel Meeting
Essure System for Permanent Birth Control

Dear Ms. Craig:

Reference is made to PMA P020014, Essure System for Permanent Birth Control and to the email dated July 15, 2015 from the Designated Federal Officer (DFO) regarding the upcoming meeting of the Obstetrics and Gynecology Devices Panel of the Medical Devices Advisory Committee scheduled for September 24, 2015.

The purpose of this submission is to provide FDA with Bayer's Final Executive Summary which provides data to support the positive benefit-risk profile of Essure.

This Executive Summary is being submitted via UPS as 5 paper binders and 30 CD packs in electronic form. The eCopy is an exact duplicate of the paper copy except for the literature references, which are included as PDF files on the eCopy.

Please let me know if you have any questions or need additional information. I can be reached by phone at (908) 200-9718 or by email at manal.morcos@bayer.com.

Bayer HealthCare Pharmaceuticals certifies that this submission has been scanned for viruses and is virus free using TREND MICRO™ Office Scan™, Program Version 10.6 or higher. For any questions regarding eCopy technical aspects of this electronic submission, please contact Maria Ferrara at (862) 404-3421 or by email at maria.ferrara@bayer.com.

September 3, 2015

Bayer HealthCare, LLC
Woman's Health.
Global Regulatory Affairs
100 Bayer Blvd
Whippany, NJ 07981

Phone: (862) 404-5729
Cell: (908) 200-9718
manal.morcos@bayer.com

Bayer HealthCare

Page 2 of 2
Respectfully,
Bayer HealthCare LLC.

Kevin Hanley
On behalf of:

A handwritten signature in black ink, appearing to read "Manal Morcos". The signature is fluid and cursive, with a long horizontal stroke extending to the right.

Manal Morcos, MS BME, MBA
Global Regulatory Affairs Director
Group Head, essure & Devices
manal.morcos@bayer.com
Phone #: (908) 200-9718

ESSURE System for Permanent Birth Control

Executive Summary

For the FDA Advisory Committee Meeting

24 September 2015

Available for public disclosure without redaction

This document is an Executive Summary of the contraceptive device “Essure[®]”, prepared for the Food and Drug Administration Meeting of the Obstetrics and Gynecology Devices Panel on September 24, 2015 to Discuss the Benefits and Risks of Essure[®]

Executive Summary

For the FDA Advisory Committee meeting for Essure

Table of Contents

Table of Contents	2
Table of Text Tables	5
Table of Text Figures	5
Abbreviations	6
1. Introduction	8
1.1 Scientific Background	10
2. Regulatory History	11
3. Device and its characteristics	12
4. Clinical Development Program	15
4.1 Studies supporting the initial Premarket Approval.....	15
4.1.1 Pre- and peri-hysterectomy studies	16
4.1.2 Phase II study (STOP 10).....	17
4.1.3 Pivotal study (STOP 2000)	17
4.1.4 Safety and efficacy data in IFU from phase II and pivotal studies.....	18
4.2 Post approval studies in the United States.....	20
4.2.1 Completed post approval studies in the United States	20
4.2.1.1 ESS205 post approval study	21
4.2.1.2 ESS305 post approval study	21
4.2.2 Ongoing post approval studies in the United States	22
4.2.2.1 Transvaginal ultrasound study (ESS305, Study 16974)	23
4.2.2.2 Novasure post approval study (Study 16975)	25
4.3 Post approval studies (observational).....	25
4.3.1 SUCCES II	25
4.3.2 Intermountain Healthcare	26
4.3.3 MarketScan.....	26
5. Training	26
6. Postmarketing Safety Reporting Overview	28
7. Topics of Interest	30
7.1 Placement rate.....	31
7.1.1 Sponsor’s clinical data	31
7.1.2 Literature	32
7.1.3 Conclusion.....	32

7.2	Confirmation test	33
7.2.1	Sponsor's clinical data	33
7.2.2	Literature review	33
7.2.3	Conclusion	34
7.3	Patient satisfaction	34
7.3.1	Sponsor's clinical data	34
7.3.2	Literature	36
7.3.3	Conclusion	37
7.4	Learning curve	37
7.4.1	Sponsor's clinical data	37
7.4.2	Literature	38
7.4.3	Conclusion	39
7.5	Contraceptive efficacy	39
7.5.1	Sponsor's clinical data	39
7.5.2	Literature	39
7.5.3	Postmarketing monitoring	41
7.5.4	Conclusion	42
7.6	Pregnancy outcome	42
7.6.1	Sponsor's clinical data	43
7.6.2	Literature	43
7.6.3	Postmarketing monitoring	45
7.6.4	Conclusion	46
7.7	Unsatisfactory device location (perforation, migration or expulsion)	46
7.7.1	Sponsor's clinical data	47
7.7.2	Literature	50
7.7.3	Postmarketing monitoring	50
7.7.4	Conclusion	51
7.8	Pain (persistent/chronic)	51
7.8.1	Sponsor's clinical data	51
7.8.2	Literature	54
7.8.3	Postmarketing monitoring	56
7.8.4	Conclusion	56
7.9	Menstrual changes	57
7.9.1	Sponsor's clinical data	57
7.9.2	Literature	60
7.9.3	Postmarketing monitoring	60
7.9.4	Conclusion	61
7.10	Headache	61
7.10.1	Sponsor's clinical data	61
7.10.2	Literature	62
7.10.3	Postmarketing monitoring	62
7.10.4	Conclusion	62

7.11	Infectious complications	63
7.11.1	Sponsor’s clinical data	63
7.11.2	Literature review	63
7.11.3	Postmarketing monitoring	64
7.11.4	Conclusion	65
7.12	Allergy reactions/hypersensitivity (to nickel)	65
7.12.1	Sponsor’s clinical data	65
7.12.2	Literature	65
7.12.3	Postmarketing monitoring	66
7.12.4	Conclusion	66
7.13	Device removal	66
7.13.1	Sponsor’s clinical data	66
7.13.2	Literature	67
7.13.3	Postmarketing monitoring	68
7.13.4	Conclusion	68
7.14	Death	68
7.14.1	Sponsor’s clinical data	69
7.14.2	Literature	69
7.14.3	Postmarketing monitoring	69
7.14.4	Conclusion	69
8.	Benefit/Risk Summary	70
9.	References	72
10.	Appendices	82
10.1	Appendix 1: Essure design changes	82
10.2	Appendix 2: Literature on Essure placement rates	83
10.3	Appendix 3: Literature on pregnancy	88
10.4	Appendix 4: Literature on post-tubal syndrome	93
10.5	Appendix 5: Menstrual changes post-Essure placement	94
10.6	Appendix 6: Instructions for Use	95
10.7	Appendix 7: Patient Information Brochure	96

Table of Text Tables

Table 4–1: Overview of studies demonstrating safety and efficacy of Essure..... 16

Table 4–2: Follow-up of women up to five years (pivotal study)..... 18

Table 4–3: Cumulative failure rates (phase II and pivotal trials combined)..... 19

Table 4–4: Adverse events reported on day of placement procedure 19

Table 4–5: Adverse events by body systems, first year of reliance ($\geq 0.5\%$), (pivotal trial) ...20

Table 4–6: Completed post approval studies 21

Table 4–7: Ongoing post approval studies..... 22

Table 4–8: Confirmation test for reliance with Essure in women with bilateral placement.... 23

Table 7–1: Essure placement rate in clinical trials 31

Table 7–2: Overall comfort level with wearing the inserts during the 1st year of follow up ...35

Table 7–3: Overall satisfaction level with wearing the inserts during the 1st year of follow up 36

Table 7–4: Placement rate in first attempt by order of procedure..... 38

Table 7–5: One-year cumulative failure rates 39

Table 7–6: Women with unsatisfactory device location with at least one insert placed
(perforation, migration or expulsion)..... 48

Table 7–7: Adverse events preventing reliance in women with bilateral placement..... 49

Table 7–8: Transvaginal ultrasound findings..... 50

Table 7–9: Pain reported at follow-up visits in phase II study..... 52

Table 7–10: Pelvic pain reported at follow-up visits (pivotal study)..... 53

Table 7–11: Adverse events of abdominal or pelvic pain and related terms 53

Table 7–12: Women with diagnoses of pain, 6 months preceding sterilization procedure 55

Table 7–13: Bleeding reported at follow-up visits in phase II study 58

Table 7–14: Changes in menstrual function reported at follow-up visits (pivotal study)..... 59

Table 7–15: Adverse events related to uterine or vaginal bleeding 60

Table 7–16: Infectious complications in Essure studies 63

Table 7–17: Essure device removal in clinical trials 67

Table 8–1: Pregnancy rates by sterilization method 70

Table of Text Figures

Figure 3–1: Essure insert..... 12

Figure 3–2: Human histology specimens after device wearing 13

Figure 6–1: Adverse event reporting over time (data through 15 JUNE 2015)..... 30

Figure 7–1: Studies based in the United States..... 32

Figure 7–2: Efficacy rates in five studies based in the United States 41

Abbreviations

AEs	Adverse events
ANSM	French National Agency for Medicine and Health Products Safety
ART	Assisted Reproductive Technology
BTL	Bilateral tubal ligation
CDC	Centers for Disease Control
CDRH	Center for Devices and Radiological Health
CI	Credible Interval for the ESS305 TVU study or Confidence Interval for other studies
CO ₂	Carbon dioxide
CREOG	Council on Resident Education for Obstetrics and Gynecology
CREST	Collaborative Review of Sterilization
CRF	Case Report Form
CSS	Clinical Sales Specialist
DESH	<u>D</u> utch <u>E</u> ssure [®] versus <u>S</u> alpingectomy for <u>H</u> ydrosalpinx Trial
FDA	Food and Drug Administration
HSG	Hysterosalpingogram
HCP	Health Care Physician
IDE	Investigational Device Exemption
IFU	Instructions for Use
ISO	International Organization for Standardization
ITT	Intent-to-treat
IVF	In-vitro fertilization
MC	Medically confirmed
MedDRA	Medical Dictionary for Regulatory Activities
Non-MC	Not medically confirmed
NSAIDs	Non-steroidal anti-inflammatory drugs
OB/Gyn	Obstetrics and Gynecology
PAC	Post alternative contraception
PAS	Post approval study
PDP	Post-device placement
PET	Polyethylene terephthalate
PIB	Patient Information Brochure

PIRs	Physician Inquiry Requests
PMA	Premarket Approval
PT	Preferred Term
SD	Standard deviation
STOP	Selective Tubal Occlusion Product
TVU	Transvaginal ultrasound
UK	United Kingdom
US	United States

This document is an Executive Summary of the contraceptive device “Essure[®]”, prepared for the Food and Drug Administration Meeting of the Obstetrics and Gynecology Devices Panel on 24 SEP 2015 to Discuss the Risks and Benefits of Essure[®]

1. Introduction

Permanent birth control is the second most common method of contraception in the United States and the most common among women 21 to 45 who have completed their families. Approximately 5.6 million women in the United States (US) rely on this method of contraception (1), which offers high efficacy without the need for ongoing patient compliance. Individual contraceptive needs and desires are variable and it is essential for women to have options.

In 2002, the US Food and Drug Administration (FDA) approved Essure[®] as a novel method of permanent birth control for women. Prior to this, if a woman desired a permanent method of fertility control, her only option was to undergo a laparoscopic tubal ligation or laparotomy/mini-laparotomy. With the approval of Essure, women were given an option that took advantage of the natural pathway to the fallopian tubes through the vagina, the cervix, and the uterine cavity using a hysteroscope. This procedure avoids surgical incisions, entry into the peritoneal cavity, and the need for general anesthesia.

The FDA announced that an advisory committee meeting for Essure will be held on 24 SEP 2015 to seek expert scientific and clinical opinion on the benefits and risks of the Essure System. In this Executive Summary, Bayer will provide data from clinical studies, literature, and post market monitoring to support the safety and efficacy and the positive benefit/risk profile of Essure. The Committee Panel Review offers all sides the opportunity to examine the data, to listen to the concerns, and to reach appropriate recommendations moving forward.

In 2002, data on safety and efficacy were reviewed by the FDA as part of the Premarket Approval (PMA) process where it was determined that all criteria for approval were met and the clinical benefits outweighed the risks. Since that time, additional company sponsored postmarketing studies and clinical trials have continued to support the positive benefit/risk profile. Dozens of external studies, both prospective and retrospective, have generated data consistent with clinical trial data both in terms of efficacy and safety. This device has been studied with a total of > 10,000 women since Essure was first developed. In addition, review of spontaneously reported safety data on over one million Essure systems distributed worldwide support a product that has for the past 13 years, offered women a safe, viable alternative to laparoscopic tubal ligation.

While Essure is highly efficacious, the data show that Essure is not without risks, including those associated with the device itself or with the placement procedure. Since approval, the Instructions for Use (IFU, Section 10.6, Appendix 6) have clearly outlined the adverse events

(AEs) experienced by some women in the clinical trials as well as other potential risks associated with the device or the procedure. Both the IFU as well as the Patient Information Brochure (PIB, Section 10.7, Appendix 7) contain comprehensive information about patient selection and counseling, directions for use, contraindications, warnings, and precautions. They also contain detailed information on the procedure itself, the need for alternative contraception for 3 months, and the importance of obtaining the Essure Confirmation Test before a woman can rely on Essure for contraception.

In addition to clear, consistent messages on the benefits and risks related to the use of Essure in the IFU and PIB, a comprehensive and mandatory physician training program, including didactics, an extensive training manual, hands-on simulation training, and proctoring of the initial cases has existed since approval. Prior to the submission of the PMA, the company performed a validation study of their proposed training program and found that physicians trained with this program achieved placement rates comparable to the clinical investigators in the clinical trials (data on file). Additional support services and enhancements were made to physician training as device modifications were introduced (see Section 10.1, Appendix 1) and as new technology became available (eg, simulation software). Training enhancements will continue to be made as appropriate.

Recently, a citizen petition regarding Essure was submitted to the FDA by a law firm that has filed several Essure product liability lawsuits naming Bayer as a defendant. Several allegations regarding the clinical trials, safety, efficacy, and training were made. The FDA closed the petition and forwarded it to the Office of Compliance within the Center for Devices and Radiological Health (CDRH) as a trade complaint. The Office of Compliance is investigating the complaint and will pursue actions as deemed necessary. As stated in its response to the former petition (submitted to CDRH on 18 AUG 2015), Bayer strongly disagrees with the allegations made in the petition.

It is vital that women continue to be offered a variety of permanent contraceptive methods in order to achieve individual reproductive objectives and avoid unintended pregnancies. While Essure is not the right method for all women, Essure is an appropriate choice for those who desire a non-hormonal, non-incisional method of permanent birth control. Ultimately, the benefit/risk analysis derived from data from the clinical development program, post approval studies, external literature and extensive postmarketing safety analysis points to a product that offers women an important option.

1.1 Scientific Background

Approximately half of the 650,000 tubal ligations that occur each year in the US are performed at a time remote from pregnancy. These “interval” tubal ligation procedures are usually performed by laparoscopy. Laparoscopic bilateral tubal ligation (BTL) is regarded as a safe and effective procedure; (2) however, it requires entry into the peritoneal cavity and general anesthesia. Many of the risks of BTL are related to these factors. The risks are increased in obese women and in those who have had previous pelvic or abdominal surgery.

The first major study to evaluate the efficacy and risks associated with BTL was the US Collaborative Review of Sterilization (CREST), a large, prospective, multicenter observational study of 10,685 women conducted by the Center for Disease Control (CDC) and Prevention that reported results on long-term failure rates for several types of sterilization procedures.(3) Overall, CREST found that the failure rate for BTL was higher than previously believed (overall failure rate of 1.3%), varied by type of sterilization procedure performed, and by age at time of procedure. Of the 9475 women in the study who underwent interval laparoscopic tubal ligation, the overall standard complication rate was 1.6 %.(4) This included an unintended major surgery rate of 0.9 %.

In 1979, the CDC began surveillance of deaths attributable to tubal ligation. From 1977 through 1981, they identified 29 deaths.(5) Utilizing medical records of women reported to the Commission on Professional and Hospital Activities in 1979 and 1980, Escobedo et al (6) projected that the most reasonable case-fatality rate for tubal sterilization is ~9 per 100,000 sterilizations if all deaths associated with the procedure are considered, and between 1 and 2 per 100,000 sterilizations when only deaths that can be attributed to sterilization per se are considered.

In 2001, Huber et al (7) performed a retrospective analysis on the database of the Swiss obstetric study group and identified 20,325 cases of laparoscopic tubal ligation performed between JAN 1993 and DEC 2001. They found that the rate of major complications was 0.1%, half of which were unintended major surgery.

Approximately half of the serious injuries that occur in conjunction with laparoscopy happen while accessing the abdominal cavity using a Veress needle and/or a trocar. In a recent review of access-related bowel injury, Jansen et al (8) found an injury rate of 0.44 per 1000 gynecological laparoscopic procedures overall. Most recently, Llarena et al (9) performed a systematic review of bowel injury in all types of gynecologic laparoscopy. Over 474,000 gynecologic laparoscopies in 90 studies published between 1972 and 2014 met eligibility criteria. The rate of bowel injury varied by procedure; the rate specific to laparoscopic sterilization was 1 in 3,333 (0.03%, 95% confidential interval 0.01-0.03%). In 2006, approximately 643,000 sterilizations were performed in the US, half of which were interval procedures.(10) This would translate to approximately 95 bowel injuries annually due to laparoscopic tubal ligation.

The use of general anesthesia, insufflation of the abdomen with carbon dioxide (CO₂) gas, and the use of steep Trendelenburg position sometimes needed to access pelvic structures also increase the risk of all laparoscopic surgeries. There may be increased risks of aspiration, respiratory dysfunction, and cardiovascular dysfunction during laparoscopy.(11)

Cottin et al (1996) reported seven cases of gas embolism during laparoscopy, 2 of which were fatal. (12) All cases occurred in women with a history of previous abdominal surgery. Their review of the literature indicated that this event is rare, with an incidence of 15/113,253 (Wadhwa et al, 1978 [13]) to 7/1194 (Hynes et al, 1992 [14]) gynecological laparoscopies.

The prevalence of permanent birth control procedures in the US and the known risks associated with laparoscopic tubal ligation described above highlight the need for additional options for women. With the approval of Essure in 2002, women gained access to a method that not only avoids entry into the abdominal cavity, but also the need for general anesthesia and its associated complications.

2. Regulatory History

The Essure System for Permanent Birth Control is a Class III medical device. Essure consists of an insert and a disposable delivery system and has undergone clinical investigations since 1996. The Essure System was initially *Conformité Européenne* (CE) Mark approved on 09 FEB 2001 and received US FDA approval on 04 NOV 2002 following the FDA PMA review pathway, the most stringent regulatory review process for devices prior to marketing. Essure was the first of its kind device and as such, the FDA referred the PMA to an outside panel of experts (advisory committee) on 22 JULY 2002, for review. The panel was asked to vote on questions that discussed: safety, efficacy, labeling, training and post approval requirements. The majority of the experts (8 of 9) on the panel agreed that the Essure System should be approved; one member abstained due to personal reasons.

Essure is commercially available in the US, Canada, several European countries, Australia and several Latin/South American and Asia Pacific countries. Initially, Essure was manufactured and marketed by Conceptus, Inc., which was acquired by Bayer on 05 JUNE 2013.

Since the PMA review and approval of Essure, the FDA has conducted 8 inspections including reviews of manufacturing, quality system processes and procedures, as well as clinical data integrity. All findings observed during the inspections were appropriately addressed by the company and subsequently closed out by the FDA.

Post approval requirements

Based on the expert opinion from the PMA Advisory Committee Panel meeting in 2002, the FDA required a post approval study (PAS) as a condition of the original Essure PMA approval. In addition to the PAS (to evaluate bilateral placement rate for newly trained physicians), the FDA required Conceptus to continue to follow all women in the phase II and

the pivotal trial through five years (see Section 4.1.4) in order to obtain longer term safety and effectiveness outcomes data.

PMA annual reports and adverse device reporting is a condition of PMA approval which has allowed for continuous interactions with the FDA throughout the life cycle of the device.

Since 2002, the FDA has reviewed and approved many post approval reports and 39 supplements to the PMA; including 15 approved PMA supplements to update the IFU (Section 10.6 [Appendix 6]) and the PIB (Section 10.7 [Appendix 7]).

3. Device and its characteristics

Device description: The Essure system is comprised of a disposable delivery system and a wound-down insert. The insert consists of a super-elastic nitinol (nickel titanium alloy) outer coil and a medical grade 316L stainless steel inner coil wrapped in polyethylene terephthalate (PET) fibers.

The inner coil is attached to the outer coil by a tin-silver alloy solder at the distal end of the insert. The wound-down (unexpanded) insert is approximately 4 cm in length and 0.8 mm in diameter. When released, the outer coil expands up to 2.0 mm in diameter, conforming itself to the varied diameters and shapes of the fallopian tube (Figure 3–1).

Figure 3–1: Essure insert


Mechanism of action: Once Essure is placed in the fallopian tube, the dynamic outer coil expands and secures the insert in place in the uterotubal junction region. Tubal occlusion is attributed to the space filling design of the device and the predictive, progressive and occlusive tissue response. The PET fibers cause tissue in-growth into and around the insert, facilitating insert retention and pregnancy prevention. During this three month period of in-growth, a woman must use an alternative form of contraception.

Histological data gathered in pre-hysterectomy clinical trials have shown that the tissue response to the Essure device is occlusive in nature, providing long term anchoring of the device.⁽¹⁵⁾ These studies demonstrated that the tissue in-growth response was predictable, occurred in all fibered specimens collected, and was localized to the inner portion of the

fallopian tube wall, with normal tubal architecture present within 5 mm of the distal end of the device. There was no evidence that fibrosis induced by the device extended into the wall of the fallopian tube or caused peritubal adhesions or serositis. The tissue response varied according to the elapsed time after placement (Figure 3–2). Acute inflammation was predominant in early specimens whereas chronic inflammation and fibrosis was more extensive in the fallopian tubes of women who had worn the device for longer than 8 weeks. The acute inflammatory and low-level chronic inflammatory responses were consistent with other devices that have used PET fibers. At 12-13 weeks post device placement, the Essure insert became completely incorporated into the fallopian tube with a proliferation of fibrotic tissue.(15) By 30 weeks, dense fibrosis replaced the tubal lumen; only scant acute inflammatory cells were present.

Figure 3–2 provides the histology of human fallopian tube.

Figure 3–2: Human histology specimens after device wearing

BU = device in intramural-isthmic portion; fiber is visible on inner coil.
 CU = cross-section of fallopian tube 5 mm distal to device location.
 Source: Valle RF et al, 2001 (15)

Device Biocompatibility

The Essure system (insert and delivery catheter) passed extensive biocompatibility testing in accordance with ISO 10993 (International Organization for Standardization). The following testing has been performed: cytotoxicity, sensitization, genotoxicity, muscle implantation, vaginal irritation, mutagenicity, subchronic toxicity, and acute systemic toxicity. Based upon biocompatibility testing and human clinical studies, the Essure device is well tolerated by the human body. In addition, all materials used in Essure are well characterized with a long history of clinical use in, and similar in form to, approved and marketed long-term implantable medical devices.

The PET fiber was chosen for this application due to its success in causing tissue in-growth in other medical applications, such as prosthetic arterial grafts, aneurysm coils, and other long-term implants since the 1940's.(16, 17, 18, 19, 20) Polyethylene terephthalate has been demonstrated to produce an immediate local tissue response characterized by macrophages, fibroblasts and collagen.(21) When used for stent-grafts in humans, the tissue reaction observed tends to be composed predominantly of macrophages and monocytes. A moderate foreign-body tissue reaction is also elicited by the PET fibers.(22) This tissue response has been noted to peak between two and three weeks, after which the tissue response slowly resolves during a 10-week period, resulting in extensive fibrosis which occludes and anchors the PET material and any associated medical device.

The unique combination of shape, memory and superelasticity properties, coupled with its proven biocompatibility and corrosion resistance, has made nitinol a suitable material for medical and dental applications since the mid 1980's. Products such as orthodontic arch wire, vena cava filters, guide wires, stents, distal protection devices, heart valves, support membranes in a catheter wall, and bone anchors are all common nitinol applications.

The nickel ions in nitinol are tightly bonded to titanium, and the nitinol used is processed so that the entire alloy surface is covered with a protective layer of titanium oxide; both of these characteristics act to minimize nickel ion release. An in vitro study was performed to determine the amount of nickel leached after the Essure insert was placed in saline solution for various lengths of time up to 6 months. The maximum nickel leaching rate was 0.14 µg/day.(23) This is far less than that for other approved nitinol implants which range from 0.42 – 8.40 µg/day.(24) In addition, it is over 2000 times less than the amount of nickel that is ingested from water and food each day (~300 µg/day).

Based on the literature research, testing data, and widespread human use of the device materials, the Essure device testing met the requirements of ISO 10993 for a surface device with permanent contact with mucosal membrane tissue and the FDA determined that the safety profile was acceptable for permanent birth control.

Design changes

Since the PMA was originally approved by the FDA in 2002, Essure has undergone several design changes to improve performance and these changes were approved by the FDA prior to marketing (Section 10.1, Appendix 1). The original approved model was the STOP device. The first US launched device was the ESS205 model which employed the same insert as the STOP device and incorporated delivery system changes. Subsequent changes to the Essure system were mainly made to the delivery catheter based on feedback gathered through direct customer input and/or post-market monitoring. For the ESS305 model, minor changes were also made to the insert to accommodate changes to the delivery system. These changes did not affect the critical design aspects of Essure and the safety and efficacy of the device. The current ESS305 model has been commercialized since 2007.

4. Clinical Development Program

Data for the phase II and pivotal studies in this document is based on the final study report with 5 years of follow-up.

The “Essure[®] Permanent Birth Control System” or “Essure[®] System” was previously known as a “Selective Tubal Occlusion Procedure” or “STOP” device or the ESS005 Model. The device is indicated for permanent birth control (female sterilization) by occlusion of the fallopian tubes.

All studies in the clinical development program of Essure are summarized in this section as follows:

- Premarket studies supporting the initial PMA are presented in Section 4.1 (Table 4–1).
- Completed post approval studies in the US are presented in Section 4.2.1 (Table 4–6).
- Ongoing post approval studies in the US are presented in Section 4.2.2 (Table 4–7).
- Post approval observational studies are presented in Section 4.3.

4.1 Studies supporting the initial Premarket Approval

For the initial PMA submission, 4 studies were included with a total of 903 women undergoing the Essure procedure (not including the 4 women without case report forms (CRFs) from Germany in STOP 01 study). Of these, 643 women (194 in the phase II study and 449 in the pivotal study) with bilateral placement who were told to rely on Essure were followed for safety and effectiveness.

Studies that supported the initial PMA approval of Essure[®] are presented in Table 4–1.

Table 4–1: Overview of studies demonstrating safety and efficacy of Essure

Study number	IDE number	Study title	Study dates / countries	Women enrolled / Women underwent placement procedure
Studies supporting the initial Premarket Approval				
Phase IA / STOP 01	IDE G960052	STOP Device Placement Feasibility Study Using Hysteroscopic Visualization.	July 1998 – Nov 2000 / Australia, Europe, Mexico, US	99 / 95
Phase IB / STOP 06	IDE G960206	Evaluation of the Safety and Principles of Operation of the Selective Tubal Occlusion Procedure (STOP) Device in Women who are scheduled to undergo a hysterectomy.	Oct 1998 – Dec 2001 / Mexico and US	63 / 63
Phase II (STOP 10)	IDE G980152	Phase II Study – A Multi-Center Clinical Trial to Evaluate the Safety and Effectiveness of the STOP Device to Prevent Pregnancy in Women Who are Seeking Permanent Contraception.	Nov 1998 – June 2000 / Australia, Belgium, Spain, USA	269 / 227
Pivotal (STOP 2000)	IDE G000055	A Multi-Center Clinical Trial to Demonstrate the Safety and Effectiveness of the STOP Device in Providing Permanent Contraception.	May 2000 – 05 Dec 2007 / Australia, Belgium, Spain, UK, US	657 / 518

IDE = Investigational Device Exemption; STOP = Selective Tubal Occlusion Procedure; UK = United Kingdom; US = United States (of America).

4.1.1 Pre- and peri-hysterectomy studies

Two studies (STOP 01 and STOP 06) were conducted in women who were scheduled to undergo hysterectomy for benign reasons in order to demonstrate the feasibility and acute safety of Essure placement and to provide support for the postulated mechanism of action and the ability of the inserts to occlude the fallopian tubes. In addition, the STOP 06 study evaluated detachment of the insert from the delivery wire, woman's tolerance of and recovery from placement procedure, stability of the insert within the tube until hysterectomy, tubal occlusion within 24 hrs to 12 weeks, local tissue response and the effect of fiber on the ability to create a local tissue response.

In STOP 01 study, 99 women were enrolled, however, 95 women were evaluated (4 women were excluded from data analysis due to the Investigator's death at a site in Germany and the CRFs from this site could not be obtained).

In STOP 06 study, 63 women underwent placement attempt with a Gamma design and 54 women had at least one device placed.

In these two studies, device handling and placement were found to be acceptable and the short-term wearing of the insert (1 to 30 weeks) was found to be acceptable with no reported side effects. The placement procedure was found to be safe with minimal post-procedure discomfort and minimal AEs and was able to be accomplished with local or regional anesthesia.

In addition, tubal occlusion based on hysterosalpingogram (HSG) occurred in 92/93 (99%) of tubes containing inserts at the time of hysterectomy, supporting the postulated mechanism of action, and the tissue in-growth reaction was predictable, localized to the insert, and did not result in adverse clinical sequelae.

4.1.2 Phase II study (STOP 10)

This was a prospective, multi-center, international phase II study of women seeking permanent contraception. Four design concepts (Beta 1, Beta 2, Beta 3 and Gamma) were clinically evaluated during the development of the Essure insert and were tested in this study. The Gamma design was ultimately selected for PMA approval.

The objectives of the study were to evaluate the woman's tolerance of and recovery from the insert placement procedure, the woman's tolerance of the placed inserts, the safety of the placement procedure, the long-term safety and stability of the placed inserts, and the effectiveness of the inserts in preventing pregnancy.

The study with the Gamma design enrolled 269 women, 227 of whom underwent an insert placement procedure. None of the women in this study became pregnant while relying on the Gamma design for contraception. One woman relying on an earlier (Beta) design did become pregnant. This was reported to FDA at the time. Specific data on device placement, confirmation testing and subject satisfaction are presented in Section 7.

Safety data: Based on the final 5-year study report (dated 06 JAN 2006), a total of 21/227 (9.3%) women reported AEs, 12/227 (5.3%) related to episodes of period pain (dysmenorrhea), ovulatory pain or changes in menstrual function. The remaining 9 reported AEs: perforation (7 [3.1%]), expulsion 1 [$<1\%$]), and other (1 [$<1\%$]) were due to unsatisfactory device location. Further information on perforations, expulsions, pain and device removals is provided in Section 7. (Note: The data described in this paragraph uses the total number of women with device placement attempt (227) as the denominator. In this study, 206 women had at least one device placed).

4.1.3 Pivotal study (STOP 2000)

The objectives of this pivotal study were to determine the bilateral placement, reliance, pregnancy, and AE rates in women undergoing Essure placement procedure. Women in this study were followed one week post-device placement (PDP); 3 months PDP; 3, 6, and 12

months post-alternative contraception (PAC); and 18 months to 5 years PAC as a part of FDA's required post approval commitment.

In total, 657 women enrolled in the study. Ninety-nine women subsequently changed their mind about participating and 40 were considered screen failures. Therefore, 518 women underwent hysteroscopy and 507 underwent attempted Essure placement. At the conclusion of this study, 364 women with bilateral devices and 2 women with unilateral devices (with unicornuate uterus or with contralateral proximal tubal occlusion on HSG) that were relying on Essure for contraception had completed their required 5-year follow-up. In addition, 20 women who were followed for safety only (ie, were not relying on Essure for contraception) completed the 5-year follow up period ([Table 4-2](#)).

Table 4-2: Follow-up of women up to five years (pivotal study)

Follow-up visit post-alternative contraception	Number completing visits to 05 DEC 2007 data cut-off date
12-month	461
24-month	435
36-month	423
48-month	403
60-month	386

Adverse events on the day of insert placement procedure occurred in 16 (3.0%) women. With the exception of one woman who had an adverse reaction to pain medication and required overnight observation, all events resolved prior to discharge. No event occurred in more than 1% of women.

Post-operative recovery was reported as "uneventful" in 313 (60%) women. The most common post-operative events were cramping and pain occurring in 31.0% and 13.5% of women respectively, followed by nausea in 9.7% of women.

Four luteal phase pregnancies were reported (pregnancies occurring prior to insert placement but not detected on the day of placement). However, none of these 4 women became pregnant while relying on Essure and all conception dates were before the insert placement procedure. No other pregnancies were reported in this study (see [Section 7.6](#)).

Overall comfort with long-term wearing of the inserts was high (see [Section 7.3](#), Patient Satisfaction). Adverse events rated by the investigator as at least possibly related to Essure during the first year of reliance are presented in [Table 4-5](#).

4.1.4 Safety and efficacy data in IFU from phase II and pivotal studies

Phase II (STOP 10) and pivotal (STOP 2000) studies form the basis for safety and efficacy data in the FDA approved IFU ([Section 10.6](#), Appendix 6). The objective of these studies was

to demonstrate safety and effectiveness of Essure system in providing permanent contraception.

In the two studies combined, the effectiveness of Essure inserts was assessed in 635 women who relied on Essure for contraception. There were no reported pregnancies (Table 4–3).

Table 4–3: Cumulative failure rates (phase II and pivotal trials combined)

One-year^c N=635	Two-year^c N=605	Three-year^c N=586	Four-year^c N=567	Five-year^c N=567
0%	0%	0%	0%	0%
95% CI, (0 - 0.10%) ^{a, b}	95% CI, 0 - 0.20%) ^{a, b}	95% CI, (0 - 0.30%) ^{a, b}	95% CI, (0 - 0.40%) ^{a, b}	95% CI, (0 - 0.50%) ^{a, b}

N = number of women; CI = confidential interval

^a 95% confidence intervals are based on a “constant-hazard” exponential failure-time model, whose parameter is determined by the total number of woman-months accumulated during the trial.

^b Combined effectiveness data obtained using Bayesian statistics.

^c The number of women “N” were considered to have completed follow-up at 1 year if patient contact occurred at ≥ 11 months, 2 years if contact occurred at ≥ 23 months, 3 years if contact occurred at ≥ 35 months, 4 years if contact occurred ≥ 47 months and 5 years if contact occurred at ≥ 59 months.

Adverse events occurring during the hysteroscopic placement procedure in these studies are summarized in the IFU (Section 10.6, Appendix 6) presented in Table 4–4. Data presented below are taken from the IFU and reflect percentage of events as a function of the number of procedures recorded.

Table 4–4: Adverse events reported on day of placement procedure

Adverse event / side effect	Phase II trial	Pivotal trial
	Number of procedures	
	N=233 (100%)	N=544 (100%)
Cramping	NA ^a	161 (29.6%)
Pain	2 (0.9%)	70 (12.9%)
Nausea/vomiting	NA ^a	59 (10.8%)
Dizziness / light headed	NA ^a	48 (8.8%)
Bleeding / spotting	NA ^a	37 (6.8%)
Other ^b	NA ^a	16 (2.9%) ^b
Vaso-vagal response	2 (0.9%)	7 (1.3%)
Hypervolemia	NA ^a	2 (0.4%)
Band detachment	3 (1.3%)	2 (0.4%)

^a Data not collected

^b Other for the pivotal trial includes: ache (3), hot/hot flashes (2), shakiness (2), uncomfortable (1), weak (1), profuse perspiration (1), bowel pain (1), sleepiness (1), skin itching (1), loss of appetite (1), bloating (1), allergic reaction to saline used for distension (1).

Table 4–5 summarizes AEs with >0.5% incidence rated as at least "possibly" related to the insert or procedure during the first year of reliance in the pivotal trial (approximately 15 months post-device placement). The most frequently reported AEs in the first year were back pain, abdominal pain/cramps and dyspareunia. Data below are taken from the IFU and are presented as the number of events as divided by the total number of subjects with at least one insert in place.

Table 4–5: Adverse events by body systems, first year of reliance (≥ 0.5%), (pivotal trial)

Adverse events by body system	Number (%)^a
Total number of women placed with at least one insert	476 (100%)
Abdominal	
Abdominal pain/abdominal cramps	18 (3.8%)
Gas/bloating	6 (1.3%)
Musculo-skeletal	
Back pain/low back pain	43 (9.0%)
Arm/leg pain	4 (0.8%)
Nervous/Psychiatric	
Headache	12 (2.5%)
Premenstrual Syndrome	4 (0.8%)
Genitourinary	
Dysmenorrhea/menstrual cramps (severe)	14 (2.9%)
Pelvic/lower abdominal pain (severe)	12 (2.5%)
Persistent increase in menstrual flow ^b	9 (1.9%) ^b
Vaginal discharge/vaginal infection	7 (1.5%)
Abnormal bleeding – timing not specified (severe)	9 (1.9%)
Menorrhagia/prolonged menses (severe)	5 (1.1%)
Dyspareunia	17 (3.6%)
Pain/discomfort - uncharacterized:	14 (2.9%)

^a Percentages reflect the number of events divided by the number of women in the trial. When numerous episodes of the same event were reported by one woman, each report was counted as a separate event. Therefore, percentages may over-represent the percentage of women who have experienced that event.

^b Eight women reported persistent decrease in menstrual flow

4.2 Post approval studies in the United States

4.2.1 Completed post approval studies in the United States

Two studies were conducted and completed after approval of Essure (Table 4–6).

Table 4–6: Completed post approval studies

Study	PMAS	Objective	Date/country	Enrolled / underwent procedure
ESS205	P020014	Post approval study for newly trained physicians	May 2003 – Oct 2005 / US	585 / 564
ESS305	P020014	ESS305 post approval study	Aug 2007 – May 2009 / US	619 / 593

PMAS = postmarketing approval study;

4.2.1.1 ESS205 post approval study

This study was intended to document the bilateral placement rate for newly trained physicians using the ESS205 system. The data were to be used to evaluate bilateral placement rate to newly trained physicians. Data collected included the following:

- Rates of successful bilateral placement of the Essure System at first attempt; and
- Identification of factors predictive of failure to achieve bilateral placement.

A total of 585 women were enrolled and 564 underwent placement attempt. Of these women, 194 placement attempts were made with the Gamma (ESS005) catheter and 370 placement attempts with the coil catheter (ESS205).

Adverse events (during and after the placement procedure) were reported in 15/585 (2.6%) women, inclusive of one woman in whom one insert perforated the left uterine fundus and embedded in the omentum, causing pain. The peritoneal portion of the insert was laparoscopically removed and each tube was banded with Falope rings. The woman reported increased pain of unknown etiology and was to be scheduled for diagnostic studies and was lost to follow-up. No unanticipated adverse device effects were reported.

4.2.1.2 ESS305 post approval study

The Essure (ESS205) system was modified to the ESS305 model based on feedback gathered through direct customer input and postmarket monitoring (see Section 10.1, Appendix 1 for design changes).

The objective of this study was to validate the ESS305 device changes in both the newly trained and experienced physicians and it is the basis for the placement rate in the approved IFU (Section 10.6, Appendix 6) from 2010 through JUNE 2015 when the IFU was modified to include the TVU/HSG confirmation test algorithm. This study initially planned to enroll 800 women but was terminated early, after enrollment of 619 women. Upon agreement from the FDA, this smaller sample size was deemed to be sufficient to demonstrate statistically significant results for the study endpoints.

A total of 625 women were enrolled and 6 women were excluded from the placement rate calculation: 5 women withdrew consent prior to placement attempt and one woman underwent placement procedure with the incorrect device (ESS205). Therefore, 619 women are included in the ITT population with the placement rate of 95.8%. Essure inserts were placed through the operating port of the hysteroscope in 612 women. Of these 612 women, 593 (96.9%) had successful bilateral placement. Data on potential differences in placement rates between newly trained and experienced physicians is presented in Section 7.4.

Device issues and malfunctions: Essure device issues were reported in 14 of the placement procedures (one procedure may have had multiple events reported). In five of these cases, the distal tip of the insert was bent. This is not a device malfunction per se since tip bending is a feature of the device (rather than cause damage when it comes into forceful contact with the uterus or the fallopian tubes). Bent tips can occur in the commercial use of Essure and may be caused by rough handling when placed into the introducer, advanced through the hysteroscope, and when manipulating the device inside the uterus.

Ten of the reported device issues represent device malfunctions. These include six insert deployment or detachment issues, one thumbwheel retraction difficulty, and three unspecified device issues. Each of these issues have been captured in Conceptus' product surveillance system and investigated under standard quality system procedures.

Safety data: Six of 619 women were reported to have AEs during and after the Essure placement procedure. One woman was hospitalized after hysteroscopy due to a uterine perforation by the hysteroscope (and not the Essure device). None of these AEs represent unanticipated adverse device effects.

4.2.2 Ongoing post approval studies in the United States

Two studies are currently ongoing (Table 4-7): the transvaginal ultrasound (TVU) study and the NovaSure PAS.

Table 4-7: Ongoing post approval studies

Study number	PMA number	Objective	Dates and countries	Women enrolled / undergoing procedure
ESS305 TVU Study (Study 16974)	P020014	Use of Transvaginal Ultrasound to Confirm Essure® Insert Placement in Women: Demonstration of Effectiveness (Essure® TVU study).	May 2011 Canada, Netherlands, Spain, US, 10 year study	620 / 597
NovaSure PAS (Study 16975)	P020014	To Evaluate the Effectiveness of Essure Post-NovaSure Radiofrequency Endometrial Ablation Procedure Following a Successful Essure Confirmation Test	Nov 2012 US	Currently enrolling (220 planned)

PMA = Premarketing approval; TVU = transvaginal ultrasound; PAS = post approval study; US = United States

4.2.2.1 Transvaginal ultrasound study (ESS305, Study 16974)

The TVU study is a currently ongoing phase IV study being conducted in order to support the use of TVU as a first line confirmation test in the US. Enrollment was completed on 11 OCT 2012. Women in this study are to be followed for 10 years following discontinuation of alternative contraception after being instructed to rely on Essure for contraception. The TVU/HSG confirmation test algorithm was approved on 29 JUNE 2015 based on the 1-year effectiveness data from this study.

The objectives of the study are to:

- Evaluate the effectiveness of the Essure procedure when the TVU/HSG algorithm is used for confirmation testing (co-primary endpoint)
- Evaluate the intent-to-treat reliance rate 3 months following TVU/HSG confirmation test algorithm (co-primary endpoint)
- Evaluate subject satisfaction with the TVU confirmation test

A total of 620 women were enrolled into the study; of these, 23 women withdrew prior to undergoing the Essure placement and 597 women underwent the Essure procedure and constitute “the intent-to-treat” (ITT) population.

All women in the ITT population (n = 597) underwent insert placement procedure. Details on placement status are provided in the topics of interest section (Section 7.1). Bilateral placement was achieved in 582 of 597 (97.5%) of women.

Overall, 547 women were told to rely on Essure based on the TVU/HSG confirmation test algorithm (Table 4–8). The reliance rate in the ITT population was 91.6% (547/597) with the posterior credible interval of 89.2% - 93.6%. The posterior probability that the reliance rate is > 86.7% is 99.99%, which exceeds the preset probability threshold of 97.5% for the study.

Table 4–8: Confirmation test for reliance with Essure in women with bilateral placement

Test	Told to rely, number (%)		
	Yes	No	All
Total	547 (100)	35 (100)	582 (100)
Transvaginal ultrasound alone	470 (85.9)	4 (11.4)	474 (81.4)
Transvaginal ultrasound followed by Hysterosalpingogram	47 (8.6)	6 (17.1)	53 (9.1)
Hysterosalpingogram alone	30 (5.5)	2 (5.7)	32 (5.5)
No test performed ^a	0	23 (65.7)	23 (4.0)

^a women withdrew prior to Transvaginal ultrasound or Hysterosalpingogram

One year effectiveness: As of 12 AUG 2014 (data cutoff for 1-yr effectiveness report), 4 pregnancies were reported, 3 pregnancies during the first year and 1 pregnancy during the second year of wearing of Essure insert (no additional pregnancies have been reported through JULY 2015). The total 1-year reliance follow-up time at data cutoff was 6217 women-months. The estimate of the monthly pregnancy rate is 3/6217 or 0.048% (95% CI: 0.01%, 0.13%) and the 1-year pregnancy rate is 0.67% (95% CI: 0.16%, 1.53%). The probability that the upper credible limit of the 1-year pregnancy rate is <2.1% is >99.99%. This is greater than the study's predefined probability cutoff of 97.5%, so this co-primary endpoint was met.

An annual report was submitted to the FDA with a cut-off date of 19 JUN 2015. As of this cut-off date, 473 subjects had a 2-year follow-up visit and 90 subjects had a 3-year follow-up visit. As of this date, no additional pregnancies have been reported and the follow-up is 16,897 women-months. Based on the life table analysis, the 2-year pregnancy rate is 0.754%.

The 4 pregnancies (none ectopic) mentioned above were reported from 4 different sites, 3 in the US, and one in the Netherlands. All 4 women were compliant with the 3 month study visit, underwent confirmation testing within the protocol designated time window, and were told to rely on Essure based upon TVU alone. A review of the root cause for Essure insert efficacy failure in these 4 cases appears to be perforation in 2 women and unsatisfactory insert location in the other 2 women. These root causes are consistent with those described with the commercial use of Essure.

The secondary efficacy objective of this interim report is to evaluate subject satisfaction with the TVU confirmation test. The analysis of subject satisfaction was reported as very satisfied or somewhat satisfied by a total of 520/528 (98.5%) women; neither satisfied or dissatisfied by 3 (0.6%) women; and somewhat dissatisfied by 3 (0.6%) women.

Safety data: As of a data cutoff date of 19 JUNE 2015 (for the recent PAS Annual Report), a total of 493 AEs were reported in this study. The severity of most (n=280) AEs was considered mild. Moderate and severe AEs were reported for 113 and 44 events, respectively. Thirteen AEs were reported as 'definitely' related to the Essure system. Hospitalization was required for 11 AEs that occurred in 8 women: pelvic inflammatory disease (PID), pleomorphic adenoma, adnexa uteri pain, intervertebral disc protrusion, abdominal pain, nausea, chest pain, pelvic pain, anxiety, cerebrovascular accident and myocardial infarction.

The most common AEs by Medical Dictionary for Regulatory Activities (MedDRA) preferred term (PT) were: pelvic pain (n=28), menorrhagia (n=23), abdominal pain (n=16), dysmenorrhea (n=15), vaginal hemorrhage (n=14), uterine hemorrhage (n=9), and ovarian cyst (n=9).

As of the data cutoff date (19 JUNE 2015), there were 28 SAEs in 20 women. Three of these events (PID, adnexa uteri pain, and pelvic pain) were categorized as “possibly related” and one event (pregnancy) as “probably related” to treatment with Essure.

4.2.2.2 Novasure post approval study (Study 16975)

The purpose of the study is to evaluate the effectiveness and safety of the Essure System when a NovaSure endometrial ablation procedure is performed following a successful Essure Confirmation Test. The primary variable is the occurrence of confirmed pregnancy at 1 year and 3 years among women relying on Essure inserts for permanent birth control when NovaSure EA is performed following a successful Essure Confirmation Test. The study began enrollment on 12 NOV 2012 and enrollment is currently ongoing. To date, no pregnancies and no new or unexpected safety findings have been reported.

4.3 Post approval studies (observational)

Several post approval *observational* studies have been conducted by Conceptus, one of which is currently ongoing (SUCCES II).⁽²⁵⁾ Three published observational studies: ESTHYME study ⁽²⁶⁾, SUCCES I ⁽²⁷⁾ and CONFORM study ⁽²⁸⁾ were conducted between 2002 and 2006.

In addition, two ongoing studies compare the rates of hysterectomy following hysteroscopic sterilization with that of laparoscopic sterilization (Intermountain Healthcare and MarketScan).

4.3.1 SUCCES II

The SUCCES II study is an ongoing multicenter, observational study being conducted in France that was initiated in 2008 by Conceptus France and which will be concluded in 2016. A total of 2,600 women were enrolled in the study. A total of 2,281 women reached the 3-months follow-up, and continued follow-up is planned for 5 years following placement through 2016.

As outlined in the study protocol, the study’s primary objective is “to assess patient’s satisfaction of the Essure procedure at 5 years.”

The secondary objectives aim to assess:

- Pain experienced when placement is performed without anesthesia and the potential risk factors of pain,
- Methods used and the interpretation of the 3-month control test and
- Effectiveness, complications and satisfaction of the Essure procedure at short-, mid- and long-term (3, 12, 24 months and 5 years).

Selected aspects of the interim results of SUCCES II, namely, the results of the placement phase conducted between 2008-2011, were published as an abstract in the Journal of Minimally Invasive Gynecology in 2012. (25) In 2,366 (93.3%) cases, the first placement attempt was successful and 97% women reported being satisfied or very satisfied.

Subsequent to assuming responsibility for AE reporting from Conceptus France, Bayer extracted all paper-based case report forms (CRFs) and transferred all data on AEs and potential AEs into its safety database, consistent with its standard reporting processes. Of the 1,059 AE cases identified, 42 were classified as reportable and disclosed to the FDA, among other regulatory authorities, including the French National Agency for Medicine and Health Products Safety (ANSM).

To date, 5 pregnancies have been reported resulting in a calculated pregnancy rate of 0.19% and a contraceptive efficacy rate of 99.8%. None of the pregnancies were ectopic.

The study data as available to date confirm the high safety and efficacy of Essure and of the Essure placement procedure. The SUCCES II study provides further confirmatory evidence of a favorable benefit-risk profile, including considerations of the alternatives in this indication (ie, laparoscopic tubal sterilization).

4.3.2 Intermountain Healthcare

This is a retrospective analysis using electronic medical record database in the US. A 10-year observational period from 01 APR 2005 – 31 MAR 2015 is planned. This is a descriptive study with the primary objective of describing hysterectomy rates in women having undergone hysteroscopic permanent birth control and women having undergone tubal ligation.

4.3.3 MarketScan

This is a retrospective cohort study using a US claims database of women undergoing hysteroscopic or laparoscopic tubal ligation procedures. The observational period is from JAN 2010 – DEC 2012. The analysis will be performed separately for individuals with and without a pre-existing diagnosis of pain prior to the procedure. The rates of post-procedural pain and surgical procedures across all groups will be examined.

5. Training

Training of physicians on the use of Essure is important to ensure its safe and appropriate use. The IFU (Section 10.6 [Appendix 6]) for Essure clearly states: “Device to be used only by physicians who are knowledgeable hysteroscopists; have read and understood the Instructions for Use and Physician Training Manual; and have successfully completed the Essure Training Program, including preceptoring in placement until competency is established, typically 5 cases.” Hysteroscopy has been a part of obstetrics/gynecology (OB/Gyn) residency training for decades. As with any procedure learned in residency, each physician must critically self-

assess their level of comfort and experience with hysteroscopy before undergoing training on Essure.

As part of the initial PMA, detailed materials were provided to the FDA describing the training program. In addition, a validation of the training program was performed in 2002, demonstrating that newly trained physicians achieved placement rates comparable to those in the clinical trials (data on file). Continuous re-evaluation, as well as new data and technological advances in simulation training, have led to updates to the program.

Specific training to use Essure, or the “Clinical Pathway,” consists of multiple steps. The first step consists of didactic review on appropriate patient selection, indications, contraindications, warning and precautions. The device itself, the clinical trials leading to approval, as well as insert placement steps and the confirmation test are also reviewed in detail. In order to assist in this stage, a Clinical Resource (“Training Manual”) is available in both printed and electronic format. Bayer last updated and revised the Clinical Resource in 2014. At the current time, didactic training involves a one-on-one discussion between the physician and a Clinical Sales Specialist (CSS). All CSS’s undergo comprehensive training that lasts for approximately three months. Bayer is in the process of also creating an interactive online training program for physicians.

After completion of didactics, physicians then move onto simulation. The EssureSim™ is a validated computer-based system that takes a physician through several Essure-based scenarios.(29, 30) Feedback is immediately available and simulated cases become progressively more difficult. A uterine model and hysteroscope may also be used to practice placing Essure. Both deployable and nondeployable demonstration devices are available. A physician completes a minimum of 5 proctored cases before being eligible to have their final case proctored and receive their Certificate of Completion. These cases must be completed with either the CSS or an Essure trained Physician Proctor in attendance. Data on the learning curve for Essure supports this approach (see Section 7.4 on Learning Curve). A physician may continue to have a CSS or a Physician Proctor present for cases beyond completion of the Clinical Pathway if so desired. Additionally, office staff of a physician who will perform the placement procedure in their office may receive training that is specific to their responsibilities.

All physicians have access to a wide variety of materials on the Essure MD website (www.EssureMD.com) and this site is updated regularly as new information becomes available. These materials include sample informed consent documents and procedure notes, checklists, patient education materials, and additional training on the Essure Confirmation Test.

Several other resources and training opportunities are available to physicians. These are summarized below:

- *Physician Inquiry Requests:* Physicians are encouraged to contact Bayer with Essure related questions or to report AEs. A dedicated Bayer Medical Information employee is available during business hours. Questions may also be escalated to the Medical Director (a board-certified OB/Gyn) as necessary or required.
- *National Consultancy Network:* Questions that are out of the scope of Bayer's ability to answer are forwarded to small group of physicians with extensive experience with Essure by the Essure Medical Director. All AEs of the case are captured and forwarded to the pharmacovigilance department; patient management decisions remain the responsibility of the querying physician.
- *Advanced Training Summits:* In conjunction with a global manufacturer and distributor of endoscopes, Bayer currently offers a 1½ day didactic, simulator and uterine model program for physicians who have completed the Clinical Pathway. These programs emphasize troubleshooting, interpretation of the HSG confirmation test, patient identification and counseling, and techniques for ensuring patient safety in the office.
- *Radiologist Training Program:* A program tailored for radiologists was developed in order to provide information and guidance on performance of the modified-HSG and appropriate interpretation of confirmation test films.
- *Residency Training:* As of January, 2013, the Council on Resident Education for OB/Gyn (CREOG) changed their guidelines to now state that all OB/Gyn residents must know how to perform hysteroscopic sterilization. In order to assist resident programs with this task and ensure that residents receive a thorough education on the procedure, a six-module resident education program was developed. This program relies on faculty input and participation and is designed to complement existing hysteroscopy training that residency programs are mandated to deliver. Access to the EssureSim™ is provided to residents as part of these programs.
- *Proctors:* Essure proctors are physicians with extensive experience with the Essure procedure. Proctors are fully trained on compliance and AE reporting as well as what their roles and responsibilities will be. Proctors may observe other physicians during or after completion of the Clinical Pathway. Physicians and appropriate staff may also go to the proctor's office to observe cases and safe placement techniques in an office setting.

6. Postmarketing Safety Reporting Overview

The safety profile of Essure is well known from clinical trials, observational research and postmarketing reporting. The amount of postmarketing reporting has increased over time which is reflective of the increasing postmarketing exposure to Essure as it can be estimated based on sold Essure kits (each containing two inserts).

Besides this, an increase of case reports since Q3 2013 has been observed. The increase coincides with Bayer's acquisition of Conceptus in JULY 2013 and with Bayer assuming responsibility for AE case handling and case reporting in OCT 2013.


Over this period, the number of medically confirmed case reports has increased only slightly, whereas the number of non-medically confirmed case reports has increased disproportionately. The difference in case numbers before and after Q3 2013 is predominantly explained by an increase in non-medically confirmed reports, mostly from the US.

Figure 6–1 presents the annualized number of individual safety (case) reports of AEs as included in the sponsor's global safety database through 15 June 2015, regardless of their relationship with Essure. Cases are presented by the year of initial receipt, date for medically confirmed events and non-medically confirmed reports separately. The latter represent cases received from patients/consumers without further confirmation by a Health Care Physician.

Several elements are likely to have influenced the observed trends globally and in particular from the US. Bayer closely monitors and continuously evaluates the safety profile of Essure. This includes monitoring of spontaneous reports of AE and associated reporting rates as well as publications from medical journals and information from general and social media. There are secular trends of increased reporting observed in postmarketing AE reporting at Bayer across the portfolio, that impacts on reporting with Essure. Besides this, the acquisition of Conceptus by Bayer led to increased public attention to Essure and an increase in interest in the product's safety profile by attorneys and their affiliates at least in part focused on potential avenues for liability, which was noticeable in social and other media, but also in reporting to Bayer call centers and websites. Bayer responded to these trends with activities that aimed at facilitating reporting from non-medically confirmed sources; for example, active online listening and outreach programs.

Also, as a matter of general policy and in the attempt to facilitate comprehensive assessment by receiving Health Authorities and regulatory bodies, Bayer favors an assessment of reportability, which may be more comprehensive than with other companies. It is, however, Bayer's understanding that Conceptus has always fulfilled regulatory reporting requirements for Essure. In parallel, reporting from social media has increased for Essure, facilitated through the availability of public mobile AE reporting applications, including the FDA released MedWatcher app, which are assisting consumers in submitting AE reports to the FDA (31). Limitations of postmarket AE reporting are well known, accepted and described in the FDA regulation.(32) Consistently, the review of the information available from postmarketing indicates that the types and rates of AEs that were received are within the expected range for the Essure inserts and the procedure associated with their placement. No new safety information has been identified to date that would indicate a change in the positive benefit/risk ratio.

Figure 6–1: Adverse event reporting over time (data through 15 JUNE 2015)

MC = medically confirmed; Non-MC = not medically confirmed.

Detailed evaluations of cumulative information available to date were conducted in support of the preparations for the advisory panel meeting. It should be noted that in parallel to taking over the AE reporting obligation, Bayer transferred all postmarketing legacy AE reports from Conceptus into the Bayer global pharmacovigilance database. While no individual follow-up reports were filed for this transfer step, all data compiled by Conceptus are now included in aggregate analyses.

The information derived from these postmarketing evaluations are consistent with the known safety profile of Essure and the information as labeled in the IFU (Section 10.6, Appendix 6). Although increases in reporting rates were observed in recent years since Conceptus was acquired by Bayer, the numerical increase can be explained by factors influencing postmarketing safety reporting. Reporting rates remain compatible in line with the data generated in clinical studies.

7. Topics of Interest

Several topics of interest were identified in preparation for this meeting. The following sections identify these topics and present data from the clinical trials, the published literature, and postmarketing monitoring. Within the literature, very few studies include a comparison group which makes it difficult to put incidence rates and relatedness into context. In addition, many studies were retrospective reviews; this needs to be considered when interpreting results. However, the data is important to consider when assessing the benefit/risk profile. In order to help minimize potential biases, we chose to focus primarily on US-based studies that were not company sponsored or funded.

7.1 Placement rate

In order to be effective, the Essure inserts must be properly placed across the uterotubal junction. There are a small number of cases where individual pelvic anatomy (eg, cervical stenosis, intracavitary lesions) precluded even attempting placement of Essure inserts. An alternate form of contraception (eg, tubal ligation) is needed in these women. In other cases, other anatomic or physiologic circumstances (eg, laterally placed tubes, tubal spasm) make a particular placement more challenging. In some of these instances, a second attempt at placement may be successful. In the small number of cases where successful bilateral placement is not achieved, a woman must be counseled about her contraception needs moving forward.

7.1.1 Sponsor's clinical data

The phase II and pivotal clinical trials for Essure were conducted with the ESS005 (STOP) device. This was replaced by the ESS205 device in 2003. Both of these devices required the physician to physically unscrew the delivery catheter from the insert at the completion of the procedure. In 2007, the ESS205 device was replaced with the ESS305, the currently marketed device. With this device, the detachment procedure is carried out by using a button and wheel on the delivery handle, both of which are activated with the thumb. In the ESS305 post approval study, the Essure bilateral placement rate improved. In addition, the time to place the device decreased. A summary of the placement data from the clinical trials is presented in [Table 7-1](#).

Table 7-1: Essure placement rate in clinical trials


Study	Women enrolled	Women undergoing hysteroscopy (ITT)	Women undergoing Essure placement attempt	Bilateral placement rate (ITT)	Device group bilateral placement rate	Number of unilateral placements
Phase II (STOP 10)	269	227	227	200/227 (88.0%)	200/227 (88%)	6
Pivotal (STOP 2000)	657	518	507	464/518 (89.6%)	464/507 (91.5%)	12
ESS205 PAS	585	585	564	525/585 (89.7%)	525/564 (93.1%)	27
ESS305 PAS	619	619	612	593/619 (95.8%)	593/612 (96.9%)	12
ESS305 TVU (Study 16974)	620	597	594	582/597 (97.5%)	582/594 (98.0%)	5

ITT = intent to treat; PAS = post approval study; TVU = transvaginal ultrasound

7.1.2 Literature

Several published studies have also investigated placement rates, with data collected prior to 2007 looking at the ESS205, and later studies looking at the ESS305 (see Figure 7–1). Additional details about these studies (eg, study design, year of data collection, number of placement attempts, comments) can be found in Section 10.2, [Appendix 2]. Overall, the results from these studies are consistent with those found in the clinical trials. Literature has also found that placement rates appear to be independent of the level of training,(33, 34) location of service (operating room verses office), anesthesia verses no anesthesia, and degree of patient obesity.(34)

Figure 7–1: Studies based in the United States

CI = confidence interval; x = number of bilateral placements; n = total number of women.
Chapa, 2015 (35); Leyser-Whalen, 2015 (40); Howard, 2013 (42); Anderson, 2013 (34); Levie, 2010 (48); Savage, 2009 (41); Shavell, 2009 (36); Shavell, 2008 (39); Levie, 2006 (33).

7.1.3 Conclusion

Placement rate is very important when physicians counsel women prior to the attempted procedure, as women need to be aware of the possibility placement may not occur. Still, with the currently available ESS305 device, both the clinical trial data as well as several external studies show a placement rate of at least 95%.

7.2 Confirmation test

Once the Essure procedure is complete, a woman must use alternate contraception and undergo a confirmation test after three months. Only after a satisfactory result on the confirmation test can a woman rely on Essure for permanent birth control. Need for a post-procedure confirmation test is also part of the algorithm for male sterilization procedures (ie, vasectomy).

The IFU summarizes the data obtained from the MAUDE database on pregnancies in the commercial setting (see Section 10.6, Appendix 6). From worldwide commercial launch in 2001 through the end of 2010, 660 pregnancies in the US were reported. In 32% of cases (n=213), patient non-compliance (eg, failure to use alternate contraception or return for Essure Confirmation Test) was identified. This points to the absolute importance of patient compliance with the confirmation test.

7.2.1 Sponsor's clinical data

In the phase 2 study, HSG was performed in the 200 women who initially had bilateral placement (and in 2 women with unilateral placement) to assess tubal occlusion at the 3-month visit. Insert location was assessed using either HSG or ultrasound and 194/200 women had satisfactory insert location and 187 women had successful bilateral occlusion at 3 months and were able to rely on Essure for contraception. Additionally, the 7 subjects with patency noted at 3 months were found to have bilateral tubal occlusion at repeat 6 month HSGs and were therefore told to rely.

In the pivotal study, of the 507 women undergoing attempted device placement, 464 had bilateral device placement; 456 underwent 3 month follow-up (HSG) to assess location of insert and occlusion and 453 were told to rely on Essure for contraception. Of the 49 remaining subjects, 31 had no devices placed, 3 were lost to follow-up prior to the 3 month visit, 3 did not have an HSG but relied (protocol violations), 2 had unsatisfactory placement detected before 3 months and 12 had unilateral placement.

In the TVU study (Study 16974), of the 597 women who underwent Essure placement attempt, 559 (93.6%) underwent confirmation testing and 547 were told to rely. Of the 38 women not undergoing confirmation testing, 23 were lost to follow up or withdrew and 15 had non-bilateral placement.

7.2.2 Literature review

In European studies, where several methods for assessing the ability to rely on Essure are available, compliance with the post-procedure confirmation test is very high (Povedano et al 2012, 96.8% [4108/4306] [37]; Ubeda et al, 2004 93% [75/81] [38]). In the US, compliance with HSG is variable. In a retrospective review of medical records in an inner city clinic population, Shavell et al (39) reported a compliance rate of 12.7% (10/79). In contrast, in a

primarily Hispanic population,(33) 92 out of 98 women (94%) completed their HSG. Additional studies found rates of 85% (243/286) in a publicly funded clinic,(40) 87% (739/850) at Kaiser Permanente Northern California,(41) and 53% (70/132) in a publicly insured population.

Several factors are associated with increased compliance with HSG. Anderson et al (34) found that women were 1.51 times more likely to undergo HSG if the Essure procedure was performed in the office setting versus the operating room. Women were 2.05 times more likely to undergo HSG if they had private insurance compared with Medicaid. Howard et al (42) found that women 35 years and older had an almost 4-fold higher odds of HSG compliance. Women who had their Essure procedure performed at a campus with a dedicated protocol to improve compliance had almost a 4-fold higher rate of compliance (odds ratio=3.67). This finding points to the fact that programs can be put into place to increase compliance.

Guiahi et al (43) reported that when a staff nurse was dedicated to tracking women's HSG results, compliance went from 78% pre-intervention to 90.9% post-intervention. Finally, Mahmud et al (44) reported that the implementation of an electronic reminder increases post-Essure HSG adherence, particularly in a resident clinic patient population. Pre-reminder compliance with HSG in a faculty practice was 90% (62/69) while compliance in the resident clinic was 47% (32/68). Post-reminder, compliance for the faculty practice remained high (19/20, 95%) while compliance in the resident clinic increased to 72% (39/54).

7.2.3 Conclusion

The Essure Confirmation Test is a vital part of the Essure procedure and must be completed in order to assess a woman's ability to rely on the Essure inserts for contraception. Some studies indicate poor compliance with the confirmation test. It appears that measures can be put into place, such as patient reminders, in order to increase women's compliance with this aspect of the procedure.

7.3 Patient satisfaction

While patient satisfaction is subjective and the tools used to assess this endpoint vary from study to study, collection of this data adds to the understanding of how women perceive a particular procedure or experience, and therefore, this is an important aspect of clinical research.

7.3.1 Sponsor's clinical data

Pivotal study (STOP 2000)

In the pivotal study, a woman's comfort with wearing STOP insert and her overall satisfaction with the STOP device were queried by the investigator at the 3-month post-device placement

visit and at 3, 6, 12, 18, 24, 36, 48 and 60 months post-alternative contraception follow-up visits. Comfort was rated as “good” to “excellent” by at least 95% women at all visits. Satisfaction was “somewhat satisfied” to “very satisfied” in at least 97% women at all visits through 5 years.

ESS305 TVU study (Study 16974)

In the ongoing ESS305 TVU study (Study 16974), both comfort and satisfaction with wearing the inserts were high at each of the 3 follow-up visits for which data is currently available: 3-months, 6-months, and 12-months.

Comfort level with the wearing of inserts was described as “excellent”, “very good” and “good” by 495 (98.4%) women at 12-month visit ([Table 7–2](#)).

Table 7–2: Overall comfort level with wearing the inserts during the 1st year of follow up

Overall comfort level ^a	3-month visit N=558 n (%)	6-month visit N=536 n (%)	12-month visit N=503 n (%)
Excellent, very good, good	552 (98.9)	518 (96.6)	495 (98.4)
Excellent	351 (62.9)	374 (69.8)	366 (72.8)
Very good	151 (27.1)	114 (21.3)	108 (21.5)
Good	50 (9.0)	30 (5.6)	21 (4.2)
Fair	4 (0.7)	8 (1.5)	4 (0.8)
Poor	0	1 (0.2)	1 (0.2)
Not reported	2 (0.4)	9 (1.7)	3 (0.6)

N, n = number of women

^a Question was asked “How does the woman rate her comfort with the wearing of inserts”?

Satisfaction with the wearing of inserts was described as very satisfied and somewhat satisfied by 495 (98.4%) women at the 12-month visit ([Table 7–3](#)).

Table 7–3: Overall satisfaction level with wearing the inserts during the 1st year of follow up

Overall satisfaction level ^a	3-month visit N=558 n (%)	6-month visit N=536 n (%)	12-month visit N=503 n (%)
Very satisfied and somewhat satisfied	547 (98.0)	513 (95.7)	495 (98.4)
Very satisfied	484 (86.7)	479 (89.4)	462 (91.8)
Somewhat satisfied	63 (11.3)	34 (6.3)	33 (6.6)
Neither dissatisfied nor satisfied	9 (1.6)	12 (2.2)	4 (0.8)
Somewhat dissatisfied	0	2 (0.4)	1 (0.2)
Not reported	2 (0.4)	9 (1.7)	3 (0.6)

^a Question was asked “How does the woman rate her overall satisfaction with the wearing of inserts”?

Similar levels of comfort and satisfaction were recorded throughout the duration of the 5 year follow-up with no trend of an increase or decrease in comfort or satisfaction over time.

7.3.2 Literature

A number of published studies also include a measure of patient satisfaction. Some studies have assessed patient satisfaction with the procedure itself. In a prospective case series of 1630 women who underwent Essure placement in Spain from JAN 2003 to JUNE 2006, 91% of women rated the procedure as a “10” (highly satisfied) on a visual analog scale, and no woman rated it below an “8.” In addition, 97% of women said they would recommend the procedure to a friend.

Additional studies included the confirmation test in the overall patient satisfaction rating as they considered this part of the overall procedure. Sinha et al (45) sent a questionnaire to 84 women post-confirmation test, which was returned by 76 (90%) women. The majority of women (70/76; 92%) were satisfied with the overall experience. Of the 6 women reporting that they were not satisfied, one stated that this was due to “extreme” pain during the procedure, three stated that this was related to their experience of the HSG, and the other two gave no reason. Most of the women (69/76; 91%) reported that they would recommend the procedure to a friend or relative. Miño et al (46) also assessed satisfaction at 3 months in 857 women. Overall patient satisfaction was rated as “very high” by 806 women (94%) and “high” by 51 (6%). None of the women reported being dissatisfied. Ubeda et al (38) performed a prospective, observational study at a private university hospital in Spain. Eighty-five cases were performed between JULY 2002 and JULY 2003. Satisfaction was characterized by good or excellent in all but 3 women (96%) in the immediate post-operative period and, at 3 months, all women said they would advise a friend or a sister to undergo the procedure.

Few studies offer long-term follow up on patient satisfaction. Anderson et al (47) performed a retrospective review of 61 women who underwent Essure placement between 2002 and 2007.

A questionnaire was sent to these women one to five years post-procedure and was returned by 50 women (81.9%). Average follow up period was 23 months (7 to 67 months). All women who returned the questionnaire expressed overall satisfaction with the procedure and all said they would recommend the procedure to others. Levie et al (48) contacted women who received Essure in the outpatient setting 13 weeks to 1 year post procedure to assess satisfaction. They were able to reach 176/209 women (84%). Within this group, the mean satisfaction score was 4.7 (standard deviation [SD] 0.71) on a scale of 1 to 5 (5 being most satisfied). In addition, 164 women (93%) reported they would have the procedure again and 173 (98%) stated that they would recommend the procedure to a friend.

In the only comparative study, Duffy et al (49) performed a cohort controlled comparative study of Essure versus laparoscopic tubal ligation in the United Kingdom (UK) and assessed satisfaction 90 days post procedure (N=89 total). Overall, satisfaction was high with both groups: 94% of the women in the Essure group and 80% of the women in the laparoscopic group were either “very” or “somewhat” satisfied at 90 days post-procedure. With regard to speed of recovery, 100% women in the Essure group were “very satisfied” verses 80% in the laparoscopic group.

7.3.3 Conclusion

Both the clinical trials and external studies demonstrate overall satisfaction (>92%) with the Essure placement procedure as well as ongoing use of the inserts. When compared to laparoscopic tubal ligation at 90 days post procedure (in one study), Essure compared very favorably with regard to overall satisfaction, speed of recovery, and choice of chosen procedure.

7.4 Learning curve

7.4.1 Sponsor’s clinical data

Data on this topic (learning curve with Essure placement) is available from two large prospective studies conducted by the sponsor, the original pivotal study and the ESS305 PAS.

Pivotal study (STOP 2000)

A learning curve analysis based on procedure time and a second analysis based on placement rates were conducted as part of this study.

Procedure time decreased significantly after the first 5 cases, and slowly thereafter for both experienced and non-experienced healthcare providers. However, there was a smaller decrease in hysteroscopy times with increased experience.

Bilateral placement rates did not improve substantially with increased experience with the device (Table 7-4).

Table 7–4: Placement rate in first attempt by order of procedure

	All	Procedure number for center			
		1 to 5	6 to 10	11 to 20	More than 20
Placement after first procedure	N	N (%)	N (%)	N (%)	N (%)
Not bilateral	43	11 (12.8%)	5 (6.8%)	11 (10.0%)	16 (7.6%)
Bilateral	438	75 (87.2%)	69 (93.2%)	99 (90.0%)	195 (92.4%)
All	481	86 (100.0%)	74 (100.0%)	110 (100.0%)	211 (100.0%)

ESS305 Post approval study

The ESS305 post approval study was designed to compare successful bilateral placement rates during the first placement attempt of newly trained versus experienced physicians. The bilateral placement rate in the experienced physician group was 342/349 (98.0%) and the bilateral placement rate in the newly trained physician group was 223/232 (96.1%). The difference in placement rate (1.9%) was not statistically significant (Fisher’s exact test $p = 0.2009$). A random effects model showed similar results ($p = 0.1855$). The conclusion of this study was that there was no significant difference in placement rates between newly trained and experienced physicians when using the ESS305 system.

Levie and Chudnoff published the results of this study in 2011.⁽⁵⁰⁾ They describe that, an average of 7.7 procedures (± 3.3 SD) were performed per physician. Of the physicians, 39 were experienced Essure users and 37 were newly trained. As noted above, there was no statistically significant difference in placement rates when comparing novices and experienced physicians. Experienced physicians, however, were approximately 3 minutes faster in performing placement procedures. The authors concluded that this study confirmed that the learning curve for the Essure procedure is steep and that once a physician performs 3 to 5 proctored cases, they can be relatively reassured that they will have high placement rates.

7.4.2 Literature

Review of the literature has provided one additional external study that specifically looked at the learning curve for placement of Essure. Janse et al⁽²⁹⁾ examined learning curves for 15 gynecologists over the period of JAN 2005 to DEC 2009 in 10 different community and teaching hospitals in the Netherlands. A total of 631 women underwent Essure placement during this time; the median number of women per gynecologist was 29 (range: 9 to 116). All procedures were performed with oral diclofenac for analgesia. The study examined procedure time, pain score during the procedure, successful bilateral placement rate, and number of complications. Only procedure time showed a statistically significant difference with experience and thus a learning curve. The authors suggested that for the other parameters, “The Essure procedure is not difficult for a gynecologist experienced in performing hysteroscopies, so immediate high success percentages can be achieved.”

7.4.3 Conclusion

Data from the clinical trials as well as external studies demonstrate that novice and experienced physicians achieve similar rates of successful placement when appropriately trained. The amount of time needed for placement does appear to decrease with experience using the product. Given this data, the training requirement for at least 5 proctored cases prior to sign off appears to be appropriate.

7.5 Contraceptive efficacy

While no method of contraception can offer 100% efficacy, Essure offers women a highly effective method of permanent birth control. Once a woman has Essure inserts successfully placed and her confirmation test confirms appropriate placement and tubal occlusion, data from both the clinical trials (see [Table 7-5](#)) and external published literature both show that Essure is highly effective.

7.5.1 Sponsor's clinical data

[Table 7-5](#) displays pregnancy data (failure rates) in the phase II study, pivotal study, and the TVU study. One additional pregnancy was reported in the ESS305 TVU study during the second year of follow-up and no pregnancies were reported during the additional 5 years of follow up in the phase II and pivotal studies.

Table 7-5: One-year cumulative failure rates

	Phase II	Pivotal	Combined (original data)	ESS305 TVU
	N = 193	N = 439	N = 632	N = 547
1 year	0% (95% CI: 0-1.53%) ^a	0% (95% CI: 0-0.68%) ^a	0% (95% CI: 0-0.47%) ^a	0.67% (95% CI: 0.16-1.53%) ^b

^a CI = Confidence interval; ^b CI = Credible interval

TVU = transvaginal ultrasound

7.5.2 Literature

Clinical trial conditions may not be reproducible in real world settings, so it is important to assess the efficacy in a variety of clinical situations. The largest US-based published study took place at Kaiser Permanente, Northern California.⁽⁴¹⁾ This was a retrospective review of 884 women, 118 physicians at 30 different facilities, with placements from JAN 2004 to DEC 2006 and follow up through DEC 2008. Initial placement attempt was successful in 850 women (96.2%). The loss to follow up prior to HSG was 13% (115 women). Bilateral occlusion was found in 687/739 women who underwent HSG. Pregnancy was identified in 8 women: one never returned for HSG, four had an HSG showing at least one patent tube, and


three had HSG interpreted as bilaterally occluded. Among the three women who conceived after the HSG was read as occluded, subsequent review suggested that none of the inserts were properly placed.

A more recent study by Anderson et al (34) was a retrospective cohort analysis of 638 women who underwent hysteroscopic sterilization between 01 JULY 2005 and 30 JUNE 2011. Approximately two-thirds of the procedures performed in the operating room and slightly less than half of the procedures performed in the office setting were performed by resident physicians in training. There was no statistical difference ($p=.10$) in the placement rate between experienced surgeons and resident physicians. In an analysis looking at all pregnancies, 3 pregnancies were reported prior to HSG confirmation, which were judged to be unrelated to the device effectiveness. There were no pregnancies reported in women after HSG confirmation of bilateral tubal occlusion, with follow-up of at least one year. In addition, there were no known pregnancies among women who did not complete the HSG.

The largest retrospective cohort study to date took place in France.(51) The purpose of this study was to compare the rates of pregnancy among women who underwent Essure hysteroscopic sterilization versus tubal ligation in France between 2006 and 2010. During the study period, French hospitals performed 109,277 tubal sterilization procedures: this included 39,169 Essure procedures and 70,108 laparoscopic tubal ligation procedures. A Cox model was performed and after adjusting for age, women with Essure became pregnant at a significantly lower rate than women with laparoscopic tubal ligation (0.36% versus 0.46%, respectively, hazard ration = 0.62 [0.40-0.96]).

Figure 7–2 summarizes five published papers and one published abstract in the US. Additional information on these studies as well as studies performed elsewhere in the world are available in Section 10.3, Appendix 3.

Figure 7–2: Efficacy rates in five studies based in the United States

CI = confidence interval

Anderson, 2013 (34); Howard, 2013 (42); Deraieu, 2012¹ (Published abstract) (52); Savage, 2009 (41); Shavell, 2009 (39); Levie, 2006 (33).

A single published paper (53) describes a potential efficacy rate of Essure of approximately 90%. Essure <90%. Rather than a prospective or retrospective study, this paper is a theoretical Markov model comparing the efficacy of Essure to tubal ligation. While the model is data based, it does not take into account all of the data supporting the efficacy of Essure, including the single largest study to date by Povedano et al (37) (N=4306) showing an efficacy of >99%. In addition, the model contains several assumptions (eg, women who receive Essure, who have an average age in their mid-thirties, are assumed to have an annual fecundity rate of 85% for 10 years if they are unable to rely on the device) that can be disputed.

7.5.3 Postmarketing monitoring

A search of the global safety database through 15 JUNE 2015 revealed 2207 reports describing pregnancy after Essure procedure. Assuming a cumulative total exposure of approximately one million devices sold, the reporting frequency of pregnancy overall is approximately 0.21%. Most reports were spontaneous reports, followed by literature as described above.

A small number of cases on intended pregnancy could be identified from the literature (n=57) describing off-label use of Essure as a treatment for hydrosalpinx prior to Assisted Reproductive Technology (ART) typically with intended fresh or frozen embryo transfers, resulting in a reporting rate of 0.006%. No increased risk for the developing fetus secondary to the presence of Essure is described for these intended pregnancies.

7.5.4 Conclusion

No method of contraception is 100% effective. Clinical trials, extensive published literature, and postmarketing monitoring of Essure all point to a product with very high efficacy. Essure is a multi-step procedure: not only must the device be placed, the woman must comply with a confirmation test 3 months later. Data from several external studies demonstrate that even when including pregnancies in women who were told not to rely on Essure for contraception and pregnancies occurring in the first 3 months, contraceptive efficacy is still high. The pregnancy rate with Essure continues to compare favorably with that of surgical sterilization and is lower than that seen with some forms of temporary contraception.

7.6 Pregnancy outcome

Pregnancy, whether planned or occurring as a consequence of contraceptive failure, may be accompanied by a number of adverse outcomes. Sporadic miscarriage is the most common complication of pregnancy, and approximately 25% of all women who become pregnant will experience a pregnancy loss.⁽⁵⁴⁾ The vast majority (approximately 80%)⁽⁵⁵⁾ are early, occurring well before 12 weeks of gestation. Sporadic miscarriage after this time complicates no more than 1 to 2% of pregnancies. Maconochie et al⁽⁵⁶⁾ found that the risk of miscarriage rises sharply after age 35, with a 75% increase for mothers aged 35 to 39 years and a 5-fold increase where the mother is aged 40 and above (relative to mothers aged 25 to 29 years).

Preterm labor, preterm premature rupture of membranes, and premature delivery are later complications of pregnancy. Preterm birth rates have been reported to range from 5% to 10% of live births in some developed countries and have shown a dramatic rise over the past 20 years. A Canadian study found a prevalence of preterm premature rupture of membranes in tertiary university healthcare centers of 2.3% and an incidence of 2.8% of all pregnancies.⁽⁵⁷⁾

The background incidence of ectopic pregnancy in the US is estimated at 1.3 – 2.1%.^(58, 59) Damage to the fallopian tubes, including prior tubal ligation, greatly increases the risk of ectopic pregnancy. The CREST study found that with some methods of tubal ligation (eg, bipolar cautery), ectopic pregnancy rates may be as high as 65% of all pregnancies that occur after these procedures.⁽⁶⁰⁾ Finally, in an overview of US background incidences⁽⁵⁸⁾ which used data from the CDC and Prevention looking at more than 67 million pregnancies, the overall rate of pregnancies resulting in fetal death was 0.5%. For pregnancies conceived with ART, the 2012 CDC national report found a stillbirth rate (death after 20 weeks gestation) of 0.6% (CDC 2012, http://www.cdc.gov/reproductivehealth/data_stats/).

When a pregnancy occurs with the Essure device in place, there are two ways to investigate the pregnancy outcome: first, by looking at unintended pregnancies in women who had Essure inserts placed for contraception, and second, by examining desired pregnancies conceived via in-vitro fertilization (IVF) or related ART procedures.

7.6.1 Sponsor's clinical data

Pivotal study (STOP 2000)

Outcome of luteal phase (pre-placement) pregnancies: Four women had insert placement in the second half of the menstrual cycle and despite negative urine pregnancy tests on the day of the procedure, these women conceived during this cycle. None continued the pregnancy: one aborted spontaneously, and the other three terminated the pregnancies (one medically and 2 surgically). All four of these women continued to use alternative contraception for the remaining 3-month time period, underwent the 3-month HSG, and were found to have bilateral occlusion and inserts in satisfactory location. All were ultimately able to rely on the Essure inserts for contraception, and are included in the "Reliant Population" for the STOP 2000 study. None conceived while relying on Essure.

One woman with successful bilateral placement underwent IVF with the inserts in situ, became pregnant and subsequently delivered a healthy baby. After her pregnancy, she continued to rely on her inserts for permanent birth control.

ESS305 TVU study (Study 16974)

Four pregnancies were reported in this study, 3 during the first year and 1 during the second year of wearing of Essure insert. Two were electively terminated and two ended in spontaneous first trimester abortions.

7.6.2 Literature

Veersema et al (61) performed a retrospective analysis of pregnancies reported with Essure in the Netherlands (2002-2010). They found 26 cases of unintended pregnancies. Of these, 17 women chose to terminate the pregnancy, 8 women had term, healthy deliveries, and the results for one are unknown.

Several additional case reports of unintended pregnancy with Essure exist. Ory et al (62) reported a conception 6 months post-HSG which resulted in a term delivery in a 30-year old. Moses et al (63) reported two cases. In the first, a 38 year old woman had 8 coils present on the right and 13 on the left. She conceived approximately 6 months after placement (she did not complete the Confirmation Test) and delivered a term pregnancy. In the second case, a 35 year old woman had documentation of devices extending from each tubal ostium into the uterine cavity on HSG. She became pregnant approximately one year later. She carried the pregnancy to 19 weeks, and chose to terminate. Hastings-Tolsma et al (64) reported a

placement with 9 coils on the right and 3 on the left. Her pregnancy occurred almost two years later and was carried to term. Finally, Arjona et al (65) reported a complicated placement which resulted in pregnancy over two years later. The pregnancy was uncomplicated.

A large proportion of pregnancies following a failed Essure procedure are terminated, making it difficult to assess pregnancy outcomes. For this reason, studies that followed pregnancies conceived with ART methods after Essure placement are useful in studying pregnancy outcomes. It must be kept in mind that women undergoing IVF may represent a very different population than women with unintended pregnancies. They may be of a different age, and they often have other concomitant issues that may adversely affect pregnancy. These are also clearly desired pregnancies. Overall, only a small fraction of pregnancy relate to intended pregnancies.

The publications describing off-label use of Essure for treatment of hydrosalpinx in the context of ART do not report an increased risk for the developing fetus secondary to the presence of Essure. No pregnancy reports from sources other than literature have been received in this context.

Pregnancy rates in a normal fertile population average approximately 20 to 25% per month and should be considered when interpreting pregnancy rates from ART. Rosenfield et al (66) published the first report of a successful IVF pregnancy and live birth in a morbidly obese woman with extensive pelvic adhesions following the placement of Essure for hydrosalpinx. Since that time, several additional studies have been published.

Arora et al (67) published a systematic review and pooled analysis on Essure for management of hydrosalpinx prior to in vitro fertilization. They included 11 studies for a total of 115 women, highlighting the fact that off-label use of Essure for the purpose of hydrosalpinx treatment is very rare. Overall, Essure was successfully placed in 96.5% of women and tubal occlusion was confirmed in 98.1% of cases. In vitro fertilization subsequent to Essure resulted in 54 pregnancies for a pregnancy rate of 39% per embryo transfer. Of the 54 pregnancies, 14 resulted in miscarriages and the remainder resulted in live births (25.9%) giving a live birth rate per embryo transfer of 29%. Legendre et al (68) conducted a national survey of 45 French hospitals and performed a retrospective analysis of 41 women with unilateral or bilateral hydrosalpinges and successful Essure placement. Mean number of visible coils at the end of the procedure was 1.61 (± 1.6); 27 women (65.9%) had fewer than 3 coils. Overall, the clinical pregnancy rate was 40.7% (22/54). Among clinical pregnancies, there were 7 early spontaneous miscarriages (31.8%) and one intrauterine death due to fetomaternal platelet alloimmunization. There were 14 live births (including 4 sets of twins) for a live birth rate per transfer of 25.9% (14/54).

Mijatovic et al (2009) reported on a case of IVF pregnancy that resulted in pre-term premature rupture of membranes, chorioamnionitis, pre-term delivery at 24 weeks, and subsequent death of the child. Hysteroscopic evaluation of the uterus post-delivery showed total tissue

encapsulation of the left device and only the tip of the right device was visible. Her delivery was attributed to an incompetent cervix and in their follow-up report, Mijatovic et al (2012) reported a subsequent near-term delivery after placement of a primary cervical cerclage.

Most recently, Ozgur et al (71) performed a retrospective study of 26 women who had Essure placed due to laparoscopic contraindications versus 76 women who had laparoscopic tubal ligation. In the Essure group, no more than 4 coils were left protruding into the uterus. Live birth rate per assisted reproductive procedure was 23.8% in the hysteroscopy group and 32.1% for the laparoscopy group. These results were deemed to be “highly comparable.” The results of the currently ongoing Dutch Essure[®] versus Salpingectomy for Hydrosalpinx (DESH) Trial will add to the ongoing body of literature.

7.6.3 Postmarketing monitoring

In order to address what impact Essure may have on pregnancy outcomes in individuals who decide to continue their pregnancy, all reported pregnancies and outcomes were reviewed for any pattern or type of abnormality which could indicate any influence of Essure on the developing pregnancy. This evaluation includes 57 intended pregnancies reported with ART and 2150 pregnancies not associated with ART. Information regarding pregnancy outcome is available for 671 of these cases while 422 cases are reported as ongoing pregnancies. Since some cases can have more than one outcome such as preterm labor and live birth, the total number of individual events will be greater than the number of cases with a pregnancy outcome. Elective termination of pregnancy was reported in 175 cases. Spontaneous miscarriage was reported in 198 cases and ectopic pregnancy in 75 cases. Considering that the average age of Essure users is towards the higher end of reproductive age at which the incidence of these events increase significantly and that the more serious adverse outcomes tend to be preferentially reported, the rates are within the expected range.

Preterm events: There were a total of 31 preterm events reported, including premature rupture of membranes, preterm labor and preterm delivery. Twelve of these are from ART procedures and 19 from unintended pregnancies.

Premature rupture of membranes: Three cases of premature rupture of membranes were reported in the medical literature from ART. In all cases the authors stated that they did not see the premature rupture of membranes as related to Essure. One case reports a healthy term delivery while no outcomes are reported for the other two.

Preterm labor: Seven cases of preterm labor were reported of which 5 reported term deliveries of healthy infants. The two other cases were due to polyhydramnios and cervical insufficiency which are well established causes of preterm labor. The outcome was not reported. Six cases are from unintended pregnancies and one is from ART.

Preterm delivery: Of the eight cases associated with ART there were 5 healthy infants, 2 neonatal deaths (at 19 and 24 weeks of gestation), and one without outcome information.

Thirteen non-ART cases included 4 births with healthy infants some of whom required transient respiratory support, 2 neonatal deaths (at 25 and 26 weeks gestation) and no outcome information for the remaining 7 cases.

Fetal demise and stillbirth: Five of the seven reported cases are derived from ART related publications. One case which was also reported as premature delivery at 19 weeks, had Essure removed and suffered another stillbirth leading to the author's conclusion that these events were independent of Essure. Two of the cases of fetal demise were reported as secondary to alloimmunization. Two additional ART cases report delivery at 21 and 27 weeks without further information. The 2 non-ART cases include delivery at 22 weeks gestation and at an undetermined gestational age respectively.

Live births: There were 217 reports of live births of which 107 were reported without further information; 82 were reported as healthy. There were 16 cases reporting transient neonatal abnormalities. Most of these were related to respiratory assistance and one fractured clavicle during delivery. Congenital anomalies were reported in 4 cases. Two cases reported genetic deletion without further specification and one case reported developmental delay. In addition to the fatalities discussed above associated with prematurity, there was one case of anencephaly.

There are very limited data on pregnancy outcomes with Essure as would be expected with a highly effective contraceptive. The ART pregnancies are very closely monitored. Published reports show the outcomes of these pregnancies to be consistent with that of the overall ART population. In summary, the reported data regarding pregnancy outcomes with Essure does not demonstrate any increase in the rate of adverse pregnancy outcomes.

7.6.4 Conclusion

The body of evidence from the literature, which is primarily reporting on closely monitored ART pregnancies, is reassuring due to the lack of reports of premature rupture of membranes and preterm labor and delivery. Similarly, all of the Essure data from postmarketing monitoring show AEs (eg, miscarriage, stillbirth, preterm delivery, fetal anomalies) to be well within the background rate for a similarly aged population.

7.7 Unsatisfactory device location (perforation, migration or expulsion)

The ideal location of the Essure insert is symmetrically spanning the interstitial segment of the fallopian tube. The location is primarily determined by the position of the Essure device at the time of deployment during the placement procedure. Unsatisfactory locations can be divided into two groups.

- Locations where the insert is not sufficiently far into the fallopian tube to cause occlusion, but rather too much of the insert is within the uterine cavity.

- Locations where the insert has advanced too far and is distal to the desired location.

The first group is detected by the presence of too many coils of the insert in the uterine cavity by either HSG or an ultrasound. Some of these inserts will spontaneously be expelled into the vagina. Those requiring intervention may be removed hysteroscopically. When women do present with symptoms they usually involve bleeding and/or dysmenorrhea which is characteristic of the uterus expelling a foreign body.

The second group represents inserts that are in the distal fallopian tube, or are perforating the fallopian tube or uterus. These cases may require a laparoscopic or laparotomy approach when insert removal is required. Embedment represents an insert that is fixed in location. This most commonly describes a device that has been deployed into the myometrium and may be extending either into the uterine cavity or into the peritoneal cavity.

7.7.1 Sponsor's clinical data

Data from the pre and peri-hysterectomy studies are presented separately as these studies were with individuals requiring hysterectomy for gynecological abnormalities rather than individuals meeting the criteria for use of Essure for permanent birth control.

STOP 01: Three subjects experienced a perforation during the procedure. One of the perforations occurred with the use of the support catheter that was associated with high number of perforations in other studies and was discontinued in 1999. Another occurred in a woman who had a prior tubal ligation (an exclusion criterion for the study). No etiology was identified for the third.

STOP 06: Three perforations were noted at the time of hysterectomy, 2 occurred in women in whom the since discontinued support catheter was used. Women who experienced the perforations reported no discomfort or difference in tolerance to the inserts than women without perforation.

[Table 7–6](#) presents the number of subjects with unsatisfactory device location in clinical trials for contraception.

Table 7–6: Women with unsatisfactory device location with at least one insert placed (perforation, migration or expulsion)

Study	Women with placement attempt, (n)	Women with perforation, n (%)	Women with expulsion, n (%)	Unsatisfactory device location, other ^c , n (%)	Total unsatisfactory location, n
Phase II (STOP 10)	227	7 ^a	1/206 (0.5%)	1/206 (0.5%)	9
Pivotal (STOP 2000)	507	5/476 (1.1%)	14/476 (2.9%)	12/476 (2.5%)	31
ESS305 TVU (Study 16974)	597	2/587 (0.3%)	3/587 (0.5%)	7/587 (1.2%)	12 ^b

TVU = transvaginal ultrasound; n = number of women

^a support catheter that has since been discontinued was used in 5 of 7 cases

^b does not include 23 subjects lost to follow up and did not undergo confirmation testing. May include patency seen on HSG.

^c includes proximal placement, distal placement and not classified.

Phase II study (STOP 10)

Expulsion: The 3-month HSG indicated a unilateral insert expulsion in one woman. In this case, at the time of placement, the Investigator suspected the insert was not placed far enough into the tube and had anticipated that the insert would be expelled.

Perforations: There were 7 women who were diagnosed with a perforation. The 3-month HSG indicated perforation in 4 women and a repeat HSG at 18 months diagnosed a perforation in a fifth woman. A sixth woman was found to have a peritubal perforation noted on gross examination following device removal for pain. During laparoscopy on a seventh woman, it was discovered that the right device had perforated the tube, although HSG demonstrated tubal occlusion.

Five of these women went on to laparoscopic sterilization and the sixth woman had a cornual resection for persistent pain. Histological evaluation in the sixth woman confirmed the device had perforated the tube just past the uterotubal junction. The seventh received a Filshie clip sterilization of the left tube after failure to place the Essure insert and began to rely on the clip and the insert in the occluded, perforated right tube for contraception.

Intraperitoneal devices: Intraperitoneal devices were successfully retrieved in 4/6 women. There were no adverse tissue reactions noted at the insert retrieval site in the peritoneal cavity in any of these women. One woman had a device left in the peritoneal cavity. She reported some pain on follow-up visits that did not seem related to the perforated device, as the reports were of rare twinges on the side with both a well-placed Essure insert and a Filshie clip. The other woman whose perforated device was not retrieved had an insert placed into the

myometrium (embedment/partial perforation). The device has remained in place with no sequelae.

Of these 7 perforations, 5 occurred with the support catheter that was discontinued due to its high rate of perforations.

Pivotal study (STOP 2000)

In the pivotal study, x-ray noted a perforation in one woman that was diagnosed on the day of the placement. Therefore this subject was followed for safety only. This woman did not undergo HSG and had a tubal ligation after Day 0. Another woman found an expelled insert ex-vivo prior to the 3-month visit.

Table 7–7 summaries AEs that initially impacted the woman’s ability to rely on the Essure insert for contraception.

Table 7–7: Adverse events preventing reliance in women with bilateral placement

Event	Number initially diagnosed	Suspected cause	Management	Ultimately affected reliance
Total	21 (4.5%)			12 (2.6%)
Expulsion	14 (3.0%)	Proximal placement (13); placement into endometrial tissue (1).	9 successfully underwent second insert placement procedure; 4 underwent laparoscopic sterilization; One woman with Crohns disease required laparotomy for management	5 (1.1%)
Perforation	4 (0.9%)	Pre-existing tubal occlusion (2); poorly identified ostium (2);	3 underwent laparoscopic sterilization; one requesting repeat insert placement procedure	4 (0.9%)
Unsatisfactory device location and perforation	1 (0.2%)	Unsatisfactory placement on one side (0 mm trailing) and difficulty placing on the perforated side.	Proximal insert removed hysteroscopically during lap sterilization, perforated side was suspicious on x-ray and confirmed during lap sterilization	1 (0.2%)
Unsatisfactory placement	2 (0.4%)	Unsatisfactory placement (1); hydro-salpinx (1)	Awaiting women’s decision to undergo lap sterilization (2)	2 (0.4%)

ESS305 TVU study (Study 16974)

Unsatisfactory device location was reported in 45/519 (8.67%) women’s TVU examination. Forty of these women underwent HSG of which 36 were told to rely (Table 7–8).

Table 7–8: Transvaginal ultrasound findings

Insert location	N	Location (if not optimal)	Right (N)	Left (N)	Both (N)
Transvaginal ultrasound was performed on Total = 528 subjects^a					
Inserts were identified on both sides with Transvaginal ultrasound = 519 subjects					
Optimal or satisfactory	474				
Unsatisfactory	45	Distal placement	12	12	2
		Proximal placement	0	2	0
		Perforation	2	1	0
		Expulsion	0	1	0
		Unclassified	4	4	2
		Other	6	7	0

^a 527 subjects with bilateral placement and 1 subject with unilateral placement

7.7.2 Literature

In cases reporting perforation, some women were asymptomatic. When symptoms were reported (most commonly abdominal pain), the onset ranged from days to months following the placement procedure. (72, 73, 74, 75, 76, 77, 78, 79, 80, 81) A 7-year retrospective study evaluated complications of tubal sterilization with Essure in 4,306 women. (37) Within this study, one woman (0.02%) had a tubal perforation discovered after pregnancy. Three women (0.06%) had the device inadvertently placed in the myometrium, and 2 women (0.04%) had asymptomatic migration of the device into the abdominal cavity. In both of these cases, a second device was placed and the migrated device left in situ. The women remained asymptomatic. A small number of case reports discuss small bowel obstruction and perforation with Essure. (82, 83, 84, 85, 86)

7.7.3 Postmarketing monitoring

A query of the global safety database on Essure through 15 JUNE 2015 revealed 4,308 reports of unsatisfactory device location. The majority of these cases (83.9%) were medically confirmed. It is acknowledged that some Essure devices when placed under hysteroscopic visualization will not be appropriately located to achieve tubal occlusion. This is the basis for conducting the Essure Confirmation Test. The reporting frequency for all unsatisfactory device locations is 0.4%.

Of the 4,308 reports, 777 were expulsions, 346 reported a missing insert during imaging or surgery and 3237 reported either a perforation/migration/embedment/dislocation or more than 1 event. The majority of events (95.6%) were from the US where the highest number of placements have occurred. Among the medically confirmed cases, most unsatisfactory device locations were diagnosed at the confirmation test (48.3%) with 17.2% detected at the placement procedure, 4.8% prior to the scheduled confirmation test, 8.2% beyond the

scheduled confirmation test at 3 months and 21.6% in which the time of diagnosis was not provided.

Devices primarily located in the uterus were removed hysteroscopically when indicated, while those located distal to the optimal location were most commonly managed with laparoscopy alone or laparoscopy in combination with additional procedures.

7.7.4 Conclusion

Unsatisfactory device placement, including perforations, are known complications of Essure and their recognition during the confirmation test with appropriate management continues to be a focus of the ongoing training and patient management programs for Essure. In the clinical trials, the incidence of all unsatisfactory device locations ranged from approximately 3.5% to 4.5%. The postmarketing data is consistent with that in the clinical trials, all of which is reflected in Sections IX, X, and XVII of the IFU (Section 10.6, Appendix 6).

7.8 Pain (persistent/chronic)

Chronic pelvic pain is a common gynecological problem with an estimated prevalence of 38 per 1000 in women aged 15 to 73 years, a rate comparable with that of asthma (37/1000) and chronic back pain (41/1000).⁽⁸⁷⁾ There is wide variation in clinical evaluation of women with chronic pelvic pain. Diagnostic laparoscopy is often carried out after referral to a gynecologist as an initial investigation to uncover pathological causes, for example, endometriosis or adhesions, but has negative results in over half the cases.

Chronic pelvic pain in association with sterilization has been described in the literature.^(88, 89) Hiemstra et al ⁽⁸⁹⁾ recommend a bio-psychosocial approach to the issue of post-sterilization pain, though this is not always possible.

7.8.1 Sponsor's clinical data

Phase II study (STOP 10)

In the phase II study, 205 women completed patient questionnaires about their procedure. The remaining 22 women who underwent procedures did not complete the questionnaire because they were bilateral failures (20 women) or had unilateral placement only (2 women). Post-procedure pain was reported in 156/205 (76%) women. The pain resolved within one day in 58% women, within 3 days in 87% and within one week in 99% of these women. Pain was resolved within 2 weeks of the procedure for the remaining 2 women (information from the 3-month follow-up visit). Thirteen of 199 (6.5%) women who completed the 6-month follow-up visit reported recurrent pain (reported on two or more visits that may or may not be consecutive) while persistent pain (reported at every visit) was not reported by any woman.

Table 7–9 summarizes the pain noted at each of the follow-up visits. Types of pelvic pain reported most commonly included dysmenorrhea, dyspareunia, and other pelvic pain (mostly ovulation pain). Of the women reporting dysmenorrhea in the first 3 months following the procedure, 11/29 (38%) reported pain in the first week following the procedure and 23/29 (79%) reported pain in the first month following the procedure. Women also reported symptoms of pain that were not related to the insert.

Table 7–9: Pain reported at follow-up visits in phase II study

Follow-up visit	Dysmenorrhea	Pelvic pain, Dyspareunia	Other pelvic pain	Other pain
3-month	29/203 (14%)	17/203 (8%)	5/203 (2%)	2/203 (<1%)
6-month	11/199 (6%)	3/199 (2%)	3/199 (2%)	1/199 (<1%)
12-month	5/196 (3%)	0	5/196 (3%)	0
18-month	2/193 (1%)	0	10/193 (5%)	2/193 (1%)
24-month ^a	8/194 (1%)	0	6/194 (3%)	6/194 (3%)
36-month ^a	7/182 (4%)	1/182 (<1%)	2/182 (1%)	1/182 (<1%)
48-month ^a	9/176 (5%)	2/176 (1%)	2/176 (1%)	2/176 (1%)
60-month ^a	4/171 (2%)	1/171 (<1%)	3/171 (2%)	2/171 (1%)

^a No data reported for some women that indicated other types of pain (N = 2 at 24 months; N = 1 at 36, 48 and 60 months).

Pivotal study (STOP 2000)

Table 7–10 summarizes various types of pelvic pain reported by women during follow-up visits.

Of the women who reported pain, one woman reported “persistent” pelvic pain at two years of follow-up. None of the women reporting at 3, 4 and 5 year follow-up indicated persistent pelvic pain of any kind. Twenty-nine women (6.1%) had at least two (recurrent) episodes of dysmenorrhea, 18 women (3.8%) reported “recurrent” dyspareunia, 14 women (3.0%) reported “recurrent” ovulatory pain, and 25 (5.3%) reported “recurrent” other pelvic pain that could not clearly be classified into a different category.

Table 7–10: Pelvic pain reported at follow-up visits (pivotal study)

Follow-up visit	Pelvic pain			
	Dysmenorrhea	Dyspareunia	Ovulatory	Other ^c
Baseline (N=518)	183 (35%)	22 (4.2%)	N/a	N/a
3-month post-device placement (N=467)	29 (6.2%)	29 (6.2%)	5 (1.1%)	32 (6.9%)
Post-alternative contraception				
3-month (N=440)	20 (4.5%)	10 (2.3%)	6 (1.4%)	26 (5.9%)
6-month (N=436)	15 (3.4%)	8 (1.8%)	3 (0.7%)	16 (3.7%)
12-month (N=460)	17 (3.7%)	15 (3.3%)	5 (1.1%)	27 (5.9%)
18-month (N=410)	14 (3.4%)	9 (2.2%)	10 (2.4%)	11 (2.7%)
24-month (N=435)	22 (5.1%)	9 (2.1%)	22 (5.1%)	13 (3.0%)
36-month (N=422) (data missing on 1 patient)	14 (3.3%)	7 (1.7%)	12 (2.8%)	6 (1.4%)
48-month (N=402)	3 (0.7%)	5 (1.2%)	6 (1.5%)	11 (2.7%)
60-month (N=386)	14 (3.6%)	8 (2.1%)	10 (2.6%)	9 (2.3%)
Recurrent (N=473) ^a	29 (6.1%)	18 (3.8%)	14 (3.0%)	25 (5.3%)
Persistent (year 1), (N=460) ^b	0	0	0	1 (0.2%)
Persistent (year 2), (N=435) ^b	0	0	0	1 (0.2%)
Persistent (year 3), (N=422) ^b	0	0	0	0
Persistent (year 4), (N=402) ^b	0	0	0	0
Persistent (year 5), (N=386) ^b	0	0	0	0

^a *Recurrent*: symptom reported at more than one visit.

^b *Persistent*: symptom reported at all prior visits.

^c *Other*: defined as pelvic pain that was not reported to be dysmenorrhea, dyspareunia or ovulatory pain.

ESS305 TVU study (Study 16974)

Table 7–11 presents AEs with abdominal, back or pelvic pain.

Table 7–11: Adverse events of abdominal or pelvic pain and related terms

MedDRA System Organ Class	MedDRA Preferred Term	Total = 597 (100%) N (%) events
Gastrointestinal disorders	Abdominal pain	16 (2.7)
	Abdominal pain, lower	1 (0.2)
Musculoskeletal and connective tissue disorders	Back pain	3 (0.5)
Reproductive system and breast disorders	Adnexa uteri pain	2 (0.3)
	Dysmenorrhea	15 (2.5)
	Dyspareunia	4 (0.7)
	Metrorrhagia	5 (0.8)
	Ovulation pain	4 (0.7)
	Pelvic pain	28 (4.7)

MedDRA = Medical Dictionary for Regulatory Activities; N = number

7.8.2 Literature

Pain associated with the Essure procedure may occur at three general time points: pain during the actual placement and immediately post-procedure, pain that occurs acutely following the procedure and pain that is persistent or “chronic pain.” Exact definitions for these time periods vary by study, however, literature exists for each of these broad categories.

Studies examining pain at time of the procedure (90, 45, 48) show that while the majority of women rate pain associated with the procedure as tolerable, a small number of women (3.1 to 17%, depending on study) rated the procedure as painful.

There are two studies comparing hysteroscopic versus laparoscopic sterilization with regard to post-procedural pain. Duffy et al (49) performed a cohort controlled comparative study in the UK. Self-assessment diaries were completed on days 7 and 90 post-procedure. Post procedure, 31% of the Essure group reported moderate or severe pain compared with 63% of the laparoscopic group (P=0.08). Syed et al (91) examined 20 cases each of laparoscopic and hysteroscopic sterilization (Staten Island, NY). Women were asked about pain one week and four weeks post-procedure on a scale of 0 (“no pain”) to 10 (“worst pain ever experienced”). Immediately post-procedure, mean pain score in the hysteroscopic group was 0.6 (median=0) and in the laparoscopic group it was 5.05 (median=6) ($t = -8.17, p < .0001$). One week post-procedure, none of the women with Essure reported pain, while the mean score in the laparoscopic group was 2.65 (median=3). Four weeks post procedure, 35% of the laparoscopic group continued to experience pain (mean=0.60, median=0).

With regard to acute pain, in Sinha et al, (45) women were surveyed three months post-Essure placement. At this time, 5/76 (6%) women reported new pain or discomfort since the procedure and 2/76 (3%) described new pain or discomfort with sexual intercourse. No further information is available on these women. Yunker et al (92) performed a retrospective cohort, single center study (Vanderbilt University, TN) of 458 women who underwent Essure placement between JAN 2005 and JUNE 2012. A pre-procedural diagnosis of chronic pain (chronic pelvic pain, chronic low back pain, chronic headache, fibromyalgia) was noted in 8.7% of women. They found the incidence of acute post-Essure pelvic pain (within 3 months) was 8.1%. Women with previous diagnoses of any chronic pain were significantly more likely to complain of acute pelvic pain (odds ratio 6.81; 95% confidential interval, 2.98 to 15.73) following the placement of Essure.

A similar study found that pre-existing pain may also affect laparoscopic procedures. Rudin et al (93) assessed potential predictive factors of post-operative pain in 59 women undergoing BTL with bipolar cautery. They found that in their multiple regression model, report of preoperative pain (eg, back, neck, shoulder, hip, leg, arm, hand, abdomen, headache) and preoperative heat perception were predictive variables for post-operative pain.

A recent analysis of the MarketScan commercial claims database revealed that women undergoing both hysteroscopic and laparoscopic sterilization frequently reported pain in the 6 months prior to their procedure (interim data on file). Claims data from JAN 2010 through DEC 2012 were analyzed. Six months of pre-index data were required for inclusion. Claims for post-partum tubal ligation or claims associated with a pregnancy in the preceding 6 weeks were excluded. Women undergoing any additional surgical procedure on the day of the sterilization procedure were also excluded. This resulted in 12,031 cases of hysteroscopic sterilization (which included both Essure and Adiana, a hysteroscopic device now off the market) and 7,286 cases of laparoscopic sterilization. Pre-procedure diagnoses are summarized in [Table 7–12](#).

Table 7–12: Women with diagnoses of pain, 6 months preceding sterilization procedure

	Hysteroscopic procedure (N = 12,031)	Laparoscopic procedure (N = 7,286)
Mean age	37 years	35.8 years
Any pain	2827 (23.5%)	1971 (27.1%)
Pelvic pain/lower abdominal pain	1515 (12.6%)	1129 (15.5%)
Low back pain	1097 (9.1%)	3724 (9.9%)
Chronic headache	529 (4.4%)	397 (5.4%)
Fibromyalgia	248 (2.1%)	150 (2.1%)

Chronic Pain: Yunker et al (92) continued to follow women in their study beyond the 3 month time period and found that at 6 months, the incidence of pain decreased from 8.1% (at 3 months) to 4.2%. As with acute pain at 3 months, women with a previous diagnosis of any chronic pain (pelvic, low back, headache, and fibromyalgia) were significantly more likely to report chronic pelvic pain at 6 months, odds ratio, 6.15 (Confidence Interval: 2.09-18.05).

Reports of pain post-procedure appear to be lower in studies outside of the US. In a large case series of 4,274 women who underwent Essure at a Spanish teaching hospital,(94) a total of seven women (0.16%) presented with chronic pain requiring insert removal. Mean time between placement and removal was 29.4 months. In all cases, the symptoms disappeared after Essure removal.

Comparative studies help place rates of pain in context. Conover et al (95) performed a retrospective review of the MarketScan claims database and compared incidence of opioid-managed pelvic pain within 12 months after hysteroscopic versus laparoscopic sterilization. They found the rates of opioid-managed pelvic pain were low overall (0.88% in the hysteroscopic group and 0.93% in the laparoscopic group) and that there was no statistical difference between the two groups.

Several studies have attempted to identify factors related to the development of pain following placement of Essure. Case reports show that perforation of the tube or uterus may

cause pain in some women (96) and removal may provide relief. In a systematic review of PubMed, Scopus, Embase, and the Cochrane library, Adelman et al (97) identified seven cases of persistent pain without evidence of perforation. All women underwent bilateral insert removal via hysteroscopy or laparoscopy. Four of the 7 women experienced complete resolution of pain, one noted some improvement in daily pain but still experienced intermittent abdominal pain, and outcome on one was unknown. In the seventh woman, chronic pelvic pain continued despite laparoscopic salpingectomy and appendectomy.

7.8.3 Postmarketing monitoring

A search of the global safety database revealed 3,539 cases containing at least one event of abdominal, pelvic and back pain (data lock point: 08 JULY 2015). Based on the total estimated exposure of 1 million devices sold, the reporting rate for abdominal, pelvic and back pain was 0.3%. A search for cases specifically reporting “chronic”, “persistent”, “consistent”, “constant” or “continuous” pain indicated that these encompassed about 10% of all cases (0.03%).

All reports were considered regardless of causality, severity, duration, temporal relationship to the insert placement process and any co-reported conditions. Each case may contain more than one reported term. Of all events, 58 were reported from interventional studies with the different models of Essure, 2,804 were spontaneous reports, 637 were from observational studies, 40 from literature and 58 from interventional studies.

About half of all cases of abdominal, pelvic and back pain were medically confirmed reports, and half were from consumers. Pain was often reported in the context of other AEs. Heavy menstrual bleeding was a common complaint in addition to pain. Medically confirmed cases more frequently reported additional events of unsatisfactory device location, including perforation, expulsion, improper placement or device malfunction. Consumer reports at times included up to 20 or 30 events describing changes ranging from menstrual cycle disturbances, to cognition, skin appearance, sexuality, mood and general well being.

Onset latency information indicated that the onset of pain may range from as early as immediately after placement to several years after placement. Most reports, report an onset within the first year after insert placement. The onset of events five years or later is very rare.

7.8.4 Conclusion

Abdominal and pelvic pain are common in women of reproductive age. Pain after Essure sterilization is an anticipated event reflected in the IFU (Section 10.6, Appendix 6) and PIB (Section 10.7, Appendix 7). Clinical trial data indicate that the highest frequency of abdominal, genital and back pain can be expected at the time of the insert placement and shortly thereafter, but women may also experience pain immediately after the insert placement period. Chronic pelvic pain in association with female sterilization has been described in the literature. The etiology often remains unclear, although improper placement

of involved devices and a history of a previous pain diagnosis have been identified as potential contributory factors.

7.9 Menstrual changes

Abnormal uterine bleeding is a broad spectrum term that includes irregular, heavy or prolonged menstrual bleeding or an altered bleeding pattern. The condition has a high prevalence in all countries; it is estimated that the worldwide prevalence of abnormal uterine bleeding is between 4% and 52% depending on the study techniques used.(98, 99)

In 1951, Williams et al (100) first hypothesized that sterilization may affect a woman's risk of abnormal bleeding (defined as increased intermenstrual bleeding and increased flow). The phrase "post tubal sterilization syndrome" was coined to describe a variety of symptoms occurring after female sterilization.(101) This may include: abnormal bleeding and/or pain, changes in sexual behavior and emotional health, exacerbation of premenstrual symptoms, significant menstrual disturbances leading to further surgery (eg, hysterectomy, tubal reanastomosis). The existence and potential mechanism of a post-tubal ligation syndrome has been debated.(101, 102) According to the ACOG Practice Bulletin (2) on the Benefits and Risks of Sterilization: "prospective studies that account for confounding factors, such as presterilization use of hormonal contraception, have found that tubal occlusion has little or no effect on menstrual patterns". The studies cited by ACOG regarding menstrual patterns/changes during various procedures of sterilization are provided in Section 10.5 (Appendix 5).

7.9.1 Sponsor's clinical data

Phase II study (STOP 10)

Post-procedure bleeding of any kind was reported in 83% of respondents. Bleeding resolved within 3 days for most women and within one week for 96% of women. When asked to compare this bleeding with bleeding during normal menses, 71% women reported it was less than their previous period, 17% reported it was the same, 8% a little more and 4% as a lot more than their previous period.

Table 7–13 provides a summary of bleeding reports at the follow-up visits.

Persistent bleeding was not reported in this study. Fifteen women (15/199 [7.5%]) reported recurrent menstrual changes, one woman (1/199 [0.5%]) reported recurrent post-coital bleeding and no women reported persistent menstrual changes.

Table 7–13: Bleeding reported at follow-up visits in phase II study

Follow-up visit	Irregular menses	Spotting	Changes in flow	Other and unspecified
3-months	not asked	not asked	not asked	not asked
6-months	3/199 (2%)	6/199 (3%)	3/199 (2%)	0
12-months	6/196 (3%)	5/196 (3%)	4/196 (2%)	0
18-months	9/193 (5%)	4/193 (2%)	5/193 (3%)	0
24-months ^a	4/194 (2%)	4/194 (2%)	10/194 (5%)	3/194 (1.5%)
36-months ^a	4/182 (2%)	3/182 (2%)	4 (182 (2%)	3/182 (1.6%)
48-months ^a	4/176 (2%)	3/176 (2%)	9/176 (5%)	3/176 (1.7%)
60-months ^a	16/171 (9%)	3/171 (2%)	11/171 (6%)	3/171 (1.8%)

^aNo data reported for some women that indicated other/unspecified bleeding (N=2 at 24 months; N=1 at 36, 48 and 60 months)

Pivotal study (STOP 2000)

Table 7–14 presents changes in menstrual function reported at follow-up visits.

Table 7–14: Changes in menstrual function reported at follow-up visits (pivotal study)

Follow-up visit	Irregular menses	Bleeding between menses	Heavier than usual menstrual flow	Less than usual menstrual flow
Baseline, N=518	9 (1.7%)	12 (2.3%)	N/A	N/A
3-months post-device placement	48 (10.3%) N=467	110 (23.6%) N=466	89 (19.2%) N=463	56 (12.1%) N=463
Post-alternative contraception				
3-months	36 (8.2%) N=440	40 (9.1%) N=440	96 (21.9%) N=439	55 (12.5%) N=439
6-months	36 (8.2%) N=437	29 (6.6%) N=437	94 (21.6%) N=435	57 (13.1%) N=435
12-months	35 (7.7%) N=455	31 (6.7%) N=460	77 (16.8%) N=458	67 (14.6%) N=458
18-months	19 (4.6%) N=410	42 (10.2%) N=411	70 (17.0%) N=411	63 (15.3%) N=411
24-months	20 (4.6%) N=435	32 (7.4%) N=435	89 (20.6%) N=432	53 (12.3%) N=432
36-months	31 (7.4%) N=420	25 (6.0%) N=420	83 (20.2%) N=411	47 (11.4%) N=411
48-months	33 (8.4%) N=393	33 (8.3%) N=396	69 (17.9%) N=386	52 (13.5%) N=386
60-months	45 (11.7%) N=386	29 (7.5%) N=386	74 (19.6%) N=377	40 (10.6%) N=377
Recurrent ^a	70 (14.8%) N=473	89 (18.8%) N=473	177 (37.5%) N=472	110 (23.3%) N=472
Persistent (year 1) ^b	3 (0.7%) N=455	2 (0.4%) N=460	7 (1.5%) N=458	12 (2.6%) N=458
Persistent (year 2) ^b	0 (0.0%) N=435	1 (0.2%) N=435	4 (0.9%) N=432	3 (0.7%) N=432
Persistent (year 3) ^b	0 (0.0%) N=420	1 (0.2%) N=420	4 (1.0%) N=411	2 (0.5%) N=411
Persistent (year 4) ^b	0 (0.0%) N=393	0 (0.0%) N=396	3 (0.8%) N=386	1 (0.3%) N=386
Persistent (year 5) ^b	0 (0.0%) N=380	0 (0.0%) N=386	2 (0.5%) N=377	0 (0.0%) N=377

^a *Recurrent*: symptom reported at more than one visit during the follow-up period (ie, symptoms do not have to be reported on consecutive visits). The denominator ("N") is the sum of all unique women who responded over the course of their follow-up period. Not all women responded at all follow-up visits.

^b *Persistent*: symptom reported at all visits during the follow-up period. The denominator ("N") is the total number of unique women responding at the latest follow-up visit.

ESS305 TVU study (Study 16974)

In this ongoing study, no data is collected asking women to compare the bleeding that they experienced during the study to their typical menstrual periods and no specific data is collected on the recurrence or persistence of such events.

A listing of all events related to uterine or vaginal bleeding is in [Table 7–15](#). Please note: women may have reported more than one such event.

Table 7–15: Adverse events related to uterine or vaginal bleeding

Adverse events	n (%)
Amenorrhea	3 (0.5%)
Dysfunctional uterine bleeding	2 (0.3%)
Hypomenorrhea	1 (0.2%)
Menorrhagia	23 (3.9%)
Menstruation delayed	1 (0.2%)
Menstruation irregular	3 (0.5%)
Polymenorrhea	1 (0.2%)
Uterine hemorrhage	9 (1.5%)
Vaginal hemorrhage	14 (2.3%)

7.9.2 Literature

Menstrual changes post-Essure have been reported in several small studies. A summary of these studies is shown in Section 10.5 (Appendix 5). The only study to compare bleeding patterns between Essure and laparoscopic sterilization looked only at short-term bleeding.⁽⁹¹⁾ The study group was a convenience sample of 40 women: 20 underwent Essure placement and 20 underwent laparoscopic tubal ligation. Bleeding was assessed during the week following the procedure. In the Essure group, 40% reported minimal spotting and 60% reported no bleeding. In the laparoscopic tubal ligation group, 25% reported minimal spotting and 75% reported light bleeding. There are no long term comparisons between the two procedures in the literature.

7.9.3 Postmarketing monitoring

A search of the global safety database through 15 JUNE 2015 revealed 2,232 cases containing at least one event of menstrual bleeding disturbances resulting in a reporting rate of 0.22%, with 0.12% from consumers/women and 0.10% were classified confirmed by Health Care Provider (HCP).

Clustering into subtypes of menstrual bleeding disturbances showed that increased or irregular menstrual bleeding (or unspecific changes in bleeding pattern) was reported at a

frequency of 0.2%, painful menstrual bleeding at 0.03%, and decreased or infrequent menstrual bleeding at 0.02%.

The pathophysiology of heavy or prolonged menstrual bleeding in Essure users outside of unsatisfactory device location is not clear. Review of the literature shows that changes in bleeding pattern were examined also in conjunction with laparoscopic tubal sterilization. Studies that considered confounding factors such as pre-sterilization bleeding pattern during hormonal contraceptive use found no link between the sterilization procedure and changes in bleeding. It can be assumed that similar confounding factors may also apply for menstrual bleeding disturbances in conjunction with Essure.

7.9.4 Conclusion

Abnormal, heavy or prolonged menstrual bleeding is a condition with a high prevalence in the fertile female age group. Other than the possibility of endometrial irritation in very proximally placed inserts, the pathophysiology behind heavy or prolonged menstrual bleeding in Essure users is not clear. Review of the literature shows that when changes in bleeding pattern were examined in conjunction with laparoscopic tubal sterilization, studies which considered confounding factors, such as pre-sterilization, bleeding pattern, or hormonal contraceptive use did not substantiate any causal relationship between sterilization and altered bleeding patterns. These confounding factors may also apply for menstrual bleeding changes with Essure as well.

7.10 Headache

Headache is a very common symptom that is reported in almost all individuals at some time in their life. The lifetime prevalence of any form of headache, migraine, and tension headache in women were 99%, 25%, and 88%, respectively. The prevalence of migraine and tension-type headache over 1-year in women were 15% and 86% respectively.(103) Chronic daily headache is reported in 5% of females.(104)

7.10.1 Sponsor's clinical data

Phase II study (STOP 10): In the phase II study 8/206 (3.9%) of women reported headache (not including migraine). An additional 2 women (1.0%) reported events of migraine.

Pivotal study (STOP 2000): Headache and related events during the 5 year follow up period have been analyzed based on data provided in the final 5-year pivotal study report. Out of 518 women with hysteroscopy performed, a total of 98 (18.9%) women reported 226 such events. Of these events, headache only was reported in 76 (14.7%) women (167 events); headache with other symptoms in 25 (4.8%) women (32 events), migraine only in 14 (2.7%) women (24 events) and migraine with other symptoms in 3 (0.6%) women (3 events). The range of events per woman was 1 – 17 (one woman reported 17 events of headache only). A minority of these

events (12 events total in 9 women, 0.7%) were considered to be at least possibly related to Essure.

ESS305 TVU study (Study 16974): In the Essure TVU study, headache was reported in 0.7% women and migraine in 0.2%.

7.10.2 Literature

Duffy et al 2005 (49) reported one woman with headache at the post-operative visit following Essure placement. No other publications were found on Essure and headache.

7.10.3 Postmarketing monitoring

A search of the global safety database through 15 JUNE 2015 revealed 773 reports describing headache in association with Essure yielding a reporting rate of 0.08%, of which 0.07% of reports were from consumers/women, while 0.01% were classified as medically confirmed.

The overwhelming majority of cases describe either “headache” or “migraine”. Over 90% of headache/migraine-related events were classified as non-serious. The reporting rate for cases containing terms suggestive of chronic or recurrent headaches was 0.008%.

Most cases in the headache/migraine case cluster are non-medically confirmed consumer cases. Many of these cases are reported in conjunction with a large range of different symptoms or events. The reporting of headache in conjunction with many other symptoms or events increases the probability of the headache being secondary to these other conditions.

In only a very small number of headache/migraine cases a targeted diagnostic intervention (head/neck imaging or neurological examination) is reported. In a few of these cases, the diagnostic intervention revealed a cause for the headache/migraine. These causes are all assessed to be unrelated to Essure (eg, brain infarction, pituitary tumor, lymphoma).

7.10.4 Conclusion

Headache and migraine are conditions with a very high prevalence, especially in the female population between the ages of 25 and 50 years of age. Thus headache is a commonly reported symptom. External triggers, as well as hormonal changes have been described as potential causal factors in these conditions. The analysis of the available postmarketing case reports of headache or migraine reported in conjunction with Essure shows a small number of serious cases and very few diagnostic investigations, suggesting that the headache symptoms in the majority of cases were not determined to be of high clinical significance. When a medical evaluation of the headache was performed, the resulting diagnoses were unrelated to Essure.

A change of contraceptive method to Essure can often entail stopping previous hormonal contraceptives which could be postulated to exert a certain influence on the occurrence rates

of headache or migraine. Apart from this change, no plausible pathophysiological link is currently known between Essure’s localized mode of action and the development of headache or migraine.

7.11 Infectious complications

Any time a body cavity is entered, there is a risk of infection. This is no different with hysteroscopy. Infectious morbidity following hysteroscopy, however, is extremely low.(105) As with all hysteroscopic procedures, Essure should not be placed in women with any evidence of active or recent upper/lower genital tract infection as stated in the IFU (Section 10.6, Appendix 6).

7.11.1 Sponsor’s clinical data

The incidence of infectious complications in the phase II, pivotal and the TVU study is summarized in Table 7–16.

Table 7–16: Infectious complications in Essure studies

Study	Number of women undergoing placement attempt	Number of infectious complications (upper genital tract) ^a	Incidence of infectious complications
Phase II (STOP 10)	227	0	0%
Pivotal (STOP 2000)	507	5 2 cases of PID, 2 of endometritis, 1 pelvic infection	0.99%
ESS305 TVU (Study 16974)	597	4 ^b 2 cases of PID 2 cases of endometritis	0.34%

TVU = transvaginal ultrasound; PID = pelvic inflammatory disease

^a may have been reported remote from the time of insert placement and therefore may have been unrelated to placement or devices.

^b Does not include one case of sepsis after surgery for uterine prolapse 2+ years after Essure placement.

7.11.2 Literature review

Adelman et al (97) performed a systemic review of complications encountered with Essure. They reviewed 594 sources, including clinical trials, previous systemic reviews, case series, and case reports. They found a single case report of a ruptured tubo-ovarian abscess 3 years after placement of Essure. Group A β-hemolytic streptococci were isolated and the woman ultimately had a total abdominal hysterectomy. Adelman et al (97) state: “While it was not clear that the Essure inserts contributed to the ascending infection, it did not prevent it as is

observed with other forms of tubal sterilization.” A second case of bilateral cornual abscesses followed endometrial ablation after Essure placement one year previously. Povedano et al (37) did a retrospective review of medical records of the 4,306 women who underwent Essure placement at their own site and found two cases of PID, both of which resolved with IV antibiotics.

7.11.3 Postmarketing monitoring

A total of 240 cases of potential pelvic infection were found, of which 214 were spontaneous reports. All but 18 cases were classified serious, and less than half of the cases (n=100) were medically confirmed. Age ranged from 20 to 47 years, with a median of 33 years in cases in which it was provided. One case was fatal secondary to toxic shock syndrome. The reported rate of infectious complications is 0.02%.

There are a total of 16 cases where documented serious infections related to Essure have been reported. Eleven of these cases had onset of symptoms within one month of the procedure. Two of these involved concomitant endometrial ablation with Essure placement one of which involved a bowel injury secondary to thermal necrosis and one Group B streptococcal infection. Of the 9 cases with Essure placement alone, two required a hysterectomy due to necrotic uterine tissue. One case involved unrecognized ongoing PID with infected tissue noted during the procedure and a second case diagnosed as severe PID had necrotic endometrium and myometrium noted on pathological examination. There was one case of tubal abscess following a difficult placement associated with perforation which required a hysterectomy. Two cases of tubal abscess were successfully managed with drainage and antibiotics. One case of PID was managed with IV antibiotics without surgery. One case of toxic shock syndrome secondary to Group A strep had a fatal outcome. The final case was a pyosalpinx in a woman with known salpingitis who had Essure placed in preparation for in-vitro fertilization. This case was managed with radiological drainage and antibiotics.

Five cases were reported one month or more following the Essure procedure. One report of PID was treated with IV antibiotics. There was one tubal abscess treated with laparoscopic drainage and IV antibiotics 4 months post placement. Two cornual abscesses diagnosed one and a half years post placement were managed by hysterectomy and a case of tuboovarian abscess was managed by salpingo-oophorectomy. Following a hysterectomy performed for pain, menorrhagia and dysmenorrhea, the surgeon dissected the fallopian tubes and diagnosed an infection.

In summary, when infections associated with Essure placement do occur, they are relatively infrequent, and generally manageable with favorable outcomes. Placement of Essure in a woman with an active genital tract infection can exacerbate the infection which may result in surgical removal of the uterus and fallopian tubes. Antibiotic treatment alone may be sufficient in women diagnosed with PID following Essure placement, although the presence of an abscess should be investigated and if found may require drainage. The number of

reports of PID following Essure is low and may indicate that when bilateral tubal occlusion is achieved, this creates a barrier to ascending genital tract infections. The cases reported with endometrial ablation performed sequentially with Essure placement indicate a higher risk of serious infectious complications. Section VI (“Warnings”) of the IFU (Section 10.6, Appendix 6) in the US has an extensive discussion regarding the risk of performing these procedures concomitantly in the warnings section.

7.11.4 Conclusion

Data indicate that the rate of serious infection following placement of Essure is low. It is important that physicians follow guidance outlined in the IFU to help prevent such occurrences.

7.12 Allergy reactions/hypersensitivity (to nickel)

Nickel allergy was a contraindication for Essure when the product was first marketed. Based on biocompatibility data (see Section 3) and lack of correlation between history of nickel allergy and clinical outcome,(23) in agreement with the FDA, this was changed to a warning in 2009. The IFU continues to state: “women who are allergic to nickel may have an allergic reaction to this device, especially those with a history of metal allergies.” Typical allergy symptoms include rash, pruritis, and hives.

7.12.1 Sponsor’s clinical data

One subject in the ESS305 TVU study (Study 16974) had an increased sensitivity to metals, reported as an AE by MedDRA PT “allergy to metal”. One case of allergic reaction to distention medium was reported in the Pivotal Study. No further information is available on this case. All other allergic reactions and similar events reported were related to medication or environmental exposure of some type.

7.12.2 Literature

Systemic allergic reactions to nitinol in Essure are also rare. To date, there have been three published case reports describing reactions necessitating surgical removal of the devices.(106)

Peter et al (107) conducted a prospective study to determine if the nitinol in Essure inserts can cause nickel sensitivity. The study included 200 women. All women received two skin patches consisting of a nickel and a control. All women underwent Essure sterilization regardless of the skin test results. Results showed that 25% of the women had a positive skin test prior to placement of Essure. No clinical symptoms of allergy were reported in any woman, regardless of skin patch results.

Such lack of concordance has been demonstrated with other nitinol-containing devices. For example, nitinol is commonly used in heart and endovascular prostheses. Despite the nickel

content, there are few reports of nickel allergy associated with these devices, with a risk estimate of 1/17,000 devices.(108)

Zurawin et al (23) reported on a case series regarding AEs of suspected nickel hypersensitivity associated with Essure, the subsequent clinical outcomes in women with these insert placements, and the correlation of nickel allergy-related AEs with positive results of nickel patch testing. Through multiple sources, 63 reports of suspected nickel hypersensitivity were identified. Twenty women had patch testing: 13 had positive and 7 had negative test results. In 9 of the 13 women who tested positive on the patch test, Essure was removed. Symptoms resolved in 4 of the 9 women, 2 did not resolve, and 2 were unknown.

7.12.3 Postmarketing monitoring

A broad query of the global safety database revealed 641 reports as per 15 JUNE 2015 where a suspicion of allergic reactions to Essure was raised by the reporting healthcare professional or consumer/women (reporting frequency, 0.06% of all placements). Most reports were spontaneous reports from the US. Ninety one women were identified, who had reported testing positive for allergy to nickel prior to or after Essure placement (reporting frequency, 0.01% of all placements). The clinical events reported in this group of 91 testing positive for nickel allergy were similar to those reported in all 641 cases.

To date, there have been no credible reports of true, immediate-type anaphylactic or fatal allergic reactions to Essure.

7.12.4 Conclusion

While nickel allergy assessed by skin patch test is a common occurrence, positive skin patch test is not a reliable predictor of clinically significant reactions to implanted nickel-containing devices. Data from the clinical trials and literature indicate these reactions are relatively rare. Postmarketing monitoring indicates that the term ‘allergy’ appears to be often used in an unspecific manner.

7.13 Device removal

Essure inserts are intended to be left in place permanently. There are no well-controlled trials that have studied removal of inserts remote from placement. Clinical judgment as to the appropriate procedure for removal must be used. Additional information on device removal was recently added to the IFU.

7.13.1 Sponsor’s clinical data

Compiled data on device removals based on data from the phase II and pivotal study reports and from the ongoing TVU (16974) study are provided in [Table 7–17](#). In addition, there are 5 subjects (1 in the phase II study and 4 in the pivotal study) who underwent device removal

at the time of salpingectomy for sterilization but for whom this information is not documented in the above mentioned study reports, although it is documented in other legacy Conceptus files.

This table does not include hysteroscopic device removals that occurred immediately at the time of placement.

Table 7–17: Essure device removal in clinical trials

	Phase II study (STOP 10)	Pivotal study (STOP 2000)	ESS305 TVU (Study 16974)
Total device removals	11	20	11
Number attempted placement	227	507	597
Laparoscopic removal	4	5 ^a	7
Other removal	2 ^c	-	-
Removal with hysterectomy	5	15 ^b	4
	3 bleeding and/or pain 2 prolapse	9 bleeding and/or pain 6 other (1 Asherman's syndrome, 1 fibroids, and 4 missing)	1 pain 1 fibroids 1 endometriosis 1 bleeding

TVU = transvaginal ultrasound

^a inserts removed laparoscopically prior to IVF in 1 woman desiring pregnancy, 4 others at time of sterilization

^b 6 women with device fragments remaining in situ after hysterectomy

^c includes one woman with devices removed via cornual resection and one woman with devices removed via laparotomy due to pain.

7.13.2 Literature

A limited number of case reports indicate that the device may be removed hysteroscopically up to 7 weeks following placement if the proximal coils are easily visible and gentle traction is applied.(109, 81, 110) Case reports show that linear salpingostomy or salpingectomy via laparoscopy or laparotomy can be used to remove the insert. In post approval published literature, salpingostomy has been described from 35 days to 8 months post-placement. Laparoscopic salpingectomy has been described from 10 weeks post-placement and up to 4 years post-placement.(110) In some cases, a cornual resection of the proximal fallopian tube may be required for removal. In these cases, women should be counseled about the risk of hysterectomy in order to achieve hemostasis.(110)

Removal technique for perforated devices or those within the peritoneal cavity will depend on the location of the device. Localization should be confirmed prior to the surgical procedure. Use of intra-operative fluoroscopy or plain film is recommended to identify the location of the

insert(s) or pieces of the insert prior to and during surgery as this may prevent future operative interventions.(97)

7.13.3 Postmarketing monitoring

A search of the global safety database through 15 JUNE 2015 revealed 1,127 reports describing removal of Essure in association with, but not exclusively due to, complaints considered to be causally related to the device. This corresponds to a reporting rate of 0.11%. Most reports were spontaneous reports from US HCPs or consumers/ women.

Forty-five percent of device removals were associated with an unsatisfactory device location such as a perforation, embedment, dislocation or expulsion. Pain and bleeding disturbances were reported in 72.6% and 29.4% of these cases, respectively, indicating that the presence of these symptoms influenced the decision for removal. Pain and heavy bleeding were the most common symptoms associated with removal. Device breakage which can occur during the placement procedure was reported in 11.7% of removals which is the appropriate management of a broken device, deployment difficulty or any device suspected of malfunction. Allergy to metals was reported in 10% of removal cases.

Complications associated with device removal were reported in less than 5% of these cases, or 0.006% of the total population. Among medically confirmed cases, device breakage with removal of all fragments was reported in 2.0% of these cases. Suspected incomplete removal was reported in 1.1% and inability to remove the device as planned in 0.6% of cases. In 0.6% of cases the device could not be located. In order to minimize this occurrence, the IFU (Section 10.6, Appendix 6) recommends the availability of intra-operative imaging procedures.

7.13.4 Conclusion

The need for removal of Essure is overall infrequent. Device removals during the placement procedure represent appropriate recognition and management of unsatisfactory device physical properties or location. Device removal is commonly associated with symptoms of pain and bleeding. The least invasive method of removal that can be safely conducted should be performed. In situations where other existing gynecological diagnoses are present, concomitant removal of the uterus with the device may be warranted.

7.14 Death

Death represents the most serious outcome of any AE. It is important to distinguish between fatal outcomes secondary to the product under evaluation, specifically Essure, from fatal outcomes due to other medical conditions in women who have decided to use Essure for contraception. We have examined all reports that could potentially reflect fatal outcomes as discussed in this section.

7.14.1 Sponsor's clinical data

Two deaths occurred in clinical trials with Essure. One woman in the pivotal study (STOP 2000) secondary to leukemia and one woman in ESS305 TVU study (Study 16974) secondary to a myocardial infarction ten days following coronary artery bypass graft surgery. None of these events occurred around the time of the Essure placement or confirmation test. Both events were considered not to be related to Essure.

7.14.2 Literature

There were no reports from the literature on death in association with Essure placement.

7.14.3 Postmarketing monitoring

The sponsor's global safety database was searched through 15 JUNE 2015 for events with a fatal outcome. This process identified 22 cases containing non-fatal events including terms such as "bleeding to death" or "I felt like dying". These were excluded as they did not report fatal outcomes. Three reports were considered as invalid cases because they did not meet the criteria for case creation. Sixteen cases were excluded from this analysis as they referred to reports on loss of pregnancies (miscarriage, stillbirth, neonatal death). These reports are included in the analysis on pregnancy outcomes.

There are 12 postmarketing reports with fatal outcomes, all from the US. In 7 of these cases a cause of death is identified whereas in 5 no cause is reported. Of the cases with no reported cause of death, 3 were from social media, one from a television website, and one from an individual who heard it from another individual who saw it on television. In all of these cases, no specifics regarding the date or cause of death or the relation to a procedure was reported.

There are 7 cases in which the cause of death is reported. Three of these involve deaths related to anesthetic complications (2 local and one general anesthesia). A woman with severe pulmonary hypertension had a cardiac arrest during removal of a carotid catheter which was fatal. She had Essure inserted in general anesthesia. A woman with sleep apnea experienced a fatal outcome in her sleep. She had taken phenergan to facilitate sleep. One case of toxic shock syndrome secondary to the transfer of Group A streptococcus into the upper genital tract during hysteroscopy was reported.

Six of these seven cases were reported by a HCP. The seventh case reported from social media describes death due to an embolism during a hysterectomy. No further information was provided.

7.14.4 Conclusion

A comparison to the deaths with tubal sterilization performed laparoscopically is helpful in putting the fatal cases reported with Essure users into perspective. As discussed in the introductory section of this document, the leading causes of death with laparoscopic

sterilization were anesthesia and infection related. In Essure users there were 3 anesthesia related fatalities and 1 infection related fatality. The projected case fatality rate for tubal sterilization (6) was projected at ~9 per 100,000 with a rate of 1-2 per 100,000 that are directly procedure related. Given the approximately 1 million procedures performed, the rates reported for Essure are within these estimates.

8. Benefit/Risk Summary

Benefit

All female permanent birth control methods are highly effective. The following table presents comparative pregnancy rates by sterilization method, as recently compiled by the American College of Obstetricians and Gynecologists.(2)

Table 8–1: Pregnancy rates by sterilization method

Method	5-year / 1,000 procedures	10-year / 1,000 procedures	Ectopic / 1,000 procedures
Postpartum partial salpingectomy	6.3	7.5	1.5
Bipolar coagulation ^a	16.5	24.8	17.1
Silicone band methods	10.0	17.7	7.3
Spring clip	31.7	36.5	8.5
Hysteroscopy (Essure) ^b	1.64	-	-
Vasectomy	11.3	-	No association

^a Secondary analysis of 5-year failure rates with bipolar coagulation performed in different decades found that failure was significantly lower in later periods, reflecting improved technique with the method: 19.5 per 1,000 procedures for 1978 – 1982 vs 6.3 per 1,000 procedures for 1985 – 1987.

^b US FDA recommended projections based on Bayesian statistical analysis.(111)

In a review of the literature, Grunert found that the tubal patency rate following incisional tubal ligation was 3.2% with a subsequent pregnancy rate of 0.8%. Grunert concluded that “although there may be a failure of absolute physical occlusion of the tubes, this cannot be directly equated with failure of the sterilization”.(112)

Clinical trials including a 5-year follow-up study did not show any pregnancies after Essure placement. Later clinical trials demonstrated low rates of unintended pregnancy.

Clinical studies in the literature confirm that hysteroscopic sterilization with the Essure device is a very effective method of contraception with reported efficacy in the literature and clinical trials of >99%.

Based on postmarketing monitoring and considering its inherent limitations, pregnancy rates with the marketed product have been consistent with rates determined in clinical trials and as found in the literature.

No contraceptive method is 100% effective and contraceptive failures with Essure have been reported in clinical studies and post-marketing. Data suggest that compliance with the requirement to undergo a confirmation test three month after the Essure placement is critical to achieve the benefits inherent to the Essure method. This is addressed in the IFU (Section 10.6, Appendix 6) and physician training.

Placement rate is very important when physicians counsel patients prior to the attempted procedure, as patients need to be aware of the possibility placement may not occur. Still, with the currently available ESS305 device, both the clinical trial data as well as several external studies show a placement rate of at least 95%.

Both the clinical trials and external studies demonstrate overall satisfaction (>92%) with both the Essure placement procedure as well as ongoing use of the inserts. When compared to laparoscopic tubal ligation at 90 days post procedure, Essure compared very favorably with regard to overall satisfaction, speed of recovery, and choice of chosen procedure.(49)

Risk

Prior to Essure, the primary option for women desiring permanent birth control at a time remote from pregnancy was laparoscopic BTL. While BTL is regarded as safe and effective,(2) it is an invasive procedure requiring entry into the peritoneal cavity and general anesthesia. Many of the risks of the procedure are related to the required surgical entry into the abdominal cavity, including trochar or Verres needle injuries to intra-abdominal structures, bleeding, sepsis, and bowel injury. The use of general anesthesia, insufflation of the abdomen with CO₂ gas, and the use of steep Trendelenburg position sometimes needed to access pelvic structures also increase the risk of all laparoscopic surgeries. There may be increased risks of aspiration, respiratory dysfunction, and cardiovascular dysfunction during laparoscopy. The ability to perform the Essure procedure without general anesthesia is therefore a clinically relevant consideration.

The main risks with Essure as described in the IFU (Section 10.6, Appendix 6) and PIB (Section 10.7, Appendix 7) include: uterine, or tubal perforation, pregnancy (including ectopic pregnancy), infection, menstrual changes, pain, allergy/hypersensitivity reactions and lack of contraceptive effect due to improper placement or expulsion/migration of the inserts.

Benefit/Risk assessment

The Essure system provides permanent birth control without invasive surgery or the need for general anesthesia, and their attendant risks. It is a highly effective form of permanent birth

control, even when taking into account the follow-up protocol required. The risk profile is well characterized and compares favorably with the risk profile of BTL.

Risks related to Essure and the placement procedure are mitigated through the detailed information presented in the IFU and PIB, a robust multi-step physician training program, and restricted distribution of Essure to only those physicians who have completed this training. In addition, there are numerous physician support programs in place. We continue to conduct postmarketing surveillance and to work closely with Health Authorities, such as the FDA, to identify and address any changes to the benefit/risk profile.

The safety profile of Essure is consistent with other methods of female permanent birth control and may present advantages for a number of women.

Data from the sponsor's extensive clinical program, comprehensive review of peer-reviewed literature, and continuous postmarketing activities support the safety and efficacy of Essure.

Overall, this assessment concludes that the benefit/risk profile of Essure remains positive.

9. References

1. National Survey of Family Growth, 2014 (<http://www.cdc.gov/nchs/nsfg.htm>).
2. ACOG Practice Bulletin: Clinical Management Guidelines for Obstetrician-Gynecologists – Benefits and Risks of Sterilization. Number 133, Feb 2013.
3. Peterson HB, Xia Z, Hughes JM, Wilcox LS, Tylor LR, Trussell J. The risk of pregnancy after tubal sterilization: findings from the U.S. Collaborative Review of Sterilization. *Am J Obstet Gynecol.* 1996;174:1161-8; discussion 1168-70.
4. Jamieson DJ, Hillis SD, Duerr A, Marchbanks PA, Costello C, Peterson HB. Complications of interval laparoscopic tubal sterilization: findings from the United States Collaborative. *Review of Sterilization. Obstet Gynecol.* 2000;96:997-1002.
5. Layde PM, Peterson HB, Dicker RC, DeStefano F, Rubin GL, Ory HW. Risk factors for complications of interval tubal sterilization by laparotomy. *Obstet Gynecol.* 1983;62:180-4.
6. Escobedo LG, Peterson HB, Grubb GS, Franks AL. Case-fatality rates for tubal sterilization in U.S. hospitals, 1979 to 1980. *Am J Obstet Gynecol.* 1989 Jan;160(1):147-50.

7. Huber AW, Mueller MD, Ghezzi F, Cromi A, Dreher E, Raio L. Tubal sterilization: complications of laparoscopy and minilaparotomy. *Eur J Obstet Gynecol Reprod Biol.* 2007 Sep;134(1):105-9.
8. Jansen FW, Kolkman W, Bakkum EA, de Kroon CD, Trimbos-Kemper CM, Trimbos JB. Complications of laparoscopy: An inquiry about closed- versus open-entry technique. *Am J Obstet Gynecol* 2004 Mar;190(3):634-8.
9. Llarena NC, Shah AB, Milad MP. Bowel injury in gynecologic laparoscopy: a systematic review. *Obstet Gynecol.* 2015 Jun;125(6):1407-17.
10. Chan LM, Westhoff CL. Tubal sterilization trends in the United States. *Fertil Steril.* 2010 Jun;94(1):1-6.
11. Ülker K, Anuk T, Bozkurt M, Karasu Y. Large bowel injuries during gynecological laparoscopy. *World J Clin Cases.* 2014 Dec 16;2(12):846-51.
12. Cottin V, Delafosse B, Viale J-P Gas embolism during laparoscopy: A report of seven cases in patients with previous abdominal surgical history. *Surgical Endoscopy* 1996; 10:166-169.
13. Wadhwa RK, McKenzie R, Wadhwa SR, Katz DL, Byers JF. Gas embolism during laparoscopy. *Anesthesiology*1978; 48:74-76.
14. Hynes SR, Marshall RL. Venous gas embolism during gynaecological laparoscopy. *Can J Anaesth* 1992;39:748-749
15. Valle RF, Carignan CS, Wright TC; STOP Prehysterectomy Investigation Group. Tissue response to the STOP microcoil transcervical permanent contraceptive device: results from a prehysterectomy study. *Fertil Steril.* 2001 Nov;76(5):974-80.
16. Greisler HP, Henderson SC, Lam TM. Basic fibroblast growth production in vitro by macrophages exposed to Dacron and polyglactin 910. *J Biomater Sci Polym Ed (Netherlands)*, 1993; 4 (5) p. 415-30.
17. Estridge TD, Feldman DS. Quantification of vascular ingrowth into Dacron velour. *J Biomater Appl (USA)*. Oct 1991, 6 (2) p. 157-69.
18. Adloff M, Arnaud JP. Surgical management of large incisional hernias by an intraperitoneal Mersilene mesh and an aponeurotic graft. *Surg Gynecol Obstet* 1987;165:204.
19. Pans A, Pierard GE. A comparison of intraperitoneal prosthesis for the repair of abdominal muscular wall defects in rats. *Eur Surg Res* 1992; 24:54.

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30. Janse JA, Goedegebuure RS, Veersema S, Broekmans FJ, Schreuder HW. Hysteroscopic sterilization using a virtual reality simulator: assessment of learning curve. *J Minim Invasive Gynecol.* 2013 Nov-Dec;20(6):775-82.
31. Bahk CY, Goshgarian M, Donahue K, Freifeld CC, Menone CM, Pierce CE, Rodriguez H, Brownstein JS, Furberg R, Dasgupta N. Increasing patient engagement in pharmacovigilance through online community outreach and mobile reporting applications: An analysis of adverse event reporting for the Essure device in the US. *Pharm Med DOI 10.1007/s40290-015-0106-6* (published, 05 Aug 2015).
32. 21CFR 803, Medical Device Reporting.
33. Levie MD and Chudnoff SG. Prospective analysis of office-based hysteroscopic sterilization. *J Minim Invasive Gynecol.* 2006 Mar-Apr;13(2):98-101.
34. Anderson TL, Yunker AC, Scheib SA, Callahan TL. Hysteroscopic sterilization success in outpatient vs office setting is not affected by patient or procedural characteristics. *J Minim Invasive Gynecol.* 2013 Nov-Dec;20(6):858-63.
35. Chapa HO. Vaginoscopy compared to traditional hysteroscopy for hysteroscopic sterilization. A randomized trial. *J Repd Med* 2015 Jan-Feb;60(1-2):43-47.
36. Shavell V, Abdallah ME, Diamond MP, Berman JM. Placement of a permanent birth control device at a university medical center. *J Repd Med* 2009 Apr; 54(4):218-222.
37. Povedano B, Arjona JE, Velasco E, Monserrat JA, Lorente J, Castelo-Branco C. Complications of hysteroscopic Essure[®] sterilisation: report on 4306 procedures performed in a single centre. *BJOG.* 2012 Jun;119(7):795-9.
38. Ubeda A, Labastida R, Dexeus S. Essure: a new device for hysteroscopic tubal sterilization in an outpatient setting. *Fertility & Sterility, Volume 82, Issue 1, July 2004, Pages 196–199.*
39. Shavell V, Abdallah ME, Diamond MP, Kmak DC, Berman JM. Post-Essure hysterosalpingography compliance in a clinic population. *J Minim Invasive Gynecol.* 2008 Jul-Aug;15(4):431-4.
40. Leyser-Whalen O, Berenson AB. Adherence to hysterosalpingogram appointments following hysteroscopic sterilization among low-income women. *Contraception.* 2013 Dec;88(6):697-9.
41. Savage UK, Masters SJ, Smid MC, Hung YY, Jacobson GF. Hysteroscopic sterilization in a large group practice: experience and effectiveness. *Obstet Gynecol.* 2009 Dec; 114(6):1227-31.

42. Howard DL, Wall J, Strickland JL. What are the factors predictive of hysterosalpingogram compliance after female sterilization by the Essure procedure in a publicly insured population? *Matern Child Health J.* 2013 Dec;17(10):1760-7.
43. Guiahi M, Goldman KN, McElhinney MM, Olson CG. Improving hysterosalpingogram confirmatory test follow-up after Essure hysteroscopic sterilization. *Contraception.* 2010 Jun;81(6):520-4.
44. Mahmud S, Pereira N, Taylor KC, Ekbladh LE. Improving adherence to hysterosalpingography after hysteroscopic sterilization using an electronic reminder. *J Minim Invasive Gynecol.* 2015 Feb;22(2):250-4.
45. Sinha D, Kalathy V, Gupta JK. The feasibility, success and patient satisfaction associated with outpatient hysteroscopic sterilisation. *BJOG.* 2007 Jun;114(6):676-83.
46. Miño M, Arjona JE, Cordon J, Pelegrin B, Povedano B, Chacon E. Success rate and patient satisfaction with the Essure sterilisation in an outpatient setting: a prospective study of 857 women. *BJOG.* 2007 Jun;114(6):763-6.
47. Anderson S, Eriksson S, Mints M. Hysteroscopic female sterilization with Essure in an outpatient setting. *Acta Obstet Gynecol Scand.* 2009;88(6):743-6.
48. Levie M, Weiss G, Kaiser B, Daif J, Chudnoff SG. Analysis of pain and satisfaction with office-based hysteroscopic sterilization. *Fertil Steril.* 2010 Sep;94(4):1189-94.
49. Duffy S, Marsh F, Rogerson L, Hudson H, Cooper K, Jack S, Hunter D, Philips G. Female sterilisation: a cohort controlled comparative study of ESSURE versus laparoscopic sterilisation. *BJOG.* 2005 Nov;112(11):1522-8.
50. Levie M; Chudnoff SG. A comparison of novice and experienced physicians performing hysteroscopic sterilization: an analysis of an FDA-mandated trial. *Fertility and Sterility.* 2011 Sep; 96(3):643-648.
51. Fernandez H. Tubal sterilization: pregnancy rates after hysteroscopic versus laparoscopic sterilization in France, 2006-2010. *Eur J Obstet Gynecol Reprod Biol.* 2014 Sep;180:133-7.
52. Deraleu and Heinlein. Contraceptive failures associated with hysteroscopic sterilization. *J Minim Inv Gynecol* 2012:S19 (S1-S35).
53. Gariepy AM, Creinin MD, Smith KJ, Xu X. Probability of pregnancy after sterilization: a comparison of hysteroscopic versus laparoscopic sterilization. *Contraception.* 2014 Aug;90(2):174-81.

54. Regan L, Rai R. Epidemiology and the medical causes of miscarriage. *Baillieres Best Pract Res Clin Obstet Gynaecol*. 2000 Oct;14(5):839-54.
55. Prine LW, MacNaughton H. Office management of early pregnancy loss. *Am Fam Physician*. 2011 Jul 1;84(1):75-82.
56. Maconochie N, Doyle P, Prior S, Simmons R. Risk factors for first trimester miscarriage--results from a UK-population-based case-control study. *BJOG*. 2007 Feb;114(2):170-86.
57. Smith G, Rafuse C, Anand N, Brennan B, Connors G, Crane J, Fraser W, Gratton R, Moutquin JM, Scott H, Schneider C, Walker M. Prevalence, management, and outcomes of preterm prelabour rupture of the membranes of women in Canada. *J Obstet Gynaecol Can*. 2005 Jun;27(6):547-53.
58. Saraiya M, Berg CJ, Shulman H, Green CA, Atrash HK. Estimates of the annual number of clinically recognized pregnancies in the United States, 1981-1991. *Am J Epidemiol*. 1999 Jun 1;149(11):1025-9.
59. Van Den Eeden SK, Shan J, Bruce C, Glasser M. Ectopic pregnancy rate and treatment utilization in a large managed care organization. *Obstet Gynecol*. 2005 May;105(5 Pt 1):1052-7.
60. Peterson HB. Sterilization. *Obstet Gynecol*. 2008;Jan 111(1):189-203.
61. Veersema S, Mijatovic V, Dreyer K, Schouten H, Schoot D, Emanuel MH, Hompes P, Brölmann H. Outcomes of pregnancies in women with hysteroscopically placed micro-inserts in situ. *J Minim Invasive Gynecol*. 2014 May-Jun;21(3):492-7.
62. Ory EM, Hines RS, Cleland WH, Rehberg JF. Pregnancy after microinsert sterilization with tubal occlusion confirmed by hysterosalpingogram. *Obstet Gynecol*. 2008 Feb;111(2 Pt 2):508-10.
63. Moses AW, Burgis JT, Bacon JL, Risinger J. Pregnancy after Essure placement: report of two cases. *Fertil Steril*. 2008 Mar;89(3):724.
64. Hastings-Tolsma M, Nodine P, Teal SB, Embry J. Pregnancy outcome after transcervical hysteroscopic sterilization. *Obstet Gynecol*. 2007 Aug;110(2 Pt 2):504-6.
65. Arjona JE1, Serrano JJ, Povedano B, Carrasco S, Castelo-Branco C. Unintended pregnancy after long-term Essure microinserts placement. *Fertil Steril*. 2010 Dec;94(7):2793-5.

66. Rosenfield RB, Stones RE, Coates A, Matteri RK, Hesla JS. Proximal occlusion of hydrosalpinx by hysteroscopic placement of microinsert before in vitro fertilization-embryo transfer. *Fertil Steril*. 2005 May;83(5):1547-50.
67. Arora P, Arora RS, Cahill D. Essure® for management of hydrosalpinx prior to in vitro fertilisation-a systematic review and pooled analysis. *BJOG*. 2014 Apr;121(5):527-36.
68. Legendre G, Moulin J, Vialard J, Ziegler DD, Fanchin R, Pouly JL, Watrelot A, Belaisch Allart J, Massin N, Fernandez H. Proximal occlusion of hydrosalpinges by Essure® before assisted reproduction techniques: a French survey. *Eur J Obstet Gynecol Reprod Biol*. 2014 Oct;181:300-4.
69. Mijatovic V, Veersema S, Emanuel MH, Schats R, Hompes PGA. Essure hysteroscopic tubal occlusion device for the treatment of hydrosalpinx prior to in vitro fertilization-embryo transfer in patients with a contraindication for laparoscopy. *Fertil Steril* 2010; 93: 1338-42.
70. Mijatovic V, Dreyer K, Emanuel MH, Schats R, Hompes PGA. Essure hydrosalpinx occlusion prior to IVF-ET as an alternative to laparoscopic salpingectomy. *Eur J Obstet Gynecol Repro Biol* 2012 vol 161 42-25.
71. Ozgur K, Bulut H, Berkkanoglu M, Coetzee K. ICSI pregnancy outcomes following hysteroscopic placement of Essure devices for hydrosalpinx in laparoscopic contraindicated patients. *Repd BioMed* 2014;29(1):113-118.
72. Gerritse MBE, Veersema S, Timmermans A, et al. Incorrect position of Essure microinserts 3 months after successful bilateral placement. *Fertil Steril*. 2009;91(3):930e1-5.
73. Borley J, Shabajee N, Tan TL. A kink is not always a perforation: Assessing Essure hysteroscopic sterilization placement. *Fertil Steril*. 2011;95:2429.e15-17.
74. Hur HC, Mansuria SM, Chen BA, et al. Laparoscopic management of hysteroscopic Essure sterilization complications: Report of 3 cases. *J Minim Invasive Gynecol*. 2008;15:362-65.
75. Moawad N, Mansuria S. Essure perforation and chronic pelvic pain. *J Minim Invasive Gynecol*. 2011;18(3):285-6.
76. Thoma V, Chua I, Garbin O, et al. Tubal perforation by Essure microinsert. *J Minim Invasive Gynecol*. 2006;13:161-63.

77. Ricci G, Restaino S, Di Lorenzo G, et al. Risk of Essure microinsert abdominal migration: Case report and review of literature. *Ther Clin Risk Manag.* 2014;10:963-68.
78. Chohan L, Kilpatrick CC. Laparoscopic removal of perforated tubal occlusion micro-inserts. *J Minim Invasive Gynecol.* 2014;21(4):534.
79. Howard DL, Christenson PJ, Strickland JL. Use of intraoperative fluoroscopy during laparotomy to identify fragments of retained. Essure microinserts: Case report. *J Minim Invasive Gynecol.* 2012;19:667-70.
80. Pyke R, Blackwood LR. Complication of the Essure implant sterilization procedure: A case report. *J GynecolSurg.* 2008;24(1):37-42.
81. Langenveld J, Veersema S, Bongers M, et al. Tubal perforation by Essure: Three different clinical presentations. *Fertil Steril.* 2008;90:2011e5-e10.
82. Al-Safi ZA, Shavell VI, Hobson DT, et al. Analysis of adverse events with Essure hysteroscopic sterilization reported to the manufacturer and user facility device experience database. *J Minim Invasive Gynecol* 2013;20(6):825-9.
83. Belotte J, Shavell VI, Awonuga AO, et al. Small bowel obstruction subsequent to Essure microinsert sterilization: a case report. *Fertil Steril* 2011;96:e4-e6.
84. Braginsky L, George ST, Locher SR. Management of perforated Essure with migration into small and large bowel mesentery. *J Minim Invasive Gynecol.* 2015;22(3):504-8.
85. Mantel H, Wijma J, Stael A. Small bowel obstruction and perforation after Essure sterilization. *Contraception* 2013;87:121-123.
86. Riley K, Beltran F, Stewart D, et al. Bowel perforation after placement of tubal occlusion contraceptive. *ObstetGynecol*, 2015.
87. Latthe P, Mignini L, Gray R, Hills R, Khan K. Factors predisposing women to chronic pelvic pain: systematic review. *BMJ.* 2006 Apr 1;332(7544):749-55.
88. Frishman GN, Brest NA. Hulka Clip application as a potential cause of chronic pelvic pain. *Contraception.* 1992 Apr;45(4):325-7.
89. Hiemstra E, Weijenborg PT, Jansen FW. Management of chronic pelvic pain additional to tubal sterilization. *J Psychosom Obstet Gynaecol.* 2008 Sep;29(3):153-6.

90. Arjona JE, Miño M, Cordon J, Povedano B, Pelegrin B, Castelo-Branco C. Satisfaction and tolerance with office hysteroscopic tubal sterilization. *Fertil Steril*. 2008 Oct;90(4):1182-6.
91. Syed R, Levy J, Childers ME. Pain associated with hysteroscopic sterilization. *JSLs*. 2007 Jan-Mar;11(1):63-5.
92. Yunker AC, Ritch JM, Robinson EF, Golish CT. Incidence and risk factors for chronic pelvic pain after hysteroscopic sterilization. *J Minim Invasive Gynecol*. 2015 Mar-Apr;22(3):390-4.
93. Rudin A, Wölner-Hanssen P, Hellbom M, Werner MU. Prediction of post-operative pain after a laparoscopic tubal ligation procedure. *Acta Anaesthesiol Scand*. 2008 Aug;52(7):938-45.
94. Arjona Berral JE, Rodríguez Jiménez B, Velasco Sánchez E, Povedano Cañizares B, Monserrat Jordan J, Lorente Gonzalez J, Castelo-Branco C. Essure[®] and chronic pelvic pain: a population-based cohort. *J Obstet Gynaecol*. 2014 Nov;34(8):712-3.
95. Conover MM, Howell JO, Wu JM, Kinlaw AC, Dasgupta N, Jonsson Funk M. Incidence of opioid-managed pelvic pain after hysteroscopic sterilization versus laparoscopic sterilization, US 2005-2012. *Pharmacoepidemiol Drug Saf*. 2015 Aug;24(8):875-84.
96. Mahmoud MS, Fridman D, Merhi ZO. Subserosal misplacement of Essure device manifested by late-onset acute pelvic pain. *Fertil Steril*. 2009 Dec;92(6):2038.
97. Adelman MR, Dassel MW, Sharp HT. Management of complications encountered with Essure hysteroscopic sterilization: a systematic review. *J Minim Invasive Gynecol*. 2014 Sep-Oct;21(5):733-43.
98. DeStefano F, Perlman JA, Peterson HB, et al. Long-term risk of menstrual disturbances after tubal sterilization. *Am J Obstet Gynecol* 1985; 152:835-41.
99. Foulkes J, Chamberlain G. Effects of sterilization on menstruation. *South Med J* 1985;78:544-7.
100. Williams et al, 1951 WILLIAMS EL, JONES HE, MERRILL RE. The subsequent course of patients sterilized by tubal ligation; a consideration of hysterectomy for sterilization. *Am J Obstet Gynecol*. 1951 Feb;61(2):423-6.
101. Gentile GP, Kaufman SC, Helbig DW. Is there any evidence for a post-tubal sterilization syndrome? *Fertil Steril*. 1998 Feb;69(2):179-86.

102. Peterson HB, Jeng G, Folger SG, Hillis SA, Marchbanks PA, Wilcox LS; U.S. Collaborative Review of Sterilization Working Group. The risk of menstrual abnormalities after tubal sterilization. U.S. Collaborative Review of Sterilization Working Group. *N Engl J Med.* 2000 Dec 7;343(23):1681-7.
103. Rasmussen 1991 Rasmussen BK, Jensen R, Schroll M, Olesen J. Epidemiology of headache in a general population - a prevalence study. *J Clin Epidemiol* 1991;44:1147-1157.
104. Scher A, Stewart WF, Liberman J, Lipton RB. Prevalence of frequent headache in a population sample. *Headache.* 1998;38:497-506.
105. Xia E, Chen F, Zhonghua Fu. Chan Ke Za Zhi. Severe complications of hysteroscopic surgeries: an analysis of 35 cases. 2001 Oct;36(10):596-9 (Abstract only; Article in Chinese).
106. Lane A, Tyson A, Thurston E. Providing re-ESSURE-ance to the nickel allergic patient considering hysteroscopic sterilization. *J Minim Invasive Gynecol.* 2015 Aug 7: S1553-4650(15).
107. Peter I, Michael V, Hugo VE. Evaluation of nickel allergy afer hysteroscopic Essure sterilization: risk or daily practice? Preliminary results. *Gynecol Surg* 2011;8(Suppl 1):S1-S225, Abstract FC.08.05.
108. Schram SE, Warshaw EM, Laumann A. Nickel hypersensitivity: a clinical review and call to action. *Int J Dermatol.* 2010 Feb;49(2):115-25.
109. Lannon BM, Lee SY. Techniques for removal of the Essure hysteroscopic tubal occlusion device. *Fertil Steril.* 2007 Aug;88(2):497.
110. Albright CM, Frishman GN, Bhagavath B. Surgical aspects of removal of Essure microinsert. *Contraception.* 2013 Sep;88(3):334-6.
111. Baskinski CM. A review of clinical data for currently approved hysteroscopic sterilization procedures. *Rev Obstet & Gynecol* 2010;3(3):101-110.
112. Grunert GM. Late tubal patency following tubal ligation. *Fertil Steril* 1981 Apr;35(4):406-8.

10. Appendices

10.1 Appendix 1: Essure design changes

Key design changes since original US FDA Premarket Approval

Device	Insert	Catheter	Clinical data	Key reason for change
STOP	Original approved design	Braided delivery catheter	Phase II (STOP10) and pivotal studies (STOP2000)	Original PMA approval
ESS205	No change	Coil Catheter	Post approval study (ESS205 PAS) conducted to validate delivery system changes	Original launched design Improved bilateral placement rate
ESS305 Currently marketed model	Similar to ESS205, removal of taper "pigtail" at proximal end of outer coil	Modified detachment and deployment mechanism; replaced "unscrew" step with second roll-back	Post approval study (ESS305 PAS) conducted to validate delivery system changes	(a) Improved bilateral placement rate (b) Improved ease of use by simplifying placement procedure (reduced the total number of operational steps required for insert release by eliminating counter-clockwise turns). (c) Improve reliability (reduce detachment / disengagement difficulties between insert and delivery system)

STOP = Selective Tubal Occlusion Product; PMA = Pre-market approval; PAS = post approval study.

10.2 Appendix 2: Literature on Essure placement rates

Essure placement rate in real-world setting

Study	Location	Design	Findings
<p>Satisfaction and tolerance with office hysteroscopic tubal sterilization. Fertil Steril. 2008 Oct;90(4):1182-6. doi: 10.1016/j.fertnstert.2007.08.007. Epub 2008 Jan 16. Arjona JE, Miño M, Cerdón J, Povedano B, Pelegrin B, Castelo-Branco C. http://www.ncbi.nlm.nih.gov/pubmed/18201703</p>	Spain	Prospective analysis of case series; Gynecology department in a teaching hospital; A total of 1,630 women, January 2003 to June 2006.	Out of a total of 1,630, fifteen women (<1%) dropped out because of a failure in the procedure (insertion success rate >99%).
<p>Hysteroscopic sterilization success in outpatient vs office setting is not affected by patient or procedural characteristics. J Minim Invasive Gynecol. 2013 Nov-Dec;20(6):858-63. doi: 10.1016/j.jmig.2013.05.020. Epub 2013 Jul 23. Anderson TL, Yunker AC, Scheib SA, Callahan TL. http://www.ncbi.nlm.nih.gov/pubmed/23891205</p>	USA (Vanderbilt, Nashville, TN)	Retrospective cohort analysis 638 women who underwent hysteroscopic sterilization between July 1, 2005 and June 30, 2011.	The overall success in bilateral placement of inserts was 95.5%, with no significant difference in so far as the setting in which the procedure was performed (office vs OR), anesthesia verses no anesthesia, surgeon experience level, obesity.
<p>Hysteroscopic female sterilization with Essure in an outpatient setting. Andersson S, Eriksson S, Mints M. Acta Obstet Gynecol Scand. 2009 88(6):743--6 http://www.ncbi.nlm.nih.gov/pubmed/19412802</p>	Sweden	61 women underwent outpatient hysteroscopic sterilization between 2002 and 2007.	Fifty-eight (95%) women were sterilized according to this method. Successful bilateral device placement was achieved in 52 women (85%) during the first attempt and in six (9.8%) during the second, overall 58/61 (95.1%). Three women declined a second attempt and underwent laparoscopic sterilization.
<p>Vaginoscopy compared to traditional hysteroscopy for hysteroscopic sterilization. A randomized trial. Chapa HO. J Reprod Med. 2015 Jan-Feb;60(1-2):43-7. http://www.ncbi.nlm.nih.gov/pubmed/25745750</p>	USA (TX)	Single center, randomized, prospective, n=90.	Successful bilateral placement at first attempt in 86/90 (95.4%) overall (including 3 subjects in whom vaginoscopy was unsuccessful and traditional hysteroscopy utilized).
<p>Initial Asian experience in hysteroscopic sterilisation using the Essure permanent birth control device. BJOG. 2005 Sep;112(9):1322-7.</p>	Singapore.	A retrospective study. 80 women seeking permanent birth control. Analysis of initial experience with the Essure permanent birth control over a	Bilateral placement at first attempt: 72 (90%); bilateral placement at second attempt 5(6.2%); overall bilateral placement (1st and 2nd attempts) 77/80 (96.2%).

Study	Location	Design	Findings
Chern B, Siow A. http://www.ncbi.nlm.nih.gov/pubmed/16101615		period of 18 months from 22 June 2001 to 22 December 2002.	
Placement of a permanent birth control device at a university medical center. J Reprod Med. 2009 Apr;54(4):218-22. Shavell VI ¹ , Abdallah ME, Diamond MP, Berman JM. http://www.ncbi.nlm.nih.gov/pubmed/19438163	USA (Detroit Medical Center, Detroit, MI)	Single center, retrospective study, 316 women, January 1, 2003, to June 30, 2007.	The rate of successful placement of the Essure was 92.1% (294/316). Of the 22 placement failures, 11 were attributed to difficulty visualizing the tubal ostia
Female sterilisation: a cohort controlled comparative study of ESSURE versus laparoscopic sterilisation. BJOG. 2005 Nov;112(11):1522-8. Duffy S, Marsh F, Rogerson L, Hudson H, Cooper K, Jack S, Hunter D, Philips G. http://www.ncbi.nlm.nih.gov/pubmed/16225573	Multicenter, UK	Cohort controlled comparative study. The day surgery and outpatient unit of three large UK hospitals. 59 women underwent Essure placement.	Bilateral placement first attempt 76% (45/59); The overall bilateral placement rate for the ESSURE device was 81% (48/59). For 3 patients the investigators were unable to enter the uterine cavity with the hysteroscope
What are the factors predictive of hysterosalpingogram compliance after female sterilization by the Essure procedure in a publicly insured population? Matern Child Health J. 2013 Dec;17(10):1760-7 Howard DL, Wall J, Strickland JL http://www.ncbi.nlm.nih.gov/pubmed/23212398	USA (2 campuses of Truman Medical Center, Kansas City, Missouri)	Multi-site, 136 patients presenting for Essure placement between January 1 2005 and December 31 2010	Placement of the Essure permanent birth control device was attempted in 136 patients, 131 of whom had successful placement (96.3 %) [includes one unilateral placement in patient with salpingectomy].
Analysis of pain and satisfaction with office-based hysteroscopic sterilization. Fert Steril. 2010 Sep4(4):1189-94 Levie M, Weiss G, Kaiser B, Daif J, Chudnoff SG. http://www.ncbi.nlm.nih.gov/pubmed/19683232	USA (NY)	From June 2003 to June 2006, 209 patients were recruited, at a faculty practice office at an inner-city urban medical center.	Successful placement occurred in 198/209 (95%)
Prospective analysis of office-based hysteroscopic sterilization. Levie MD, Chudnoff SG J Minim Invasive Gynecol. 2006;13(2):98-101. http://www.ncbi.nlm.nih.gov/pubmed/16527710	USA, NY	Prospective, longitudinal analysis, of 102 women at university outpatient office.	Bilateral placement of the device was successful in 98 (96%) of 102 patients; 92 had f/u HSG, 1 pregnancy (luteal phase pregnancy-pregnant at procedure or shortly after)
A comparison of novice and experienced physicians performing hysteroscopic sterilization: an analysis of an FDA-mandated trial. Levie M; Chudnoff SG Fertility and Sterility. 2011 Sep; 96(3):643-648.e1. http://www.ncbi.nlm.nih.gov/pubmed/?term=21782168	Multicenter USA (Subanalysis of ESS305 PAS)	Multicenter prospective cohort study. Seventy-six sites throughout the U.S. comprising community hospitals, teaching institutions, surgery centers, and office-based practices. A total of 578 women seeking	N=578, A total of 625 patients were recruited with 578 eligible cases performed by 37 newly trained and 39 experienced physicians: 562/578 (97.2%) overall had successful placement, with 98% successful placement rate versus 96.1% for experienced versus novice physicians,

Study	Location	Design	Findings
<p>Success rate and patient satisfaction with the Essure sterilisation in an outpatient setting: a prospective study of 857 women. Miño M, Arjona JE, Cordón J, Pelegrin B, Povedano B, Chacon E. BJOG 2007 Jun;114(6):763-6 http://www.ncbi.nlm.nih.gov/pubmed/?term=17516970</p>	Spain	hysteroscopic sterilization. Single center, retrospective, 857 women with Essure system, inserted in an outpatient setting, January 2003 to January 2005.	respectively. At the first attempt, 96.5% of insertions (827 women) were successful, which included 812 bilateral insertions and 15 unilateral (14 of whom had previous salpingectomy and one had a unicornuate uterus). Bilateral success was achieved at the second attempt in 18 women (2.1%). In the 12 women in whom insertion was not achieved, five were bilateral and seven were unilateral failures. Overall bilateral placement rate 845/857=98.6%.
<p>Complications of hysteroscopic Essure® sterilisation: report on 4306 procedures performed in a single centre. Povedano B, Arjona JE, Velasco E, Monserrat JA, Lorente J, Castelo-Branco C. BJOG 2012. Jun;119(7):795-9 http://www.ncbi.nlm.nih.gov/pubmed/22360159</p>	Spain	Retrospective 7-year study, Office hysteroscopic unit in a teaching hospital, 4306 women sterilized from March 2003 and June 2010.	Of the 4306 women undergoing hysteroscopic Essure® sterilisation, 4242 (98.5%) were successfully completed, of which 4075 were successfully sited at one attempt (96%), without the need for a second-stage procedure at a later date.
<p>The feasibility, success and patient satisfaction associated with outpatient hysteroscopic sterilisation. D. Sinha, V. Kalathy, J.K. Gupta BJOG. 2007 Jun;114(6):676-83. http://www.ncbi.nlm.nih.gov/pubmed/17516957</p>	UK	Prospective cohort study, large teaching hospital, 112 women enrolled August 2002 and June 2006.	Successful bilateral tubal placement of the Essure microinserts was achieved in 103/112 (92%, 95% confidence interval 85-96%) women.
<p>Post-Essure hysterosalpingography compliance in a clinic population. Shavell VI, Abdallah ME, Diamond MP, Kmak DC, Berman JM. J Minim Invasive Gynecol. 2008 Jul-Aug;15(4):431-4. doi: 10.1016/j.jmig.2008.03.004. Epub 2008 Apr 28. http://www.ncbi.nlm.nih.gov/pubmed/18588851</p>	US (Wayne State, Detroit)	83 women, Retrospective chart review, University teaching hospital, placement attempted during January 2003 through June 2007.	Placement of the Essure permanent birth control device was attempted in 83 patients, of which 79 were successfully completed (95.2%).
<p>Hysteroscopic sterilization in a large group practice: experience and effectiveness. Savage UK, Masters SJ, Smid MC, Hung YY, Jacobson GF Obstet Gynecol. 2009 Dec;114(6):1227-31. doi: 10.1097/AOG.0b013e3181c2a10d.</p>	North California, Kaiser Permanente (USA)	Retrospective review, 884 women, 118 physicians at 30 different facilities, from January 2004 to December 2006, f/u through December 31, 2008	The initial placement attempt was successful in 850 (96.2%), failed in 31 (3.5%), and could not be determined in three (0.3%). 13 patients underwent two attempts (<i>eventual bilateral placement rate not given, just initial</i>).

Study	Location	Design	Findings
http://www.ncbi.nlm.nih.gov/pubmed/19935023			
Essure: a new device for hysteroscopic tubal sterilization in an outpatient setting A. Ubeda, R. Labastida, S. Dexeus Fertility & Sterility, Volume 82, Issue 1, July 2004, Pages 196–199 http://www.sciencedirect.com/science/article/pii/S0015028204006843	Barcelona, Spain	Prospective, observational study, Private university hospital, 85 women, enrolled July 2002 - July 2003.	Successful placement was achieved in 81 (95%) patients, including one case of two-step procedure (placement of devices in two consecutive menstrual cycles because of unilateral tubal spasm), one unicornuate uterus, and a woman in whom only one device could be placed and whose hysterosalpingogram revealed contralateral tubal obstruction.
Hysteroscopic permanent tubal sterilization using a nitinol-dacron intratubal device without anaesthesia in the outpatient setting: procedure feasibility and effectiveness. Litta P1, Hum Reprod. 2005 Dec;20(12):3419-22. Epub 2005 Aug 5. http://www.ncbi.nlm.nih.gov/pubmed/?term=16085664	Italy	Prospective study of 36 consecutive cases of outpatient hysteroscopic tubal sterilization using a nitinol-dacron intratubal device without anaesthesia. Successful device placement was assessed after 3 months by hysterosalpingography and diagnostic hysteroscopy.	Successful bilateral placement was obtained in 32 patients (88.9%); in one (2.8%) the placement was monolateral; and in three (8.3%) the procedure failed.
Confirmation of Essure placement using transvaginal ultrasound. Veersema S1, J Minim Invasive Gynecol. 2011 Mar-Apr;18(2):164-8. doi: 10.1016/j.jmig.2010.10.010. Epub 2011 Feb 1. http://www.ncbi.nlm.nih.gov/pubmed/21282075	Netherlands	Prospective multicenter cohort study Outpatient departments of 4 teaching hospitals in the Netherlands. 1145 who underwent hysteroscopic sterilization using the Essure device between March 2005 and December 2007.	Of 1145 hysteroscopic sterilization procedures, bilateral placement was successful in 1034 (90.3%), unilateral placement was successful in 13 patients (1.1%), and bilateral placement was successful after the second attempt in 25 patients (2.2%). The overall successful placement rate was 93.6% (1072 of 1145 intentions to treat).
3D ultrasound to assess the position of tubal sterilization microinserts. Legendre G1, Hum Reprod. 2011 Oct;26(10):2683-9. doi: 10.1093/humrep/der242. Epub 2011 Aug 16. http://www.ncbi.nlm.nih.gov/pubmed/21846730	France	Retrospective observational study examined the case records of 311 women who underwent hysteroscopic sterilization from October 2002 through October 2008.	The insertion procedure was completed in 94.2% patients.
Learning curve of hysteroscopic placement of tubal sterilization microinserts in 15 gynecologists in the Netherlands.	Netherlands	Prospective multicenter study, Ten community (teaching) hospitals in the Netherlands. A	Excluding the patients lost to follow-up evaluation, a success rate of 80.0% for bilateral placement at first attempt was

Study	Location	Design	Findings
Janse JA1, Fertil Steril. 2013 Sep;100(3):755-60. doi: 10.1016/j.fertnstert.2013.05.019. Epub 2013 Jun 12. http://www.ncbi.nlm.nih.gov/pubmed/?term=23768989		total of 631 women who underwent permanent sterilization by tubal microinserts.	observed. If successful unilateral placement at first attempt is included (ptx w prior tubectomy), a success rate of 82.2% was achieved. Including the successes at second attempt, an overall success rate of 89.5% was achieved
Oral analgesia vs intravenous conscious sedation during Essure Micro-Insert sterilization procedure: randomized, double-blind, controlled trial. Thiel JA1J Minim Invasive Gynecol. 2011 Jan-Feb;18(1):108-11. doi:10.1016/j.jmig.2010.10.006 http://www.ncbi.nlm.nih.gov/pubmed/21195962	Canada	Randomized, double-blind, placebo-controlled trial. Tertiary care ambulatory women's clinic. 87 women of reproductive age women requesting permanent sterilization.	Bilateral placement of the Essure device was successful in 84 patients (97%).
Incorrect position of Essure microinserts 3 months after successful bilateral placement. Gerritse MB, Fertil Steril. 2009 Mar;91(3):930.e1-5. doi: 10.1016/j.fertnstert.2008.08.013. Epub 2008 Oct 21. http://www.ncbi.nlm.nih.gov/pubmed/?term=18945426	Netherlands	100 patients who underwent hysteroscopic sterilization using Essure between December 2003 and June 2004.	Bilateral placement of Essure microinserts in one session was successful in 93 women (93%).
Hysteroscopic tubal sterilization with Essure intratubal devices: a case-control prospective with inert local anesthesia or without anesthesia P. Lopes, E. Gibon, T. Linet, et al. Eur J Obstet Gynecol Reprod Biol, 138 (2008), pp. 199–203 http://www.ncbi.nlm.nih.gov/pubmed/?term=17822834	France	Between February 2002 and May 2005, one operator performed 140 sterilization procedures in this prospective study: the first 70 were performed using local anesthesia and the following 70 began without administration of anesthesia.	N=140, successful bilateral micro-insert placement in the first 70 cases,utilizing paracervical block, was 82.8% and did not differ significantly from the next 70 cases, without anesthesia (91.4%). The first 70 patients had hysteroscopic sterilization under anesthesia: the success rate for placement was 58/70 (82.8%). Hysteroscopies for the next 70 patients, planned without anesthesia, had a placement success rate of 91.4% (64/70). Overall: 87.1% (122/140).
Efficacy of Essure hysteroscopic sterilization--5 years follow up of 1200 women Ríos-Castillo JE Gynecol Endocrinol. 2013 Jun;29(6):580-2. doi: 10.3109/09513590.2013.777419. Epub 2013 Apr 5 http://www.ncbi.nlm.nih.gov/pubmed/?term=23557170	Spain	Outpatient hysteroscopy facility in a large teaching hospital. One-thousand three-hundred and twenty-one women who underwent hysteroscopic sterilization with Essure device, between January 2003 and May 2005.	Satisfactory insertion was accomplished in the first attempt in 1166 women (97.16%). After the second attempt, successful insertion rate rise to 98.6%. (n = 31, 2.6%).

10.3 Appendix 3: Literature on pregnancy

Studies reporting the efficacy of Essure

Study	Location	Design	Findings
<p>Complications of hysteroscopic Essure(®) sterilisation: report on 4306 procedures performed in a single centre. Povedano B, Arjona JE, Velasco E, Monserrat JA, Lorente J, Castelo-Branco C. BJOG. 2012 Jun;119(7):795-9. http://www.ncbi.nlm.nih.gov/pubmed/22360159</p>	Spain	Retrospective 7-year study, Office hysteroscopic unit in a teaching hospital, 4306 women sterilized from March 2003 and June 2010.	A total of seven women (0.16%), three before and four after the 3-month follow-up, became pregnant. The pregnancy rate after assessment for satisfactory placement was 0.09%.
<p>Success rate and patient satisfaction with the Essure sterilisation in an outpatient setting: a prospective study of 857 women. Miño M, Arjona JE, Cordon J, Pelegrin B, Povedano B, Chacon E. BJOG. 2007 Jun;114(6):763-6. http://www.ncbi.nlm.nih.gov/pubmed/?term=17516970</p>	Spain	Single center, retrospective, 857 women with Essure system inserted in an outpatient setting, January 2003 to January 2005.	1 luteal phase pregnancy (elective termination) [1/857 = 0.1%].
<p>Prospective analysis of office-based hysteroscopic sterilization. J Minim Invasive Gynecol. 2006;13(2):98–101. Levie MD, Chudnoff SG. http://www.ncbi.nlm.nih.gov/pubmed/16527710</p>	USA, NY	Prospective, longitudinal analysis, university out-patient office.	Bilateral placement of the device was successful in 98 (96%) of 102 patients; 92 had f/u HSG, 1 pregnancy (luteal phase pregnancy-pregnant at procedure or shortly after – SAB at 8w).
<p>Hysteroscopic sterilization in a large group practice: experience and effectiveness. Savage UK, Masters SJ, Smid MC, Hung YY, Jacobson GF Obstet Gynecol. 2009 Dec; 114(6):1227-31. http://www.ncbi.nlm.nih.gov/pubmed/?term=19935023</p>	Kaiser Permanente, Northern California	Retrospective review, 884 women, 118 physicians at 30 different facilities, from January 2004 to December 2006, f/u through December 31, 2008.	8 pregnancies identified: one never returned for hysterosalpingography, four had a hysterosalpingogram with at least one patent tube, and three had the hysterosalpingogram interpreted as bilaterally occluded.
<p>What are the factors predictive of hysterosalpingogram compliance after female sterilization by the Essure procedure in a publicly insured population? Howard DL, Wall J, Strickland JL Matern Child Health J. 2013 Dec;17(10):1760-7. http://www.ncbi.nlm.nih.gov/pubmed/23212398</p>	2 campuses of Truman Medical Center, Kansas City, Missouri	Multi-site, 136 patients presenting for Essure placement between January 1 2005 and December 31 2010.	No documented pregnancies after the Essure procedure.

Study	Location	Design	Findings
<p>Hysteroscopic sterilization success in outpatient vs office setting is not affected by patient or procedural characteristics. Anderson TL, Yunker AC, Scheib SA, Callahan TL. J Minim Invasive Gynecol. 2013 Nov-Dec;20(6):858-63 http://www.ncbi.nlm.nih.gov/pubmed/23891205</p>	Vanderbilt (Nashville, TN)	Retrospective cohort analysis 638 women who underwent hysteroscopic sterilization between July 1, 2005 and June 30, 2011	3 Pregnancies were reported prior to HSG confirmation (ITT=3/638=0.5%); these were unrelated to device effectiveness. No pregnancies were reported in patients after HSG confirmation of bilateral tubal occlusion, with follow-up of at least 1 year. In addition, there were no known pregnancies among patients who did not complete HSG assessment.
<p>Placement of a permanent birth control device at a university medical center. J Reprod Med. 2009 Apr;54(4):218-22. Shavell VI¹, Abdallah ME, Diamond MP, Berman JM. http://www.ncbi.nlm.nih.gov/pubmed/19438163</p>	Detroit Medical Center (Detroit, MI)	Single center, retrospective study, 316 women, January 1, 2003 to June 30, 2007	3 pregnancies, giving a post-Essure pregnancy rate of 0.95%
<p>Hysteroscopic female sterilization with Essure in an outpatient setting. Anderson S, Eriksson S, Mints M. Acta Obstet Gynecol Scand. 2009;88(6):743-6 http://www.ncbi.nlm.nih.gov/pubmed/19412802</p>	Sweden	61 women underwent outpatient hysteroscopic sterilization between 2002 and 2007	No pregnancies reported
<p>Initial Asian experience in hysteroscopic sterilisation using the Essure permanent birth control device. BJOG. 2005 Sep;112(9):1322-7. Chern B, Siow A. http://www.ncbi.nlm.nih.gov/pubmed/16101615</p>	Singapore	A retrospective study. 80 women seeking permanent birth control. Analysis of initial experience with the Essure permanent birth control over a period of 18 months from 22 June 2001 to 22 December 2002. [<i>probably ESS205</i>]	As of 22 June 2003, 77 women had accumulated 1218 months of wearing the Essure device. All 67 women on regular follow up and with radiological evidence of tubal occlusion did not report any pregnancy. The remaining 10 women who have defaulted follow up were contacted through telephone and there were no pregnancies reported as well in this group.
<p>Pregnancies after hysteroscopic sterilization: a systematic review. Contraception. 2013 May;87(5):539-48. doi: 10.1016/j.contraception.2012.08.006. Epub 2012 Oct 4. Cleary TP, Tepper NK, Cwiak C, Whiteman MK, Jamieson DJ, Marchbanks PA, Curtis KM. http://www.ncbi.nlm.nih.gov/pubmed/23040124</p>	Review	22 studies of women who underwent Essure [®] placement	Eleven articles that documented bilateral tubal occlusion with hysterosalpingogram (HSG) or placement with X-ray or ultrasound following Essure [®] placement did not report any pregnancies with follow-up ranging from 7 months to 7 years. The remaining 11 articles identified 102 reported pregnancies. Eighteen of these pregnancies occurred prior to the 3-month period required before imaging for contraceptive reliability. Two articles did not report what follow-up imaging was performed among women after Essure [®] placement; one of

Study	Location	Design	Findings
			these articles reported three pregnancies.
Female sterilisation: a cohort controlled comparative study of ESSURE versus laparoscopic sterilisation. Duffy S, Marsh F, Rogerson L, Hudson H, Cooper K, Jack S, Hunter D, Philips G. BJOG. 2005 Nov;112(11):1522-8. http://www.ncbi.nlm.nih.gov/pubmed/16225573	UK	Cohort controlled comparative study. The day surgery and outpatient unit of three large UK hospitals. 59 women underwent Essure placement	The overall bilateral placement rate for the ESSURE device was 81% (48/59). One ESSURE patient became pregnant during the study (1/48 = 2.1%). (Author comments "X-ray showed that the left device was suspiciously sited. Patient failed to attend on two occasions for HSG. It is likely that she had an immediate expulsion post-placement of the left-sided device").
Satisfaction and tolerance with office hysteroscopic tubal sterilization. Arjona JE, Fertil Steril. 2008 Oct;90(4):1182-6. doi: 10.1016/j.fertnstert.2007.08.007. Epub 2008 Jan 16. http://www.ncbi.nlm.nih.gov/pubmed/?term=18201703	Spain	Prospective analysis of case series 1615 patients Up to 3.5years follow-up/analysis	3 pregnancies (0.2%) [all occurred within 90d]
Unintended pregnancies after Essure sterilization in the Netherlands. Veersema S, Fertil Steril. 2010 Jan;93(1):35-8. doi: 10.1016/j.fertnstert.2008.10.005. Epub 2008 Nov 21. http://www.ncbi.nlm.nih.gov/pubmed/19027109	Netherlands	Retrospective, national, multicenter ~ 6000 sterilizations (estimated from distributor) Up to 5 ~ follow-up	Ten cases of unintended pregnancies after Essure sterilization out of an estimated 6000 hysteroscopic sterilizations were performed in 45 hospitals in the Netherlands from August 2002 to May 2008 as estimated from the data from the Dutch distributor (<i>Author comments "Most pregnancies occurred in patients with only one device placement and bilateral occlusion on HSG. Other cases included misinterpretation of HSG, undetected abnormal device position by ultrasound, one undetected preprocedure pregnancy, and two patient failures to follow up with the physician advice"</i>).
Essure® permanent birth control effectiveness: a seven-year survey. Jost S1, Eur J Obstet Gynecol Reprod Biol. 2013 Jun;168(2):134-7. doi: 10.1016/j.ejogrb.2012.12.042. Epub 2013 Jan 30. http://www.ncbi.nlm.nih.gov/pubmed/?term=23375900	France	Retrospective National Survey 53,003 kits sold. Up to 7.5 follow-up/analysis	58 pregnancies (0.1%).
Contraceptive Failures Associated with Hysteroscopic Sterilization. Deraleu & Heinlein, Journal of Minimally Invasive Gynecology S19 (2012) S1–S35 http://www.jmig.org/article/S1553-4650(12)00398-	USA, CA	Retrospective study, Kaiser Permanente Northern California database 2,621 women, 1-8 year follow-up/analysis	22 pregnancies (0.8%)

Study	Location	Design	Findings
6/abstract Tubal sterilization: pregnancy rates after hysteroscopic versus laparoscopic sterilization in France, 2006-2010. Fernandez H, Eur J Obstet Gynecol Reprod Biol. 2014 Sep;180:133-7. doi: 10.1016/j.ejogrb.2014.04.043. Epub 2014 May 14. http://www.ncbi.nlm.nih.gov/pubmed/?term=24993770	France	Retrospective cohort study To compare the rates of pregnancy among women who underwent Essure hysteroscopic sterilization versus tubal ligation in France between 2006 and 2010	During the study period, French hospitals performed 109,277 tubal sterilization procedures: 39,169 Essure sterilizations and 70,108 laparoscopic tubal ligations A Cox model has been performed. Following sterilization, after adjustment on age Essure patients became pregnant at a significantly lower rate than laparoscopic ligation patients 0.36% versus 0.46%, respectively (HR=0.62 (040-096).
Efficacy of Essure hysteroscopic sterilization--5 years follow up of 1200 women Rios-Castillo JE Gynecol Endocrinol. 2013 Jun;29(6):580-2. doi: 10.3109/09513590.2013.777419. Epub 2013 Apr 5 http://www.ncbi.nlm.nih.gov/pubmed/?term=rios-castillo+essure	Spain	Retrospective analysis of case series	Three pregnancies had been reported after 5 years follow up, which implies an overall absolute rate of 0.25%. This represents a Pearl index of 0.05 after 72,000 months of surveillance. All of them occurred in the first year of use of the microinsert. There has been no unintended pregnancy in the next 4 years
Definitive contraception with Essure device: Single institutional experience on 517 procedures Aparicio-Rodríguez-Miñón P Ginecol Obstet Mex. 2015 Jan;83(1):16-22 http://www.ncbi.nlm.nih.gov/pubmed/26016312	Spain	Retrospective descriptive study, 517 procedures	7 (1.35%) unintended pregnancies were observed.
Hysteroscopic tubal sterilization with Essure intratubal devices: a case-control prospective with inert local anesthesia or without anesthesia P. Lopes, E. Gibon, T. Linet, et al. Eur J Obstet Gynecol Reprod Biol, 138 (2008), pp. 199–203 http://www.scopus.com/record/display.url?eid=2-s2.0-44249108785&origin=inward&txGid=A2B0BAFB03F13FF7328FB6DC3CE46984.I0QkqbljGqqLQ4Nw7dqZ4A%3a8	France	Prospective case-controlled cohort study, 140 women	No pregnancy was reported in our study. The meanfollow-up for the first group of patients (with local anesthesia) was 48 months (39–58 months). The mean follow-up for the second group was shorter: 30 months (19–39 months) because it was a before-and-after study.
Essure: a new device for hysteroscopic tubal sterilization in an outpatient setting A. Ubeda, R. Labastida, S. Dexeus Fertility & Sterility, Volume 82, Issue 1, July 2004, Pages 196–199 http://www.sciencedirect.com/science/article/pii/S0015028204006843	Spain	Prospective, observational study, Private university hospital, 85 women, enrolled July 2002 - July 2003	None of the patients has become pregnant after device placement.

Study	Location	Design	Findings
3D ultrasound to assess the position of tubal sterilization microinserts. Legendre G1, Hum Reprod. 2011 Oct;26(10):2683-9. doi: 10.1093/humrep/der242. Epub 2011 Aug 16. http://www.ncbi.nlm.nih.gov/pubmed/21846730	France	Retrospective observational study examined the case records of 311 women who underwent hysteroscopic sterilization from October 2002 through October 2008. Imaging with 3D-US or pelvic X-radiography or both was performed 3 months after the procedure to verify device position. Hysterosalpingography (HSG) was performed when a bilateral procedure was not completed because of a history of salpingectomy or blocked tube, when doubt persisted after 3D-US or pelvic radiography, or for comparative purposes in a prospective study. The positions seen on 3D-US were classified in four categories according to a specific scale we devised.	2 pregnancies (live birth, elective termination). Neither pregnancy nor early expulsion occurred when 3D-US found that the devices were correctly placed.
Confirmation of Essure placement using transvaginal ultrasound. Veersema S1, J Minim Invasive Gynecol. 2011 Mar-Apr;18(2):164-8. doi: 10.1016/j.jmig.2010.10.010. Epub 2011 Feb 1. http://www.ncbi.nlm.nih.gov/pubmed/21282075	Netherlands	Prospective multicenter cohort study. Outpatient departments of 4 teaching hospitals in the Netherlands. 1145 women who underwent hysteroscopic sterilization using the Essure device between March 2005 and December 2007	The cumulative pregnancy rate after 18 months was 3.85 per thousand women.

10.4 Appendix 4: Literature on post-tubal syndrome

Studies cited in ACOG “*Practice Bulletin on the Benefits and Risks of Sterilization*” regarding menstrual patterns/changes during various procedures of sterilization:

- Bhiwandiwalla PP et al. (Level II-3): Menstrual pattern changes following laparoscopic sterilization with different occlusion techniques: a review of 10,004 cases. *Am J Obstet Gynecol* 1983; 145:684-94. <http://www.ncbi.nlm.nih.gov/pubmed/6219585>.
- DeStefano F et al. (Level II-2): Long-term risk of menstrual disturbances after tubal sterilization. *Am J Obstet Gynecol* 1985; 152:835-41. <http://www.ncbi.nlm.nih.gov/pubmed/4040707>.
- Foulkes J, Chamberlain G. (Level II-3): Effects of sterilization on menstruation. *South Med J* 1985;78:544-7. <http://www.ncbi.nlm.nih.gov/pubmed/3992301>.
- Rivera R et al. (Level I): Menstrual patterns and progesterone circulating levels following different procedures of tubal occlusion. *Contraception* 1989;40:157-69. <http://www.ncbi.nlm.nih.gov/pubmed/2758840>.
- Sahwi S et al. (level I): Changes in menstrual blood loss after four methods of female tubal sterilization. *Contraception* 1989;40:387-98. <http://www.ncbi.nlm.nih.gov/pubmed/2582767>.
- Thranov I et al. (Level I): Hormonal and menstrual changes after laparoscopic sterilization by Falope-rings or Filshie-clips. *Fertil Steril* 1992;57:751-5. <http://www.ncbi.nlm.nih.gov/pubmed/1532561>.
- Rulin MC et al. (Level II-2): Long Term effect of tubal sterilization on menstrual indices and pelvic pain. *Obstet Gynecol* 1993;82:118-21. <http://www.ncbi.nlm.nih.gov/pubmed/8515910>.

10.5 Appendix 5: Menstrual changes post-Essure placement

Studies with data on menstrual changes post-Essure placement

Study	Location / Design	Findings
Kerin JF, Munday D, Ritossa M, Rosen D. Tissue encapsulation of the proximal Essure micro-insert from the uterine cavity following hysteroscopic sterilization. <i>J Minimally Inv Gynecol.</i> Mar-Apr 2007;14(2):202-204. http://www.sciencedirect.com/science/article/pii/S1553465006005474#	US, multicenter / Retrospective observational study of 545 women who underwent Essure placement, 1997-2005	Of the 545 placements, 20 (3.7%) required a subsequent diagnostic second look hysteroscopy for abnormal and/or heavy bleeding. Mean interval to second look, 19.73 months (SD 11.74). Findings included 2 endocervical polyps, 3 endometrial polyps, 1 submucosal myoma, 1 case benign hyperplasia. No evidence of visual epithelial abnormalities, tissue inflammation or abnormal vascular patterns.
Sinha D, Kalathy V, Gupta JK, Clark TJ. The feasibility, success and patient satisfaction associated with outpatient hysteroscopic sterilization. <i>BJOG</i> 2007 Jun;114(6):676-683. http://www.ncbi.nlm.nih.gov/pubmed/17516957	UK / Prospective cohort study, 112 consecutive Essure placements	Successful placement 103/112 (92%); HSG in 84/112 (75%); 100% tests satisfactory. 76/84 (90%) returned questionnaire. At 3 months: 20/76 (26%) reported persistent change in menstrual pattern, 14/76 (18%) reported heavier, 2/76 (3%) reported irregular, 1/76 (1%) reported intermenstrual spotting, 3/76 reported amenorrhea.
Miño M, Arjona JE, Cordon J, Pelegrin B, Povedano B, Chacon E. Success rate and patient satisfaction with the Essure sterilisation in an outpatient setting: a prospective study of 857 women. <i>BJOG.</i> 2007 Jun;114(6):763-6. http://www.ncbi.nlm.nih.gov/pubmed/?term=17516970 .	Spain / Single center, prospective study of 857 women undergoing Essure	Successful placement 857/863(99%); 3 month follow up 100%. On survey, no woman reported any change to volume or pattern of menstruation.
Chern B and Siow A. Initial Asian experience in hysteroscopic sterilization using the Essure permanent birth control device. <i>BJOG: an International J Obstet Gynecol.</i> 2005;Sep(112):1322-1327. http://onlinelibrary.wiley.com/doi/10.1111/j.1471-0528.2005.00436.x/full).	Singapore / Retrospective study, June 2001-December 2002, 80 women undergoing Essure placement	Successful placement: 77/80 (96.2%); 3 month confirmation test 67/77 (87%), 100% satisfactory; at 3 month follow up, no patient complained of menstrual disturbance.
Ubeda A, Labastida R, Dexeus S. Essure: a new device for hysteroscopic tubal sterilization in an outpatient setting. <i>Fertility & Sterility</i> , Volume 82, Issue 1, July 2004, Pages 196–199. http://www.ncbi.nlm.nih.gov/pubmed/15237011 .	Spain / Prospective, observational study 85 women, July 2002 – July 2003	Successful placement: 81/85 (95%); 3 month confirmation test 75/81 (93%), 100% satisfactory. No patients reported change to menstrual pattern.
Anderson S, Eriksson S, Mints M. <i>Acta Obstet Gynecol Scand.</i> 2009;88(6):743-6. Hysteroscopic female sterilization with Essure in an outpatient setting. http://www.ncbi.nlm.nih.gov/pubmed/19412802 .	Sweden / Retrospective review of 61 women who had Essure placed 2002-2007. Questionnaires sent 1-5 years post-procedure	Successful placement: 58/61 (95%); 3 month confirmation test 100%; 57/58 (98%) confirmation tests satisfactory. Questionnaire sent to all women; 50/61 (82%) returned. 9/50 (18%) reported heavier periods (none sought treatment) and 8/50 (16%) reported lighter periods.

US = United States; UK = United Kingdom; SD = standard deviation; HSG = hysterosalpingogram.

10.6 Appendix 6: Instructions for Use

ESSURE[®]
permanent birth control
INSTRUCTIONS FOR USE

essure[®]
permanent birth control

Instructions for use	English	USA only ⚠	1 - 6
Instrucciones de uso	Español		7 - 10
Mode d'emploi	Français		11 - 14
Gebruiksaanwijzing	Nederlands		15 - 18
Gebrauchsanweisung	Deutsch		19 - 22
Brugsanvisning	Dansk		23 - 26
Instructions for use	English	Outside the USA only ⚠	27 - 30
Istruzioni per l'uso	Italiano		31 - 34
Bruksanvisning	Norsk		35 - 38
Instruções de utilização	Português		39 - 42
Bruksanvisning	Svenska		43 - 46
Käyttöohjeet	Suomi		47 - 50
Οδηγίες Χρήσης	Ελληνικά		51 - 54

INSTRUCTIONS FOR USE

IMPORTANT

- Caution: Federal law restricts this device to sale by or on the order of a physician. Device to be used only by physicians who are knowledgeable hysteroscopists; have read and understood the Instructions for Use and Physician Training Manual; and have successfully completed the **Essure**[®] training program, including preceptoring in placement until competency is established, typically 5 cases.

IMPORTANT

- Patient should not rely on **Essure** inserts for contraception until an **Essure** Confirmation Test [modified hysterosalpingogram (HSG)] demonstrates bilateral tubal occlusion **and** satisfactory location of inserts.
- If **Essure** inserts cannot be placed bilaterally, then the patient should not rely on **Essure** inserts for contraception. The effectiveness of the **Essure** procedure when unilateral insert placement occurs has not been evaluated.
- This product does not protect against HIV infection or other sexually transmitted infections.

I. OVERVIEW OF ESSURE PROCEDURE AND PRINCIPLES OF OPERATION

Step 1: **Essure** insert placement procedure ("procedure").

Step 2: Patient must remain on alternative contraception until **Essure** Confirmation Test (modified HSG).

Step 3: **Essure** Confirmation Test (modified HSG) must show fallopian tube satisfactory insert location and occlusion before the patient can rely on **Essure** for contraception.

Using a transvaginal approach, one flexible **Essure** insert ("insert") is placed in the proximal portion of each fallopian tube lumen. The insert expands upon release to conform to and acutely anchor in the tubal lumen. Subsequently, the insert elicits a benign tissue in-growth that permanently occludes the fallopian tube, resulting in contraception.

II. DEVICE DESCRIPTION

The **Essure** system ("system") is comprised of a disposable delivery system and a wound-down insert. A disposable introducer is also provided to facilitate delivery system entry into the hysteroscope operating channel. The insert consists of a super-elastic Nitinol outer coil and a stainless steel inner coil wrapped in polyethylene terephthalate (PET) fibers. The wound-down insert is approximately 4 cm in length and 0.8 mm in diameter (**Figure 1a**). When released, the outer coil expands up to 2.0 mm in diameter, conforming itself to the varied diameters and shapes of the fallopian tube (**Figure 1b**).

Figure 1a
Essure Insert
Wound-down configuration, attached to the delivery system
(NOT TO SCALE)


Figure 1b
Essure Insert
Expanded configuration, detached from the delivery system
(NOT TO SCALE)

The disposable delivery system, (shown in **Figure 2**), consists of a delivery wire, a release catheter, a delivery catheter and a delivery handle.

NOTE: The delivery wire and the release catheter are not visible in **Figure 2**.

Figure 2
Essure System
Showing detail of placement procedure symbols.
(NOT TO SCALE)

The wound-down insert is attached to a Nitinol delivery wire and sheathed by a flexible delivery catheter. A black positioning marker on the delivery catheter aids in proper insert placement.

The delivery handle controls the device delivery and release mechanism. The thumbwheel on the delivery handle retracts the delivery catheter. The button allows the physician to change the function of the thumbwheel from retracting the delivery catheter to deploying the outer coils. The delivery wire is detached from the insert by continuing to rotate the thumbwheel. To remind the physician of these placement procedure steps, symbols are located on the delivery handle (refer to **Figure 2**).

The **DryFlow**[®] introducer ("introducer") (**Figure 3**) helps facilitate entry and advancement of the insert during insertion into the hysteroscope, also minimizing fluid back splash. (Additional introducers are available.)

Figure 3
DryFlow Introducer
(NOT TO SCALE)


III. MECHANISM OF ACTION

A. Placement at Utero-Tubal Junction

Proper placement occurs when the insert spans the utero-tubal junction (UTJ) (**Figure 5**). This location is distal enough to avoid expulsion, yet proximal enough to visualize trailing coils to show placement.

Figure 4:
Ideal Essure Insert Placement

B. Dynamic Anchoring

The insert is a dynamic and flexible spring-like device. The outer coil expands upon deployment, conforms to and pushes against the fallopian tube wall, acutely anchoring the insert in the UTJ.

C. Tubal Occlusion and Tissue In-Growth

Tubal occlusion is attributed to the space filling design of the device and the benign occlusive tissue response. PET fiber causes tissue in-growth into and around the insert, facilitating insert retention and pregnancy prevention.

IV. INDICATIONS FOR USE

Essure is indicated for women who desire permanent birth control (female sterilization) by bilateral occlusion of the fallopian tubes.

V. CONTRAINDICATIONS

Contraindicated for patients who:

- Are uncertain about ending fertility.
- Can have only one insert placed (including contralateral proximal tubal occlusion or suspected unicornuate uterus).
- Have previously undergone a tubal ligation.
- Are pregnant or suspect pregnancy.
- Have delivered or terminated a pregnancy less than 6 weeks prior to the **Essure** procedure.
- Have an active or recent upper or lower pelvic infection.
- Have a known allergy to contrast media.

VI. WARNINGS

- Pregnancies (including ectopic pregnancies) have been reported among women with inserts in place. Some of these pregnancies were due to patient non-compliance, which included failure to:
 - use alternate contraception during the 3-month "waiting period" prior to **Essure** Confirmation Test (modified HSG);
 - return for the **Essure** Confirmation Test (modified HSG) to determine if the inserts are in the correct location and tubal occlusion is present; and
 - use alternate contraception or undergo sterilization by another method if the **Essure** Confirmation Test (modified HSG) reveals tubal patency. In this case, the clinician should inform the patient of the **Essure** Confirmation Test (modified HSG) finding and counsel her not to rely on the **Essure** System for contraception.

Therefore, it is critical that clinicians properly counsel patients regarding the risk of pregnancy (including ectopic pregnancy) attributable to non-compliance during all stages of the **Essure** procedure.

- Physicians performing the **Essure** procedure must adhere to the **Essure** Confirmation Test (modified HSG) protocol in these Instructions for Use. The protocol for interpretation of the **Essure** Confirmation Test (modified HSG) is different from a standard HSG for infertility. In addition to patient non-compliance, incorrect interpretation of the **Essure** Confirmation Test (modified HSG) has led to pregnancy.
- The **Essure** procedure should be considered irreversible. Safety or effectiveness of reversal surgery is unknown.
- Effectiveness rates for the **Essure** procedure are based on patients who had bilateral placement. Little effectiveness data exist for unilateral insert placement in patients with unicornuate uteri or contralateral proximal tubal occlusion (PTO).
- Do not attempt hysteroscopic **Essure** insert removal once placed unless 18 or more trailing coils are seen inside the uterine cavity. Attempted removal with less than 18 trailing coils may result in fractured insert, fallopian tube perforation or other injury.
- To reduce risk of hypervolemia, terminate procedure if distension fluid deficit exceeds 1500cc or hysteroscopic time exceeds 20 minutes. Excess fluid deficit may signal uterine or tubal perforation. If noted, discontinue procedure and evaluate patient for possible perforation.
- The **Essure** insert may conduct energy and cause patient injury if contacted by an active electro-surgical device. Avoid electro-surgery on uterine cornua and proximal fallopian tubes without visualizing inserts. During Laparoscopically Assisted Vaginal Hysterectomy, do not place instruments more proximal than the ampullary portion of the tube.
- Avoid direct contact between inserts and monopolar radio frequency (RF) energy instruments as bench studies suggest damage to surrounding tissue should monopolar RF instruments contact inserts.
- Patients undergoing immunosuppressive therapy (e.g. systemic corticosteroids or chemotherapy) are discouraged from undergoing the **Essure** procedure because the immunosuppressant may lead to decreased tissue in-growth.
- To reduce risk of uterine perforation, terminate procedure if excessive force is required to achieve cervical dilation.
- Never attempt to advance **Essure** insert(s) against excessive resistance. If tubal or uterine perforation occurs or is suspected, discontinue procedure and work-up patient for possible complications related to perforation, including hypervolemia. A false positive HSG and pregnancy have been associated with tubal perforation by insert in the literature; evaluate **Essure** Confirmation Test for perforation if excessive resistance is experienced during procedure. Only 1.8% (12/682) of clinical trial patients had device related perforations. If necessary, retrieval of perforating inserts requires surgery.
- Some **Essure** patients have reported pelvic pain that may be device related. If device removal is indicated, this will require surgery. (refer to section XVIII. INSERT REMOVAL)
- Effects, including risks, of **Essure** inserts on in vitro fertilization (IVF) have not been evaluated. Risks to the patient, fetus and continuation of pregnancy are unknown.
- The **Essure** micro-insert includes nickel-titanium alloy, which is generally considered safe. However, *in vitro* testing has demonstrated that nickel is released from this device. Patients who are allergic to nickel may have an allergic reaction to this device, especially those with a history of metal allergies. In addition, some patients may develop an allergy to nickel if this device is implanted. Typical allergy symptoms reported for this device include rash, pruritus, and hives.

Essure and Global Endometrial Ablation

- DO NOT perform the **Essure** procedure concomitantly with endometrial ablation. Ablation causes intrauterine synechiae which can compromise (i.e. prevent) the **Essure** Confirmation Test (modified HSG). Women with inadequate confirmation tests cannot rely on **Essure** for contraception. Bench and clinical studies demonstrate balloon thermal (THERMACHOICE® Uterine Balloon System) and hydro-thermal (HTA® System) endometrial ablation can be safely and effectively performed with **Essure** inserts in place. However, balloon thermal and hydro-thermal endometrial ablation should only be performed after the **Essure** Confirmation Test (modified HSG). Bench and clinical studies have been conducted which demonstrate that bipolar radio frequency (RF) NovaSure® Impedance Controlled Endometrial Ablation System endometrial ablation of the uterus can be safely performed with **Essure** inserts in place. However, thermal injury to the proximal portion of the fibrotic in-growth that causes tubal occlusion may occur. It is unknown whether partial thermal injury will interfere with tubal occlusion. Contraception rates following NovaSure with **Essure** inserts in place are under investigation. Bipolar RF ablation device may contact inserts. As a result of contact, heat from a bipolar RF device may be propagated along the insert. This could cause bowel or bladder injury if there is an unrecognized tubal perforation and part of the insert lies outside of the tubal serosa. Therefore, do not perform bipolar RF ablation in patients who may have undiagnosed perforation e.g., patients whose **Essure** procedure was difficult or atypical. Performing intrauterine procedures without direct visualization may result in trailing coils of insert being ensnared in another device. When device is withdrawn, the insert may be removed and tubal patency be restored. In one study, the risk of insert being ensnared by a bipolar RF endometrial ablation array was approximately 3%. Note: All trademarks are property of their respective companies.
- Safety of cryo-ablation, laser-ablation, or microwave ablation with **Essure** inserts in place is unknown and little data exist. Microwave energy near metallic implants may pose risk of serious patient injury; therefore avoid use of microwave endometrial ablation devices near inserts.

VII. PRECAUTIONS

- Women undergoing sterilization at a younger age are at greater risk of regretting their decision.
- Perform the **Essure** procedure during early proliferative phase of the menstrual cycle to increase ostia visualization and prevent placement in a patient with an undiagnosed (luteal phase) pregnancy. Women with menstrual cycles shorter than 28 days should undergo careful ovulation day calculations.
- As with any tubal sterilization procedure, performing endometrial ablation following **Essure** insert placement may increase the risk of post-ablation tubal sterilization syndrome.
- Safety and effectiveness of **Essure** is not established in patients under 21 or over 45 years old, nor in patients who delivered or terminated a pregnancy less than 8-12 weeks before procedure.
- Do not advance the **Essure** delivery system once the end of the black positioning marker has reached the tubal ostium as unsatisfactory insert placement could result.
- Use caution and avoid the **Essure** inserts when undertaking blind intrauterine procedures as disturbing the inserts could interrupt their ability to prevent pregnancy. Direct visualization of inserts during intrauterine procedures is optimal. Insert retention and location should be verified by hysteroscopy, X-ray, or ultrasound following intrauterine procedures. There could be risks associated with intrauterine procedures and the presence of inserts not currently identified.
- Have appropriate equipment, medication, staff, and training in place to handle emergency situations such as vaso-vagal response.
- Uterine or fallopian tube anomalies may make it difficult to place **Essure** inserts. Identify and assess both tubal ostia prior to insert placement. Do not attempt to place an insert in one tubal ostium unless there is a reasonable expectation that the contralateral tube is accessible and patent.
- Do not advance the **Essure** system if the patient is experiencing extraordinary pain or discomfort. Terminate the procedure and work-up patient for possible perforation.
- Never attempt to re-sterilize the **Essure** system as it is single use only. Resterilization may adversely affect device function or cause patient injury.
- Utilize eye protection as back splash of hysteroscopic distention fluid (saline) into face can occur.
- Use an introducer to avoid insert tip damage.
- Keep the operating channel of the hysteroscope open to avoid damage to insert or introducer.
- Do not place more than one **Essure** insert in a single fallopian tube.
- Do not use **Essure** if sterile package is open or damaged, or if insert is damaged.
- MRI Information.** The **Essure** insert was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

Non-clinical testing demonstrated that the **Essure** insert is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the **Essure** insert produced the following temperature rise during MRI performed for 15-min in the 3-Tesla (3-Tesla/128-MHz, Excite, Software G3.0-052B, General Electric Healthcare, Milwaukee, WI) MR system: Highest temperature change +1.7°C

Therefore, the MRI-related heating experiments for the **Essure** insert at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 3.0-W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.8-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.7°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the **Essure** insert. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Dimensions: Wound-down and expanded length: 4-cm
Expanded diameter: 1.5 to 2.0-mm

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	173-mm ²	53-mm ²	621-mm ²	277-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

VIII. ADVERSE EVENTS

A. Patient Population

From November 1998 - June 2001, a total of 745 women underwent the **Essure** procedure in two clinical trials that evaluated safety and effectiveness (227 Phase II; 518 Pivotal¹). Placement of at least one insert was achieved in 682 women (206 Phase II study; 476 Pivotal). If bilateral placement was not initially achieved, some women underwent additional procedure(s).

B. Observed Adverse Events

Adverse events resulting from the placement procedure are detailed in Table 1.

Table 1
Adverse events reported day of placement procedure in Phase II & Pivotal Trials

Adverse Event/ Side Effect	Phase II Trial		Pivotal Trial	
	Number (N=233 procedures)	Percent	Number (N=544 procedures)	Percent
Cramping	**	**	161	29.6%
Pain	2	0.9%	70	12.9%
Nausea/vomiting	**	**	59	10.8%
Dizziness/light headed	**	**	48	8.8%
Bleeding/spotting	**	**	37	6.8%
Other	**	**	16*	2.9%
Vaso-vagal response	2	0.9%	7	1.3%
Hypervolemia	**	**	2	0.4%
Band Detachment	3	1.3%	2	0.4%

*Includes: ache (3), hot/hot flashes (2), shakiness (2), uncomfortable (1), weak (1), profuse perspiration (1), bowel pain (1), sleepiness (1), skin itching (1), loss of appetite (1), bloating (1), allergic reaction to saline used for distension (1).

**Data not collected

During and immediately following the procedure, the majority of participants experienced mild to moderate pain. The majority of participants experienced spotting for an average of 3 days after the procedure. Pain was managed with oral non-steroidal anti-inflammatory drugs (NSAIDs) or oral narcotic pain reliever.

Table 2 summarizes adverse events rated as “possibly” related to the insert or procedure during the first year of reliance in the Pivotal trial (approximately 15 months post-device placement). Percentages reflect the number of *events* divided by the number of *participants* in the trial. When numerous episodes of the same event were reported by one participant, each report was counted as a separate event. Therefore, percentages may over-represent the percentage of *women* who have experienced that event.

Table 2
Pivotal Trial
Adverse Events by Body Systems, First Year of Reliance*
(N=476 patients implanted with at least one insert)

Adverse Events by Body System	Number	Percent
Abdominal:		
Abdominal pain/abdominal cramps	18	3.8%
Gas/bloating	6	1.3%
Musculo-skeletal:		
Back pain/low back pain	43	9.0%
Arm/leg pain	4	0.8%
Nervous/Psychiatric:		
Headache	12	2.5%
Premenstrual Syndrome	4	0.8%
Genitourinary:		
Dysmenorrhea/menstrual cramps (severe)	14	2.9%
Pelvic/lower abdominal pain (severe)	12	2.5%
Persistent increase in menstrual flow	9**	1.9%
Vaginal discharge/vaginal infection	7	1.5%
Abnormal bleeding - timing not specified (severe)	9	1.9%
Menorrhagia/prolonged menses (severe)	5	1.1%
Dyspareunia	17	3.6%
Pain/discomfort - uncharacterized:	14	2.9%

* Only events occurring in ≥ 0.5% are reported

** Eight women reported persistent decrease in menstrual flow

In the Phase II trial, 12/206 (5.8%) women with at least one insert reported episodes of period pain, ovulatory pain, or changes in menstrual function.

C. Potential Adverse Events Not Observed in Clinical Studies

The following adverse events were not experienced by Phase II/Pivotal clinical trial participants prior to marketing but are still possible and/or have occurred in the commercial setting:

- Pregnancy and ectopic pregnancy in women relying on **Essure**² inserts.
- Perforation of internal bodily structures other than the uterus and fallopian tube.
- Adnexal infection/salpingitis.
- Adverse events associated with the **Essure** Confirmation Test (modified HSG).
- Pregnancy due to uterine or fallopian tube procedures causing failure of insert.
- Adverse events associated with surgery attempting to reverse the procedure, pregnancy following a reversal, or an IVF procedure.
- Adverse events associated with gynecologic surgical procedures (e.g. endometrial ablation).

D. Adverse Event Reporting

Report any adverse event involving the **Essure** system to **Bayer HealthCare LLC** immediately by calling 888-84 BAYER (888-842-2937).

IX. CLINICAL STUDIES

A. Purpose of the Study, Study Design, Primary Endpoints

Two clinical trials were conducted (Phase II Trial; Pivotal Trial) prior to marketing to demonstrate safety and effectiveness of the **Essure** system in providing permanent contraception. Additionally, two studies were performed after pre-marketing approval to evaluate rates of bilateral **Essure** placement, however, the study using the ESS305 is the only post approval study presented in this label.

1. Phase II Trial

Phase II was a prospective, multi-center, single-arm, non-randomized, international study which evaluated:

- Participant's tolerance of, and recovery from, procedure
- Safety of the procedure
- Participant's tolerance of implanted inserts
- Long-term safety and stability of implanted inserts
- Effectiveness of the inserts in preventing pregnancy

2. Pivotal Trial

The Pivotal trial was a prospective, multi-center, single-arm, non-randomized, international study which used the U.S. Collaborative Review of Sterilization (CREST study) as a qualitative benchmark. The study included the following primary endpoints:

- Prevention of pregnancy
- Safety of insert procedure
- Safety of insert wearing

Secondary endpoints included:

- Participant satisfaction with procedure
- Participant satisfaction with insert wearing
- Bilateral device placement rate
- Profile development for appropriate procedure candidates

¹ Pivotal trial: 657 women initially enrolled; 518 underwent the procedure; 99 changed their minds about participating; 23 did not meet the inclusion criteria and were terminated from study; 17 failed screening tests.

² One Phase II trial participant received a prior device design that was discontinued in 1998 (the Beta design of the STOP device) and became pregnant after nearly two years of reliance. That pregnancy is not included in the effectiveness rate calculations, since that design was not subject of the Premarket Approval Application (PMA) supporting **Essure** approval.

3. ESS305 Post Approval Study

The ESS305 Post Approval Study was a prospective, multi-center, single-arm, non-randomized United States study intended to document the bilateral placement rate using the ESS305 delivery system. The main primary endpoint for the study was successful bilateral placement rate at first attempt.

B. Patients Studied

Table 3

Age Distribution (Combined Data from Pivotal Study and Phase II Study); Average age: 33

<28 years old	28-33 years old	≥34 years old
14%	40%	46%

Table 4
Patient Demographics

	Phase II and Pivotal Studies Combined N=745
Race*	
White/Caucasian	428
Latin	31
Black	24
Other	9
Gravidity	Mean=2.91 (0 – 11)
Parity	Mean=2.23 (0 – 6)
Body Mass Index (BMI) (kg/m²)	Mean=27 (16 – 57)

* Data from Pivotal study only; race not collected in Phase II study.

- The Phase II and Pivotal trials combined consisted of 664 participants in whom bilateral insert placement was achieved after one or more attempts (200 Phase II study; 464 Pivotal). Tables 6 and 7 present patient demographic information. Participants were between 21 and 45 years of age and seeking permanent contraception. All participants had at least one live birth, regular menstrual cycles and were willing to use alternative contraception for the three months following the procedure.
- The study population of the ESS305 Post Approval Study consisted of 581 women in whom insert placement was attempted using the ESS305 delivery system. A total of 70 investigators performed the procedures at 70 US sites. All study participants were between 21 and 51 years of age and were seeking permanent contraception prior to enrollment.

C. Methods

All Phase II and Pivotal Trial participants were screened for eligibility. Medical history, physical examination and required laboratory tests were performed.

Insert placement was attempted in each fallopian tube. In the Pivotal trial, pelvic X-rays were performed within 24 hours of placement to serve as a baseline evaluation of insert location. Participants used alternative contraception for three months following the procedure.

An **Essure** Confirmation Test (modified HSG) was performed three months post procedure to evaluate insert location and fallopian tube occlusion. If bilateral placement and occlusion were satisfactory, participants discontinued alternative contraception and relied on the inserts for contraception.

D. Results

As of the final 5-year follow-up data extracts (Phase II-January 6, 2006; Pivotal-December 5, 2007), 643 trial participants with bilateral placement (194 Phase II; 449 Pivotal) contributed 35,633 months of follow-up time with zero pregnancies reported. Adverse events reported in Pivotal and Phase II clinical studies conducted prior to marketing are provided in Section VIII B above, and events by study are below. Tables 5A, 5B and 6 present principal safety and effectiveness results.

Table 5A
Insert Reliance Rates
in the Phase II and Pivotal Clinical Studies

Outcome	Phase II and Pivotal Trials N=745	
	Number	Percent
Reliance Rate*: <i>Among women with bilateral placement</i>	643/664	97%

*The reliance rate is the number of women who relied on **Essure** for contraception divided by the number of women with bilateral insert placement. In the Phase II trial, the following adverse events prevented reliance: Perforation (7/206; 3.4%, including one patient that that relied for 31 months before laparotomy and cornual resection due to pain, the other six never relied); Expulsion (1/206; 0.5%); Unsatisfactory insert location (1/206; 0.5%); Initial tubal patency (7/200; 3.5%) was found at the 3-month **Essure** Confirmation Test (modified HSG), however, all had tubal occlusion at a 6-month repeat **Essure** Confirmation Test (modified HSG). In the Pivotal trial, the following adverse events prevented reliance: Perforation (5/476; 1.1%); Expulsion (14/476; 2.9%, nine out of the fourteen underwent a successful second placement procedure; Unsatisfactory insert location (3/476; 0.6%); Initial tubal patency (16/456; 3.5%) was found at the 3-month **Essure** Confirmation Test (modified HSG), however, all had tubal occlusion at a 6 or 7-month repeat **Essure** Confirmation Test (modified HSG).

Table 5B
Insert Placement Rate at first attempt in the Commercial Setting
Using the Essure System ESS305

Placement Status	ESS305 Post Approval Study	
	Number	Percent
Bilateral Placement*:	593/612	96.9%
Intent-to-Treat Bilateral Placement**:	593/619	95.8%
Unilateral Placement or No devices placed***:	26/619	4.2%

*Bilateral placement rate excludes seven patients in whom insert placement was not attempted.

** Intent-to-treat bilateral placement rate includes all participants who underwent hysteroscopy regardless of whether insert placement was attempted.

*** Unilateral placement occurred in 10 participants; 6 participants had no placement in either fallopian tube; insert placement was not attempted in 7 participants. Tubal obstruction and visualization problems due to uterine anatomy were the most common reasons for placement failure.

Table 6
Effectiveness Results as of December 2007

CUMULATIVE FAILURE RATES Phase II and Pivotal Trials Combined				
One-Year ^C	Two-Year ^C	Three-Year ^C	Four-Year ^C	Five-Year ^C
0% N=635 (95% CI 0 – 0.10%) ^{A, B}	0% N=605 (95% CI 0 – 0.20%) ^{A, B}	0% N=586 (95% CI 0 – 0.30%) ^{A, B}	0% N=567 (95% CI 0 – 0.40%) ^{A, B}	0% N=567 (95% CI 0 – 0.50%) ^{A, B}

^A 95% confidence intervals are based on a "constant-hazard" exponential failure-time model, whose parameter is determined by the total number of woman-months accumulated during the trial.

^B Combined effectiveness data obtained using Bayesian statistics.

^C The number of women "N" were considered to have completed follow-up at 1 year if patient contact occurred at ≥ 11 months, 2 years if contact occurred at ≥ 23 months, 3 years if contact occurred at ≥ 35 months, 4 years if contact occurred ≥ 47 months and 5 years if contact occurred at ≥ 59 months.

No pregnancies were reported in the clinical trials. However, no method of contraception is 100% effective and pregnancies have occurred in the commercial setting.

Effectiveness rates for the **Essure** procedure compare quite favorably to the effectiveness rates for other methods of tubal sterilization at these time points.

X. ESSURE EFFECTIVENESS IN THE COMMERCIAL SETTING

In the commercial setting, unintended pregnancies have been reported in women who have worn the inserts.

Table 7 summarizes the reasons for pregnancy from reports received by Bayer HealthCare LLC and additional reports from the published scientific literature.

Table 7: Summary of Pregnancies Reported in Commercial Use of Essure*

Potential Contributing Factor	United States		Outside the United States**		Total	
	n	% of US causes	n	% of OUS causes	n	%
Patient Non-compliance (e.g., failure to use alternate contraception or return for Essure Confirmation Test)	213	32%	16	18%	229	31%
Perforation*** / #	91	14%	4	5%	95	13%
Unsatisfactory Placement***	32	5%	13	15%	45	6%
Physician Non-compliance	22	3%	13	15%	35	5%
Pregnant at time of Placement (Luteal)	26	4%	6	7%	32	4%
Inadequate Confirmation Test***	28	4%	0	0%	28	4%
Expulsion***	20	3%	4	5%	24	3%
Tubal Patency***	19	3%	1	1%	20	3%
Insufficient Information to determine	209	32%	31	35%	240	32%
Total	660		88		748****	

*Table includes pregnancy reports received directly by Bayer HealthCare LLC, recorded in the FDA MAUDE database and reported in the scientific literature; data reported to FDA in PMA Annual Reports. Pregnancies in **Essure** patients may be underreported.

Outside of the United States, the **Essure Confirmation Test may be an x-ray or transvaginal ultrasound; device location alone, not occlusion, is primarily used to determine whether the patient may rely on **Essure**.

*** Most of these pregnancies are due to misinterpreted **Essure** Confirmation Tests. Please note that many misinterpretations are due to the fact that occlusion is seen on the HSG films even though the insert is not properly located.

****Number of pregnancies reported from worldwide commercial launch in 2001 through end of 2010. 497,306 **Essure** kits sold during this time. Note that an accurate pregnancy rate is difficult to obtain as the number of devices actually implanted is not known.

*The causal association cannot be established between the perforation and the pregnancy. However, perforations have been identified in pregnant women who were relying on **Essure** for contraception.

The majority of unintended pregnancies are preventable. Most unintended pregnancies are related to patient non-compliance and physician misinterpretation of the **Essure** Confirmation Test. In order to ensure maximum contraceptive effectiveness by **Essure**, the physician should ensure that the patient is properly counseled in accordance with Section XI. It is also important to evaluate both insert location and occlusion carefully before telling the patient that she may rely on **Essure** for contraception.

References

- Adelman MR, Dassel MW, Sharp HT. Management of complications encountered with **Essure** hysteroscopic sterilization: A systematic review. *JMIG*, 2014; 21(5): 733-743.
- Albright CM, Frishman GN, Bhagavath B. Surgical aspects of removal of **Essure** microinsert. *Contraception*, 2013. 88: 334-336.
- Arjona, J.E., et al., *Satisfaction and tolerance with office hysteroscopic tubal sterilization*. *Fertil Steril*, 2008. 90(4): p. 1182-6.
- Arjona, J.E., et al., *Unintended pregnancy after long-term Essure microinserts placement*. *Fertil Steril* 2010. 94(7): p. 2793-5.
- Connor, V.F., *Essure: a review six years later*. *J Minim Invasive Gynecol*, 2009. 16(3): p. 282-90.
- Duffy, S., et al., *Female sterilisation: a cohort controlled comparative study of Essure versus laparoscopic sterilisation*. *BJOG*, 2005. 112(11): p. 1522-8.
- Grosdemouge I., et al., *[Essure implants for tubal sterilisation in France]*[Article in French]. *Gynecol Obstet Fertil*, 2009. 37(5): p. 389-95.
- Guiahi, M. et al., *Retrospective Analysis of Hysterosalpingogram Confirmatory Test Follow-Up after Essure Hysteroscopic Sterilization; 4-year experience in a community setting*. *J Minim Invasive Gynecol*, 2008. 15:S78.
- Hastings-Tolsma, M., P. Nodine, S.B. Teal, and J. Embry, *Pregnancy outcome after transcervical hysteroscopic sterilization*. *Obstet Gynecol*, 2007. 110: p. 504.
- Hur H-C, Mansuria SM, Chen BA, Lee, TT. Laparoscopic management of hysteroscopic **Essure** sterilization complications: Report of 3 cases. *JMIG*, 2008; 15(3): 362-365.
- Langenveld, J., et al., *Tubal perforation by Essure: three different clinical presentations*. *Fertil Steril*, 2008. 90(5): p. 2011 e5-10.
- Lannon BM, Lee, S-Y. Techniques for removal of the **Essure** hysteroscopic tubal occlusion device. *Fertility and Sterility*, 2007; 88(2): 497-498.
- Levie, M.D. and S.G. Chudnoff, *Prospective analysis of office-based hysteroscopic sterilization*. *J Minim Invasive Gynecol*, 2006. 13(2): p. 98-1.
- Mino, M., et al., *Success rate and patient satisfaction with the Essure sterilisation in an outpatient setting: a prospective study of 857 women*. *BJOG*, 2007. 114(6): p. 763-6.
- Moses, A.W., et al., *Pregnancy after Essure placement: report of two cases*. *Fertil Steril*, 2008. 89(3): p. 724 e9-11.
- Ory, E.M., et al., *Pregnancy after microinsert sterilization with tubal occlusion confirmed by hysterosalpingogram*. *Obstet Gynecol*, 2008. 111 (2 Pt 2): p. 508-10
- Ploteau S. and P. Lopes, *Pregnancy after hysteroscopic tubal sterilization despite two hysterosalpingograms showing bilateral occlusion*. *Eur J Obstet Gynecol*, 2009. 147(2): p. 238-9.
- Savage UK, et al., *Hysteroscopic sterilization in a large group practice: experience and effectiveness*. *Obstet Gynecol*, 2009. 114(6): p. 1227-31.
- Shavell, V.I., et al., *Post-Essure hysterosalpingography compliance in a clinic population*. *J Minim Invasive Gynecol*, 2008. 15(4): p. 431-4.
- Veersema, S., et al., *Unintended pregnancies after Essure sterilization in the Netherlands*. *Fertil Steril*, 2010. 93(1): p. 35-8.

Table 8 provides estimates of the percent of women likely to become pregnant while using a particular contraceptive method for one year. These estimates are based on a variety of studies.

Table 8
Pregnancy Rates for Birth Control Methods
(For One Year of Use)

Method	Typical Use Rate of Pregnancy
Sterilization:	
Male Sterilization	0.15%
Female Sterilization	0.5%
Hormonal Methods:	
Implant (<i>Norplant™</i> and <i>Norplant™ 2</i>)	0.05%
Hormone Shot (<i>Depo-Provera™</i>)	3%
Combined Pill (<i>Estrogen/Progestin</i>) and Progestin-only Pill	8%
NuvaRing	8%
Ortho Evra	8%
Intrauterine Devices (IUDs):	
Copper T (<i>ParaGard</i>)	0.8%
LNG-IUS (<i>Mirena</i>)	0.2%
Barrier Methods:	
Male Latex Condom ¹	15%
Diaphragm ²	16%
Female Condom	21%
Spermicide: (gel, foam, suppository, film)	
	29%
Natural Methods:	
Withdrawal	27%
Natural Family Planning (calendar, temperature, cervical mucus)	25%
No Method:	85%

¹ Used Without Spermicide ² Used With Spermicide

Data adapted from Trussell J. Contraceptive efficacy. In Hatcher RA, Trussell J, Nelson AL, Cates W, Stewart FH, Kowal D. Contraceptive Technology: Nineteenth Revised Edition. New York NY: Ardent Media, 2007.

XI. PATIENT SELECTION AND COUNSELING

Consider risks and benefits as described in Sections VIII and IX.

A. Patient selection factors should include:

- Certainty about desire to end fertility
- Evaluation for pelvic infection, cervicitis, undiagnosed vaginal bleeding, anatomical variants and/or uterine pathology that may make patient unsuitable for procedure

The decision to undergo treatment is at patient discretion, following physician counseling and informed consent.

IMPORTANT: Counsel patients that this product does not protect against either HIV infection or other sexually transmitted infections.

B. Patient counseling should include the following:

- Patient must use alternative contraception (except an IUD or IUS) for three months post-procedure until insert placement and fallopian tube occlusion is confirmed by the **Essure** Confirmation Test (modified HSG). Notify patient that alternative contraception is critical due to a theoretical increased risk of ectopic pregnancy. Ensure patient is supplied with the most effective means of contraception for which she is a candidate during this time frame.
- The Patient Information Booklet details risks associated with placement and wearing of the inserts.
- The procedure is permanent, and irreversible.
- Like all birth control methods, there is a risk of pregnancy.
- Bilateral insert placement may not be successful. Discuss a management plan with the patient in the event that bilateral placement is not achieved.

XII. HOW SUPPLIED

Each **Essure** kit contains:

- Two (2) **Essure** systems
- Two (2) **DryFlow** Introducers
- One (1) Instructions for Use (IFU)
- One (1) Patient Identification Card

The **Essure** system is available in one size only.

STERILE: Each **Essure** system is supplied sterile. Do not reuse or resterilize. Inspect each package and do not use if damaged.

STORAGE: Store in a cool, dry place.

XIII. PHYSICIAN TRAINING MANUAL

Refer to the **Essure** Physician Training Manual for detailed information not included in this IFU.

XIV. DIRECTIONS FOR USE

A. Prior to Insert Placement Procedure

1. Perform procedure during early proliferative phase of the menstrual cycle to increase ostia visualization and prevent placement in a patient with an undiagnosed (luteal phase) pregnancy. Women with menstrual cycles shorter than 28 days should undergo careful ovulation day calculations. Insert placement should not be performed during menstruation.
2. Administer a pregnancy test within 24 hours prior to the procedure.
3. Administration of nonsteroidal anti-inflammatory drugs (NSAIDs) is strongly recommended 1-2 hours pre-procedure as this was shown to increase likelihood of placement success in clinical trials. To reduce anxiety, an anxiolytic agent may be given 30 minutes before the procedure.

B. Essure Insert Placement Procedure

1. Use universal precautions and sterile technique.
2. Check all equipment for damage; ensure there are no missing parts.
3. Place patient in the lithotomy position. Use under buttocks drape with fluid control pouch for fluid management.
4. Introduce speculum to allow cervical access. Prep cervix with betadine or other antibacterial solution according to standard practice.
5. Administer anesthesia, preferably local block with or without IV sedation.
6. Insert sterile hysteroscope, with camera and operating channel (≥ 5 French), through cervix into uterine cavity, and remove speculum. Do not perform cervical dilation unless necessary; if necessary, dilate only enough for hysteroscope insertion.
7. Distend uterus with a pre-warmed physiologic saline solution. Maintain uterine distension throughout procedure for optimal visualization. Use gravity feed to minimize tubal spasm and over-distension. Monitor fluids to reduce the risk of hypervolemia; terminate procedure if fluid deficit exceeds 1500cc or hysteroscopic time exceeds 20 minutes.
8. Identify and assess both tubal ostia prior to insert placement. Do not attempt placement in one tubal ostium unless expectation of contralateral tubal patency exists.
9. Insert introducer through sealing cap on the operating channel of the hysteroscope. Keep the channel open to avoid damage to the device and introducer.
10. Insert delivery catheter through the introducer; avoid bending the insert tip. Under direct visualization, advance catheter through the operating channel into the proximal fallopian tube with gentle, constant forward movement to prevent tubal spasm. If excessive resistance occurs (i.e. catheter does not advance toward tubal ostium and/or catheter bends or flexes excessively), terminate procedure to avoid uterine perforation or placement into a false passage.


Figure 5: Insert the introducer through sealing cap on the hysteroscope operating channel, then place delivery catheter through the introducer.

11. Align the proximal end of the black positioning marker with the ostium. Proper alignment leads to ideal placement for the insert, with the outer coil spanning the uterotubal junction. Do not rotate the thumbwheel until the marker is properly aligned.


Figure 6: Advance until black positioning marker at tubal ostium, indicating proper position.

12. Stabilize the delivery system handle against the hysteroscope to prevent inadvertent forward movement. Refer to the Physician Training Manual for techniques for stabilizing the handle.


Figure 7: Stabilize handle against camera head or some other fixed object to prevent inadvertent forward movement of the **Essure** system

Before proceeding with the **Essure** procedure, recall that two distinct operations will take place. The first is retraction of the delivery catheter away from the insert, prior to actual detachment of the insert. Full retraction is accomplished by rotating the thumbwheel to the point where you cannot rotate the thumbwheel any further. Actual detachment is accomplished after retraction by pressing the handle button and then continuing to rotate the thumbwheel. Only after detachment of the insert has occurred can you remove the delivery system.


Figure 8: Rotate thumbwheel to retract delivery catheter, exposing wound down insert

13. Rotate the thumbwheel on the handle toward you until the wheel no longer rotates (corresponds to the symbol ② on the delivery system handle). The delivery catheter and black positioning marker will move away from the tubal ostium and disappear into the operating channel, exposing approximately 1 cm of the wound-down insert.

NOTE: The thumbwheel cannot be reversed, therefore, ensure that the proximal edge of the black positioning marker is aligned at the ostium before rotating thumbwheel.

14. Stop and check proper positioning, which corresponds to the symbol ② on the delivery system handle. Confirm placement of the gold marker band just outside the ostium, and confirm visualization of the distal tip of the green release catheter. If more than 1 cm of the insert is visible in the uterus, then the insert should be repositioned by moving the entire system further into the tube, if possible.

STOP and Check


Figure 9: Visualize gold band at ostium and green release catheter.

15. Press the button on the delivery handle which corresponds to the symbol ③ on the handle button. This enables the thumbwheel to be rotated further for insert deployment.

NOTE: DO NOT PUSH THE BUTTON until delivery system is in correct position for insert placement.


Figure 10: Press button to enable the thumbwheel to rotate again.

16. Rotate the thumbwheel once more (symbol R^4) until the thumbwheel cannot turn any further. When the thumbwheel cannot be rotated any further and the expanded outer coils are visible, remove the delivery system.

Note: Hold the valved introducer in place during removal of the delivery system as it may also be inadvertently withdrawn. If the introducer is removed, replace with a new introducer provided in the **Essure** system packaging.


Figure 11: Rotate thumbwheel to allow outer coils to expand and release the insert from the delivery system

17. Repeat procedure in contralateral fallopian tube.

C. Assessment of Insert After Deployment

1. Assess insert position immediately after deployment. Inserts showing 0-17 trailing coils should be left in place and evaluated via **Essure** Confirmation Test (modified HSG). Ideally, 3 to 8 expanded outer coils should be trailing into the uterus.


Figure 12: Ideal placement of insert

2. If no trailing coils are visible, examine delivery system upon removal from the hysteroscope. Refer to images below to determine if the insert has been deployed from the delivery system. **IMPORTANT:** If insert was inadvertently deployed in the uterine cavity and not in the tube, remove from uterus and attempt another placement.


Figure 13: Delivery systems showing absence of insert after deployment (top) and with insert attached (bottom)

WARNING: Do not attempt insert removal hysteroscopically unless 18 or more coils of the Essure insert are trailing into the uterine cavity. Removal of insert may not be possible; attempted removal of inserts having fewer than 18 trailing coils may cause insert to fracture or patient injury.

3. If insert removal is indicated, perform removal immediately after failed placement as follows:
- As necessary, administer analgesia/anesthesia to reduce or prevent patient discomfort.
 - Introduce a grasping instrument through hysteroscope operating channel.
 - Grasp both the outer and inner coils of the insert together.
 - Withdraw the grasping instrument and hysteroscope simultaneously; the insert may stretch or elongate. Do not pull insert through the operating channel.
4. If insert removal is successful, repeat the **Essure** insert placement procedure. If removal is not accomplished, leave insert in place. If insert is not completely removed, do not place additional inserts. Take a diagnostic X-ray to determine if insert fragment remains *in vivo*. If fragment remains, refer to section XVIII Insert Removal.
5. Document procedural concerns. Review during **Essure** Confirmation Test (modified HSG). Note possible perforations due to:
- excessive or sudden loss of resistance;
 - inability to visualize coils;
 - problems with identification of tubal ostium;
 - poor distension;
 - poor illumination;
 - or poor visualization secondary to endometrial debris.

D. Post Procedure Instructions and Patient Follow-up Requirements

- Ensure patient uses alternative contraception until the **Essure** Confirmation Test (modified HSG). Counsel patient on theoretical increased risk of ectopic pregnancy.
- Schedule patient for an **Essure** Confirmation Test (modified HSG) three months following the procedure to evaluate insert location and fallopian tube occlusion.
- Provide patient with the Patient Identification Card and instruct her to show it to her physicians.

XV. MANAGEMENT OF CASES WITH UNSUCCESSFUL INSERT PLACEMENT

In the event of placement failure, inform patient permanent contraception is not complete.

Counsel patient on undergoing a second procedure, especially if unilateral placement was achieved. In the Pivotal trial, 83% of those who underwent a second procedure achieved bilateral placement. Before a second placement attempt, determine tubal patency by an **Essure** Confirmation Test (modified HSG). Schedule after patient's next menses. If second attempt fails, success with subsequent attempts is unlikely.

If one insert is left *in vivo*, counsel patient to not rely on the insert for contraception. Do not remove a unilaterally placed insert unless the patient experiences an adverse event(s) due to its presence. Three months after bilateral placement, follow-up patient with a second **Essure** Confirmation Test (modified HSG) to verify location and tubal occlusion.

If the patient chooses laparoscopic sterilization, clip or coagulate both fallopian tubes distal or proximal to the insert. Do not perform clipping or coagulation adjacent to or over the insert.

XVI. PERFORMING AND EVALUATING ESSURE CONFIRMATION TEST (MODIFIED HSG) THREE MONTHS POST-INSERT PLACEMENT

The **Essure** Confirmation Test (modified HSG) is performed three months post-placement to evaluate insert location and fallopian tube occlusion. Follow the summary instructions below for performing and evaluating the **Essure** Confirmation Test (modified HSG). For the full protocol, please visit **EssureMD.com** or refer to the Physician Training Manual.

A. Performing the Essure Confirmation Test (modified HSG)

Guidelines:

- Obtain good cornual filling so that uterine cavity silhouette is clearly seen.
- Place fluoroscopy beam as close to A/P projection as possible.
- Do not dilate cervix unless necessary; if dilation occurs, maintain a good cervical seal.
- Downward traction on cervical tenaculum may be required for midpositional uteri. Remove speculum prior to fluoroscopy for best visualization of uterine anatomy.
- Take a minimum of 6 radiographs to assess insert location and tubal occlusion.


Figure 14: Corresponding radiographic view of the **Essure** insert

Radiograph 1 – “Scout Film” - Uterus and inserts without contrast

Radiograph 2 – Minimal Fill of Cavity - uterus and inserts with small amount of contrast

Radiograph 3 – Partial Fill of Cavity - uterus and inserts when nearly full of contrast


Radiograph 4 – Total Fill of Cavity - uterus and inserts when the cornua is distended by contrast

Radiographs 5 & 6 - Magnifications of uterine cornua - insert within the fallopian tube with right (5) and left (6) cornua

CAUTION: Avoid excessive intrauterine pressure beyond Radiograph 4 to avoid undue patient discomfort and vaso-vagal reaction.

NOTE: Immediately repeat **Essure** Confirmation Test (modified HSG) if insert cannot be located or position is unclear.

Figure 15

Radiograph 1: Scout Film

Figure 16

Radiograph 4: Total Fill of Cavity

B. Assessing Insert Location

During evaluation, note four “markers” at each end of the inner and outer coils. Note that the distal markers are fixed in relation to one another, but the proximal markers may move or seem stretched because of the flexibility of the outer coil. Ideal insert location is when the inner coil crosses the utero-tubal junction.

Assess insert location:

- Expulsion or proximal placement: Insert is not present, or $\geq 50\%$ of inner coil trailing into uterine cavity.
- Satisfactory placement: Distal end of inner coil is within the tube, with $< 50\%$ of inner coil trailing into the uterine cavity, or proximal end of inner coil ≤ 30 mm into tube from where contrast fills cornua.
- Distal placement or perforation: Insert is in tube but proximal end of the inner coil > 30 mm distal from where contrast fills cornua or the insert is completely or partially perforated.

C. Assessing Tubal Occlusion

Determine whether the contrast is visible beyond the insert and note any degree of proximal tubal filling even if tube is occluded.

Assess tubal occlusion:

- Satisfactory occlusion: Tube is occluded at the cornua.
- Satisfactory occlusion: Contrast seen within tube but not past distal end of outer coil.
- Unsatisfactory occlusion: Contrast seen past the distal end of the insert or in the peritoneal cavity.

D. Assessing Ability to Rely

If location and tubal occlusion are both rated satisfactory, instruct patient to discontinue alternative contraception.

If location is unsatisfactory, instruct patient to not rely on the inserts for contraception.

If location is satisfactory but occlusion is unsatisfactory, instruct patient to remain on alternative contraception. Repeat the **Essure** Confirmation Test (modified HSG) in three months. If occlusion is still unsatisfactory, instruct patient to not rely on the inserts for contraception.

XVII. MANAGEMENT OF INSERT EXPULSION OR UNSATISFACTORY INSERT LOCATION

1. Insert is too proximal (≥50% of the inner coil trails into the uterine cavity): Do not advise patient to rely on insert, continue alternative contraception or consider incisional sterilization.
2. Insert is too distal (proximal end of the inner coil is > 30mm from cornua): If tube is patent counsel patient on repeat placement procedure. If tube is occluded, advise patient on potential for false positive diagnosis of occlusion. Also counsel patient on incisional sterilization or remaining on alternative contraception.
3. Unilateral or bilateral insert expulsion: If tube is patent, counsel patient on repeat placement procedure. If tube(s) where insert was expelled is occluded, advise patient on potential for false positive diagnosis of occlusion. Also counsel patient on incisional sterilization or remaining on alternative contraception.
4. Unilateral or bilateral insert perforation: If tube is patent, counsel patient on repeat placement procedure. If tube is occluded, advise patient on potential for false positive diagnosis of occlusion. Also counsel patient on incisional sterilization or remaining on alternative contraception.
5. If a patient opts for incisional sterilization, clip or coagulate both fallopian tubes distal or proximal to the insert even if any remaining insert is in a satisfactory location. Do not perform clipping or coagulation adjacent to or over the insert. An attempt should be made to retrieve a perforated insert if the physician believes it can be done safely, however insert retrieval may not be possible. Refer to section XVIII. INSERT REMOVAL.

XVIII. INSERT REMOVAL

WARNING: Essure inserts are intended to be left in place permanently. Do not remove insert(s) unless patient is experiencing an adverse event(s) associated with its presence, or if removal is demanded. If insert removal is indicated, patient should be counseled on unknown risks as techniques for insert removal post-placement have not been evaluated in clinical studies.


For all surgical device removal procedures, care should be taken to avoid transecting the insert during removal. Avoid use of any instrument that is likely to result in fragmentation of the insert. Application of electrocautery to the outer coil should be avoided. During removal, grasping both the inner and outer coils together may help prevent excessive stretching of the outer coil, which could result in fragmentation.

Location of **Essure** inserts should be confirmed through imaging prior to any attempted surgical removal.

Limited case reports describe hysteroscopic insert removal up to seven weeks following placement. In these cases the proximal coils were visible within the uterine cavity and were easily removed with gentle traction. When hysteroscopic removal is not feasible, linear salpingotomy or salpingectomy via laparoscopy or laparotomy can be used to remove an insert within the tube. In some cases, a cornual resection of the proximal fallopian tube may be required for removal. In these cases patients should be counseled about the risk of hysterectomy in order to achieve hemostasis.

1. To perform a linear salpingotomy, make a small incision (approximately 2 cm in length) along the antimesenteric border of the fallopian tube overlying the insert. Use of vasoconstrictive agents is at the discretion of the operating surgeon. The insert needs to be exposed and may need to be freed from the surrounding tissue prior to grasping the coils. Once the insert is exposed, a grasping instrument may be used to extract the insert using gentle traction. Removal may be along with, or independent of, an incisional sterilization procedure.
2. When removing insert via salpingectomy, the location of the proximal and distal portions of the insert within the fallopian tube should be reconfirmed intraoperatively by palpation and/or imaging prior to transecting the tube. The insert may be exposed and visualized via salpingotomy prior to transection or removal of the fallopian tube.
3. The technique for removal of an insert that has perforated the uterus or tube or is within the peritoneal cavity will depend on the location of the insert. Localization should be assessed with imaging prior to the surgical procedure and confirmed intraoperatively. Availability of intraoperative fluoroscopy and/or intraoperative x-ray is recommended to identify the location of the insert or fragments of the insert during surgery.

XIX. SYMBOL LEGEND

	Sterilized using ethylene oxide		Use by		Authorized European Representative
	Batch code		Keep away from heat		Device complies with European Directive
	Do not reuse		MR Conditional		Keep dry
	Catalog number		Do not use if package is open or damaged		Content
	Attention, See Instructions for Use				

Manufactured by:
Bayer HealthCare LLC
 1011 McCarthy Blvd.
 Milpitas, CA 95035
 USA
 essure.com

© 2002 Bayer HealthCare LLC
 All Rights Reserved
 Printed in USA

Essure® and **DryFlow®** are trademarks or registered trademarks of Bayer Essure Inc. in the United States and/or other countries. All other trademarks are the property of their respective owners.

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
INSTRUCTIONS FOR USE

I. Description

The **Essure**[®] Permanent Birth Control System is comprised of several components. The **Essure** micro-insert, a dynamically expanding micro-insert, is attached to a delivery wire and a release catheter. The entire assembly is sheathed within a delivery catheter. This system, (shown in Figure 1), is attached to a handle that facilitates micro-insert delivery and deployment. A valved **DryFlow**[®] Introducer is also provided with the **Essure** system. It is intended to help protect the **Essure** micro-insert as it is being passed through the rubber port of the hysteroscope working channel.

Figure 1
Essure Delivery System
Showing detail of placement procedure symbols
(NOT TO SCALE)

II. Mechanism of action

Under hysteroscopic visualisation, the **Essure** system delivers an **Essure** micro-insert to the proximal section of the fallopian tube lumen. When the **Essure** micro-insert expands on release, it acutely anchors itself in the fallopian tube. Subsequently, the micro-insert elicits an intended benign tissue response, resulting in tissue in-growth into the micro-insert that anchors the micro-insert firmly into the fallopian tube. This benign tissue response is local, fibrotic and occlusive in nature.

Each **Essure** system is sterilised using ethylene oxide and is supplied sterile for single use only. Do not reuse or resterilise. Resterilisation may adversely affect proper mechanical function and could result in patient injury.

III. Indications for use

The **Essure** system is intended for use as a tubal occlusion micro-insert for purposes of permanent contraception.

IV. Contraindications for use

- Patient uncertainty about her desire to end fertility.
- Pregnancy or suspected pregnancy.
- Delivery or termination of a second trimester pregnancy less than 6 weeks before **Essure** micro-insert placement.
- Active or recent pelvic infection.
- Untreated acute cervicitis.
- Unexplained or severe vaginal bleeding.
- Gynaecological malignancy (suspected or known).
- Known abnormal uterine cavity or fallopian tubes that makes visualisation of the tubal ostia and/or cannulation of the proximal fallopian tube difficult or impossible.
- Allergy to contrast media (a hysterosalpingogram may be required three months post-micro-insert placement).
- Patient currently taking corticosteroids.

V. Warnings

- The **Essure** procedure should only be performed by skilled hysteroscopists who have completed the Bayer Healthcare LLC training programme for this procedure.
- Persons allergic to nickel titanium may suffer an allergic reaction to the micro-insert.
- Do not use the **Essure** system if the package is open or damaged. Do not use if the micro-insert is damaged.
- When introducing the **Essure** micro-insert into the fallopian tube, never advance the micro-insert(s) against excessive resistance.
- Do not continue to advance the **Essure** system once the positioning marker on the catheter has reached the tubal ostium. Advancement beyond this point could result in unsatisfactory micro-insert placement or tubal/uterine perforation.
- If a tubal perforation occurs or is suspected, do not continue with the **Essure** micro-insert placement attempt. A very small percentage of women in the **Essure** clinical trials (1.8% or 12/682 patients) were identified as having device related tubal perforations. Retrieval of perforating micro-inserts, if necessary, will require laparoscopy or other surgical methods.
- If **Essure** micro-insert placement attempts are not successful after 10 minutes of attempted cannulation per tube, the case should be discontinued and potentially rescheduled.
- Once the micro-insert has been placed (i.e., detached from the delivery wire), micro-insert removal should not be attempted hysteroscopically unless 18 or more coils of the **Essure** micro-insert are trailing into the uterine cavity. Removal of such a micro-insert should be attempted immediately following the placement. However, removal may not be possible.
- The patient must use alternative contraception until an x-ray performed three months post-micro-insert placement demonstrates satisfactory micro-insert location.
- Patients who undergo placement of the **Essure** micro-insert may, in future years, be offered intrauterine therapies that utilise electrical energy. It is recommended that electrocautery be avoided in surgical procedures undertaken on the uterine cornua and fallopian tubes. All other procedures in the pelvis should avoid the use of electrocautery within 4 cm of the micro-insert. Due to the presence of the **Essure** micro-inserts, there may be risks associated with such procedures that at this time have not been identified.
- Any intrauterine procedure such as endometrial biopsy, D&C, hysteroscopy (diagnostic or operative) including endometrial ablation, could interrupt the ability of the micro-inserts to prevent pregnancy. In addition, the presence of the **Essure** micro-inserts could entail risks associated with such procedures that at this time have not been identified.
- Bench and clinical studies demonstrated that endometrial ablation of the uterus can be performed safely and effectively with the GYNECARE THERMACHOICE[®] Uterine Balloon System, the Hologic NovaSure[™] Endometrial Ablation System and the Boston Scientific Hydro ThermAblator[™] immediately following **Essure** micro-insert placement. No specific studies have been conducted to evaluate **Essure** micro-insert expulsion or contraception rates following combined **Essure** and endometrial ablation procedures.
- Patients may decide, in future years, to undergo in vitro fertilisation (IVF) to become pregnant. The effects of the **Essure** micro-inserts on the success of IVF are unknown. If pregnancy is achieved, the risks of the micro-insert to the patient, to the foetus and to the continuation of the pregnancy are unknown.

* Trademark of ETHICON, INC.

** Trademark of Hologic, Inc.

*** Trademark of Boston Scientific Corporation

VI. Precautions

- Whenever possible, micro-insert placement should be performed during days 7-14 of the menstrual cycle (where day 1 represents the first day of bleeding) in order to enhance visualisation of the fallopian tube ostia and decrease the potential for micro-insert placement in a patient with an undiagnosed pregnancy.
- Unusual uterine anatomy may make it difficult to place the **Essure** micro-inserts.
- In order to reduce the risk of uterine perforation, the procedure should be discontinued if excessive force is required to achieve cervical dilatation.
- Both tubal ostia should be identified and assessed hysteroscopically prior to proceeding to **Essure** micro-

insert placement. No attempt should be made to place a micro-insert in one tubal ostium unless there is a reasonable expectation that the opposite tube is accessible and patent.

- Performing endometrial ablation immediately following placement of **Essure** micro-inserts may increase the risk of post-ablation tubal sterilisation syndrome, a rare condition that has been reported in women with a history of tubal sterilisation who undergo endometrial ablation
- Do not advance the **Essure** system if the patient is experiencing extraordinary pain or discomfort.
- Store the **Essure** system in a cool, dry place.

VII. Possible adverse effects

A. Pregnancy

There is a risk of pregnancy, ectopic pregnancy, and risks associated with the treatment for both. If the patient conceives and chooses to continue an intrauterine pregnancy, she should be informed that the risks of the micro-insert to the patient, to the foetus and to the continuation of the pregnancy are unknown.

B. Risks associated with the micro-insert placement procedure

- Local anaesthesia, oral analgesia/sedation, regional anaesthesia (i.e., spinal, epidural), oral or conscious (intravenous) sedation or general anaesthesia may be administered to the patient to prevent or reduce discomfort. Regardless of the type of anaesthesia, patients may not be able to resume normal activities for 12-24 hours following the procedure.
- Pain, cramping and vaginal bleeding may occur during and following the micro-insert placement procedure. Typically, these incidents are tolerable, transient and successfully treated with medication.
- During and/or directly following the micro-insert placement procedure, there is the risk that the patient will experience nausea or vomiting. This is expected to be transient and may be treated with medication as required.
- Patients may experience fainting or vasovagal response on the day of the procedure.
- There is a risk of perforation or dissection of the fallopian tube or uterine cornua. Bleeding and scarring may result from such a perforation or dissection; however, treatment is typically not required.
- There is a risk of uterine perforation by the hysteroscope, **Essure** system or other instruments used during the procedure with possible injury to the bowel, bladder and major blood vessels. Surgical intervention may be required, but is unlikely, if such injury were to occur. To reduce the risk of uterine perforation, the procedure should be terminated if excessive force is required to achieve cervical dilatation.
- There is a risk that the **Essure** micro-insert may be inadvertently placed into the myometrium of the uterus and not into the fallopian tube lumen. If one micro-insert has already been properly placed in one fallopian tube, in addition to inadvertent placement into the myometrium, the physician may attempt to place a third micro-insert to complete the procedure. If bilateral fallopian tube placement is not achieved, this may result in the patient having one micro-insert in the fallopian tube and/or one micro-insert in the myometrium that cannot be relied upon for contraception. Placement of the micro-insert in the myometrium may result in post-operative pain or other adverse event. If surgical removal of the micro-insert(s) is required, salpingectomy or hysterectomy may be required.
- There is a risk that the **Essure** micro-insert may be placed too distally in the fallopian tube. If removal of the micro-insert is necessary, surgery (laparoscopy or laparotomy) will be required.
- There is a risk that the **Essure** micro-insert may be placed too proximally in the fallopian tube. If 18 or more coils of the **Essure** micro-insert are visible at the time of placement, an immediate attempt should be made to remove the micro-insert (see section XII, **Essure** micro-insert removal). If micro-insert removal is attempted, there is a possibility that the removal will not be successful or that the **Essure** micro-insert may break, leaving a fragment of the micro-insert *in vivo*. If micro-insert removal is attempted and/or achieved, there is also a possibility that the patient may experience increased pain, cramping and bleeding during and following the **Essure** micro-insert placement procedure.
- There is a risk that the **Essure** micro-insert may perforate through the tubal wall or uterine cornua, which could result in the micro-insert being released into the peritoneal cavity. Post-operative pain and/or menstrual disturbance or other adverse event may occur as a result. If the patient elects to undergo incisional sterilisation or other surgical intervention, micro-insert retrieval from the peritoneal cavity may be attempted if the physician believes it is safe to do so. However, micro-insert retrieval may not be possible if the micro-insert cannot be visualised or accessed by the physician.
- There is a risk that **Essure** micro-insert placement will only be achieved in one fallopian tube. If this occurs, patients may be left with one micro-insert *in vivo* that cannot be relied upon for permanent contraception.
- There is a risk that **Essure** micro-insert placement will not be possible in either fallopian tube.
- There is a minimal risk of excess fluid absorption of the physiologic saline fluid used for distension of the uterus, to perform the hysteroscopic procedure.
- As with all invasive procedures, the micro-insert placement procedure can cause an infection. An infection could cause damage to the uterus, fallopian tubes or pelvic cavity. This could require antibiotic therapy, or rarely, hospitalisation or surgery, including hysterectomy.

C. Risks associated with **Essure** micro-insert wearing

- There is a risk that the **Essure** micro-insert could move out of the fallopian tubes. This movement could be expulsion (movement out of the fallopian tube and into the uterine cavity/cervix/vagina or out of the body) or migration (movement to the distal fallopian tube or out of the fallopian tube and into the peritoneal cavity). Additional x-rays may be required to identify the location of the micro-insert(s), and surgery may be required to remove the micro-insert(s). Device movement could result in pregnancy, ectopic pregnancy and/or pain/menstrual disturbance or other adverse events.
- As with currently available methods of mechanical permanent contraception (i.e., clips, rings), if the **Essure** micro-insert is to be removed, surgery will be required. Further, it is possible that surgical removal of the fallopian tubes (salpingectomy) and uterus (hysterectomy) may be required.
- Abdominal/pelvic pain and cramping may occur. Pain and cramping may be a more likely occurrence during the menstrual period, during and after sexual intercourse or with other physical activity.
- Intermenstrual bleeding or heavier than normal menstrual bleeding may be experienced.
- Occasionally, a woman may regret her decision to undergo permanent contraception and experience mild depression or other emotional disturbances as a result.

D. Risks associated with follow-up procedures

- There is the risk of radiation associated with the pelvic x-ray that is required three months following micro-insert placement to evaluate micro-insert location. There may also be a need for an **Essure** Confirmation Test (HSG). There are approximately .033 rads in the fluoroscopic portion (< 30 seconds) of a hysterosalpingogram procedure. As a point of comparison, radiation exposure from a barium enema is 0.85 rads which is higher than the required **Essure** Confirmation Test (HSG). The amount of radiation exposure from one pelvic x-ray is about the same as the amount an individual would receive from one year of natural background radiation.
- The following additional risks are associated with the **Essure** Confirmation Test (HSG) procedure if needed: vasovagal response; infection, which may require antibiotic treatment and in rare cases could require hospitalisation; intravasation; perforation of the uterus; uterine cramping and/or bleeding; pain or discomfort; allergic reaction to latex. Latex exposure has been reported to be associated with anaphylactic reactions in rare cases, which may lead to death.
- The use of contrast media, used to perform the **Essure** Confirmation Test (HSG), has been associated with allergic reaction in some patients. Allergic reaction can result in hives or difficulty breathing. In some individuals, an anaphylactic response may occur which may lead to death.

E. Risks associated with potential future procedures

- Patients who undergo placement of the **Essure** micro-insert may, in future years, be offered intrauterine therapies that utilise electrical energy. It is recommended that electrocautery be avoided in surgical procedures undertaken on the uterine cornua and fallopian tubes. All other procedures in the pelvis should avoid the use of electrocautery within 4 cm of the micro-insert. Due to the presence of the **Essure** micro-inserts, there may be risks associated with such procedures that, at this time, have not been identified.
- Any intrauterine procedure such as endometrial biopsy, D&C, hysteroscopy (diagnostic or operative) including endometrial ablation, could interrupt the ability of the micro-inserts to prevent pregnancy. In addition, the presence of the **Essure** micro-inserts could entail risks associated with such procedures that, at this time, have not been identified.
- Patients may decide, in future years, to undergo in vitro fertilisation (IVF) to become pregnant. The effects of the **Essure** micro-inserts on the success of IVF are unknown. If pregnancy is achieved, the risks of the micro-insert to the patient, to the foetus and to the continuation of a pregnancy are unknown.
- The **Essure** micro-inserts are MR-safe and radiopaque. The **Essure** micro-inserts are also MR-compatible, except for pelvic imaging where they may cause some artefacts.
- **There is the potential that unknown risks exist.**

VIII. Directions for use


A. Prior to micro-insert placement procedure

1. Micro-insert placement should be performed during days 7-14 of the menstrual cycle (where day 1 represents the first day of bleeding), in order to enhance visualisation of the fallopian tube ostia and decrease the potential for micro-insert placement in a patient with an undiagnosed pregnancy.
2. A pregnancy test administered by the physician or designee, should be conducted within 24 hours prior to or immediately preceding the micro-insert placement procedure.
3. Administration of a non-steroidal anti-inflammatory drug (NSAID) such as Indocid (orally or via suppository) is strongly recommended one to two hours before the micro-insert placement procedure, since clinical trial data demonstrate that the use of NSAIDs significantly increase the likelihood of placement success. If using only a paracervical block, Diazepam (PO) or a similar agent may also be offered 30 minutes prior to the procedure to reduce anxiety.

B. Essure micro-insert placement procedure


The **Essure** micro-insert placement procedure can be performed in an ambulatory or day surgery setting. Sterile techniques should be used during the micro-insert placement procedure. The amount of time required to complete the micro-insert placement procedure should not exceed 30 minutes.

1. Place the patient in the lithotomy position.
2. Introduce a speculum into the vagina to allow access to the cervix. Prep the cervix with betadine or other suitable antibacterial solution according to standard practice.
3. Local anaesthesia is the preferred method for implantation of the micro-inserts. A paracervical block may be administered. Midazolam (IV), or a similar agent, may also be administered to prevent or reduce discomfort if needed.
4. Insert a sterile hysteroscope, with attached camera and operating channel (≥ 5 French), through the cervix into the uterine cavity. If necessary, perform cervical dilation to allow insertion. In order to prevent uterine perforation, the procedure should be discontinued if excessive force is required to achieve cervical dilatation.
5. Uterine cavity distension should be accomplished with a physiologic saline infusion through the working channel of the hysteroscope. It is strongly recommended that the saline solution be pre-warmed to body temperature and introduced under gravity feed to minimise spasm of the fallopian tubes. Excellent uterine distension must be achieved and maintained throughout the procedure. Standard fluid monitoring procedures should be followed throughout the procedure. The fallopian tube ostia should be identified by hysteroscopic visualisation.
6. Both tubal ostia should be identified and assessed hysteroscopically prior to proceeding to **Essure** micro-insert placement. No attempt should be made to place a micro-insert in one tubal ostium unless there is a reasonable expectation that the opposite tube is patent.
7. Once the fallopian tube ostia have been identified, insert the introducer through the sealing cap on the hysteroscope working channel. The operating channel stopcock should remain in the open position (the device and/or introducer can be damaged if the stopcock closes on either device). Place the **Essure** delivery system through the introducer and advance through the operating channel of the hysteroscope. If undamaged from the first micro-insert placement, the valved introducer may remain in the operating channel throughout the **Essure** procedure.


Insert the introducer through the sealing cap on the hysteroscope working channel, then place the **Essure** delivery system through the introducer.

8. Advance the **Essure** delivery system into the proximal fallopian tube with slow, steady movement to prevent tubal spasm. Advance the delivery system until the positioning marker on the delivery catheter reaches the fallopian tube ostium. This visual marker indicates that the **Essure** micro-insert is spanning the distal intramural to proximal isthmic segments of the fallopian tube, with the outer coil spanning the uterotubal junction. This is the ideal placement for the **Essure** micro-insert.

Advance until the black positioning marker is at tubal ostium. This is a visual indicator for proper position for deployment.


9. Proper concentric alignment of the delivery catheter with the tubal lumen is suggested by the ability to advance the catheter under direct visualisation without undue resistance. Resistance to advancement is usually apparent in two ways: 1) the black marker on the outside surface of the catheter is seen not to advance forward towards the tubal ostium, and/or 2) the delivery catheter bends or flexes excessively, thus preventing the physician from applying forward pressure on the catheter assembly. When such resistance to forward motion of the catheter is observed, no further attempts should be made to place the micro-insert in order to avoid the possibility of uterine perforation or inadvertently placing the micro-insert in the uterine musculature rather than within the tubal lumen. A follow-up **Essure** Confirmation Test (HSG) should be undertaken to determine tubal patency.
10. If it is not possible to advance the catheter to the positioning marker after several minutes, a perfusion test with a patency catheter may be employed, if it has not already been utilised, to determine tubal patency. If the tube is blocked or the catheter cannot be advanced to the positioning marker, the case should be terminated. If micro-insert placement is not successful after 10 minutes of attempted cannulation per tube, the case should be terminated.
11. When the delivery catheter has been advanced to the positioning marker, deploy the micro-insert. To do so, first stabilise the handle of the **Essure** micro-insert against the hysteroscope camera or some other fixed object, to prevent inadvertent forward movement of the **Essure** system during retraction of the delivery catheter.

Stabilise handle against camera head or some other fixed object to prevent inadvertent forward movement of the **Essure** system

12. Being certain that the black positioning marker is at the fallopian tube ostium, rotate the thumbwheel on the handle toward you until the wheel no longer rotates. This operation corresponds to the symbol ② on the delivery system handle. This facilitates withdrawal of the delivery catheter. You will see the black positioning marker move away from the tubal ostium (towards the hysteroscope) and disappear into the operating channel. Withdrawal of the delivery catheter exposes the wound-down **Essure** micro-insert. Approximately 1 cm of the micro-insert (wound-down coils) should appear trailing into the uterus when the delivery catheter is withdrawn.


Rotate Thumbwheel

Rotate thumbwheel to retract catheter


13. To confirm proper positioning, place gold marker band just outside the ostium, which corresponds to the symbol ② on the delivery system handle. Visualisation of the gold band just outside the ostium, as well as visualisation of the distal tip of the green release catheter, will confirm proper positioning. If more than 1 cm of the micro-insert is visible in the uterus, then the micro-insert should be repositioned by moving the entire system further into the tube, if possible, before proceeding to the next step.

STOP and Check


Visualise gold band at ostium

14. Press the button on the delivery handle to enable the thumbwheel to be further rotated, which corresponds to the symbol ③ on the handle button.


Press button to enable thumbwheel to rotate again.

15. Rotate the thumbwheel toward you to deploy the outer coil of the micro-insert, which corresponds to the symbol ④ on the delivery system handle. Continue to rotate the thumbwheel until it stops rotating. When the thumbwheel cannot be rotated any further and the expanded outer coils are visible, withdraw the system.

Rotate thumbwheel to deploy the outer coil of the micro-insert

16. The position of the deployed **Essure** micro-insert will be assessed under hysteroscopic visualisation. There should ideally be 3 to 8 expanded outer coils of the **Essure** micro-insert trailing into the uterus.

Expanded outer coils of the **Essure** micro-insert trailing into the uterus indicates ideal placement

17. If the physician is dissatisfied with the micro-insert's placement based on the hysteroscopic view, or suspects tubal or uterine perforation, the micro-insert(s) should be left in place and evaluated via pelvic X-ray or HSG three months post device placement.

WARNING: AFTER THE MICRO-INSERT HAS BEEN PLACED AND RELEASED INTO THE FALLOPIAN TUBE, DO NOT ATTEMPT TO REMOVE THE MICRO-INSERT HYSTEROSCOPICALLY UNLESS 18 OR MORE COILS OF THE ESSURE MICRO-INSERT ARE TRAILING IN THE UTERINE CAVITY. Removal of such a micro-insert should be attempted immediately during the placement attempt. However, removal may not be possible (see section XIII, **Essure** Micro-insert Removal). If the micro-insert was inadvertently deployed in the uterine cavity and not into the tube, the micro-insert should be removed from the uterus and another attempt made at micro-insert placement in the tube.


18. Repeat the **Essure** micro-insert placement procedure in the contralateral fallopian tube.
19. Record the length of the micro-insert trailing into the uterine cavity, noting any issues with identifying or confirming either tubal ostium or any concerns regarding potential perforation. These should be noted in patient records for subsequent reference when reviewing the **Essure** Confirmation Test (See Section IX - *Essure Confirmation Test* below).
20. **Remind the patient to use an alternative form of contraception (except an IUD) for the first 3 months following the micro-insert placement procedure.**
21. Schedule the patient for an **Essure** Confirmation Test three months following the **Essure** micro-insert placement procedure to evaluate micro-insert retention and location.

IX. Essure Confirmation Test

- A. An **Essure** Confirmation Test should be performed three months after micro-insert placement to evaluate micro-insert retention and location. The **Essure** Confirmation Tests (transvaginal ultrasound (TVU), pelvic x-ray or hysterosalpingogram (HSG)) should be performed only by an experienced gynecologist, ultrasonographer and/or radiologist trained on the respective **Essure** confirmation test protocol. A detailed protocol with images and test performance tips is provided with training; additional copies may be obtained by downloading a copy from essure.com.
- B. For the first-line confirmation test, either a pelvic x-ray or a TVU may be performed three months after an uncomplicated bilateral micro-insert placement procedure.
 1. X-ray and TVU should not be used as the **Essure** Confirmation Test under the following circumstances:
 - a) Difficult placement procedure including one or more of the following:
 - (1) Concern at the time of placement of possible perforation due to excessive force required for micro-insert delivery and/or a sudden loss of resistance.
 - (2) Difficulty identifying the tubal ostia during placement due to anatomical variation or technical factors such as poor distention, suboptimal lighting or endometrial debris.
 - (3) Surgeon is uncertain about placement.
 - b) Procedure time > 15 minutes (scope in-scope out).
 - c) Placement with zero or > 8 trailing coils
 - d) Unusual post-operative pain, transient or persistent, or onset at some later time post procedure, without any other identifiable cause.
 2. If X-ray or ultrasound is not indicated, patient must proceed to an HSG to evaluate micro-insert location and tubal occlusion. Trans-abdominal ultrasound cannot be substituted for TVU. If X-ray or ultrasound evaluation is equivocal or unsatisfactory, patient must proceed to an HSG to evaluate micro-insert location and tubal occlusion.


C. Transvaginal Ultrasound

1. A minimum of three images must be obtained and retained for documentation:
 - a) A coronal or oblique coronal view demonstrating a portion of each micro-insert in the cornua labeled "scout image":


Bilateral micro-inserts are identified in this transverse (coronal / oblique coronal) view.

- b) A coronal or oblique coronal image of the linear axis of the left micro-insert including the proximal end crossing the myometrium in the cornua (interstitial portion of the fallopian tube) or in contact with the utero-serosal tubal junction and labeled "left."
 - c) A coronal or oblique coronal image of the linear axis of the right micro-insert crossing the myometrium in the cornua (interstitial portion of the fallopian tube) or in contact with the utero-serosal tubal junction and labeled "right."
 - d) All three images should be captured on film and placed in the subject's medical record to document satisfactory micro-insert retention and location.
2. Classification of Micro-insert Location
 - a) Micro-insert identification: In a single scout image, a portion of each micro-insert must be visualized in the cornua in the coronal or oblique coronal view to ensure bilateral placement and reduce the risk of duplicate imaging of the same micro-insert. The linear axis of the micro-inserts should appear relatively symmetric.
 - b) Optimal Location
Micro-insert location is optimal when the proximal end of the micro-insert is in contact with the uterine cavity or endometrium, and the linear axis is within the myometrium in the cornua (interstitial portion of the fallopian tube) and can be visualized at or crossing the utero-serosal tubal junction (USTJ). The portion of the micro-insert located in the fallopian tube may or may not be visualized. The linear axis of the micro-insert must be visualized to confirm it is not coiled or elongated.

Optimal Location

- c) Satisfactory Location
Micro-insert location is satisfactory when the proximal end of the micro-insert is distal to the endometrium, however the linear axis is within the myometrium in the cornua (interstitial portion of the fallopian tube) and can be visualized at or crossing the utero-serosal tubal junction (USTJ). The portion of the micro-insert located in the fallopian tube may or may not be visualized. The linear axis of the micro-insert must be visualized to confirm it is not coiled or elongated.

Satisfactory location


- d) Unsatisfactory Location
 - (1) Micro-insert location is unsatisfactory if a portion of each micro-insert cannot be visualized in the cornua in the coronal or oblique coronal view in one scout image.

- (2) Expulsion is suspected if one or both micro-inserts are not identified in the cornua in a coronal view in a single scout image.
 - (3) Distal placement is suspected if the proximal end of the micro-insert is not located in the myometrium in the cornua (interstitial portion of the fallopian tube), and not crossing or in contact with the USTJ.
 - (4) Proximal placement is suspected if greater than 50% or the majority of the micro-insert is visualized in the uterine cavity or if the linear axis of the micro-insert(s) is visualized in the midline sagittal view.
 - (5) Perforation is suspected if the linear axis of one or both micro-inserts are parallel to the endometrial stripe in the sagittal view, or if the linear axis of a micro-insert is visualized crossing the myometrium in the midline sagittal view.
 - (6) Unclassified position: If the linear axis of a micro-insert cannot be identified, suggesting it is coiled, bent or elongated, micro-insert location is considered unsatisfactory. If the surrounding soft tissue cannot be clearly defined, position is considered unsatisfactory.
3. If ultrasound evaluation is equivocal or unsatisfactory, patient must proceed to an HSG to evaluate micro-insert location and tubal occlusion.

D. Pelvic X-ray

1. Capture an image of the uterus with both Essure micro-inserts clearly seen. The lie and curvature of the micro-inserts should be noted.

Corresponding radiographic view of the Essure micro-insert


2. Evaluate pelvic X-ray as follows:
 - a) Satisfactory: Micro-inserts appear to be in the tubal lumen and spanning the uterotubal junction, and appear relatively symmetrical. Patients whose X-rays are determined to be "satisfactory" may begin to rely on the Essure micro-insert for contraception.
 - b) Suspicious: One or both of the micro-inserts appear to be distal or proximal to optimal position, or may be partially or completely perforated through the tube, and/or appear relatively asymmetrical. Patients whose X-rays are determined to be "suspicious" should be instructed to continue alternative contraception and undergo an HSG.
 - c) Unsatisfactory: Obvious intraperitoneal micro-insert location or expulsion.
3. If x-ray evaluation is equivocal or unsatisfactory; or micro-insert location is suspicious, patient must proceed to an HSG to evaluate micro-insert location and tubal occlusion.

E. Performing and Evaluating modified HSGs

1. The HSG is performed to further evaluate Essure micro-insert location and fallopian tube occlusion if necessary based on the X-ray or ultrasound findings. Follow the instructions below for performing and evaluating the HSG.
2. Performing the HSG - Guidelines:
 - a) Obtain good cornual filling so that the uterine cavity silhouette is clearly seen.
 - b) Place fluoroscopy beam as close to A/P projection as possible.
 - c) Do not dilate cervix unless necessary; if dilation occurs, maintain a good cervical seal.
 - d) Downward traction on cervical tenaculum may be required in for midpositional uteri. Remove speculum prior to fluoroscopy for best visualization of uterine anatomy.
 - e) Take a minimum of six radiographs to assess micro-insert location and tubal occlusion.
 - (1) Radiograph 1 - "Scout Film" - Uterus and micro-inserts without contrast.
 - (2) Radiograph 2 - Minimal Fill of the Cavity - Uterus and micro-inserts with small amount of contrast.
 - (3) Radiograph 3 - Partial Fill of the Cavity - Uterus and micro-inserts when nearly full of contrast.
 - (4) Radiograph 4 - Total Fill of Cavity - Uterus and micro-inserts when the cornua is distended by contrast.
 - (5) Radiographs 5 & 6 - Magnifications of uterine cornua - Micro-insert within the fallopian tube with right (5) and left (6) cornua.

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3. Assessing Micro-insert Location
 - a) During evaluation, note four "markers" at each end of the inner and outer coils). Note that the distal markers are fixed in relation to one another, but the proximal markers may move or seem stretched because of the flexibility of the outer coil. Ideal micro-insert location is when the inner coil crosses the utero-tubal junction.
 - b) Assess micro-insert location:
 - (1) Expulsion or proximal placement: Micro-insert is not present or $\geq 50\%$ of inner coil trailing into the uterine cavity.
 - (2) Satisfactory placement: Distal end of inner coil is within the tube, with < 50% of inner coil trailing into the uterine cavity or proximal end of inner coil ≤ 30 mm into the tube from where contrast fills cornua.
 - (3) Distal placement or perforation: Micro-insert is in tube but proximal end of inner coil > 30 mm distal from where contrast fills cornua or the micro-insert is completely or partially perforated.
4. Assessing Tubal Occlusion
 - a) Determine whether the contrast is visible beyond the micro-insert and note any degree of proximal tubal filling even if the tube is occluded.
 - b) Assess tubal occlusion:
 - (1) Satisfactory occlusion: Tube is occluded at the cornua.
 - (2) Satisfactory occlusion: Contrast seen within tube but not past distal end of outer coil.
 - (3) Unsatisfactory occlusion: Contrast seen past the distal end of the micro-insert or in the peritoneal cavity.
5. Assessing Ability to Rely
 - a) If location and tubal occlusion are both rated satisfactory, instruct patient to discontinue alternative contraception.

- b) If location is unsatisfactory, instruct patient to not rely on the micro-inserts for contraception.
- c) If location is satisfactory but occlusion is unsatisfactory, instruct patient to remain on alternative contraception. Repeat the HSG in three months. If occlusion is still unsatisfactory, instruct patient to not rely on the micro-inserts for contraception.

X. Management of Unsatisfactory Micro-insert Location (UML)

A. Unsatisfactory micro-insert location diagnosed by hysterosalpingogram

1. Proximal Location: more than 50% of the length of the inner coil of the micro-insert(s) is trailing into the uterus.
2. Distal Location: the micro-insert(s) is in the fallopian tube but the proximal end of the inner coil is more than 30 mm from the contrast filling the uterine cornua.
3. Complete micro-insert(s) expulsion: micro-insert(s) absent from the body
4. Perforation: micro-insert(s) partially or fully perforated.
5. Intraperitoneal micro-insert(s) location: micro-insert(s) obviously outside the fallopian tube(s).

B. Management of micro-insert expulsion or unsatisfactory micro-insert location

1. Bilateral micro-insert expulsion with bilateral occlusion: the patient should be counselled about the option to have incisional sterilisation or to rely on her bilateral PTO for contraception, in light of the potential for a false positive diagnosis of tubal occlusion by **Essure** Confirmation Test (HSG).
2. Bilateral micro-insert expulsion with occlusion in one tube and patency in contralateral tube: the patient may be considered for an additional micro-insert procedure to replace the micro-insert in the tube that is patent, so that she may be able to rely on one **Essure** micro-insert and contralateral PTO for contraception. The patient should be counselled regarding this option, in light of the potential for a false positive diagnosis of tubal occlusion by **Essure** Confirmation Test (HSG). She should also be counselled about the option to have incisional sterilisation.
3. Unilateral micro-insert expulsion or unsatisfactory unilateral micro-insert location (in myometrium or intraperitoneal cavity), with contralateral micro-insert in a satisfactory location: if the **Essure** Confirmation Test (HSG) demonstrates tubal blockage in the tube from where the micro-insert was expelled or where the micro-insert should have been placed, the patient may rely on the satisfactorily located micro-insert and the contralateral PTO, in light of the potential for a false positive diagnosis of tubal occlusion by **Essure** Confirmation Test (HSG). She should also be counselled regarding the option to undergo incisional sterilisation.
4. Unsatisfactory unilateral micro-insert location (in myometrium or intraperitoneal cavity) with contralateral micro-insert in a satisfactory location: if the **Essure** Confirmation Test (HSG) demonstrates tubal patency in the tube that should have been placed with a micro-insert, the patient may be offered the opportunity to return for an additional micro-insert placement procedure to re-attempt placement. She should also be counselled regarding the option to undergo incisional sterilisation.
5. Unilateral micro-insert expulsion; unsatisfactory unilateral micro-insert location (in myometrium or intraperitoneal cavity); unsatisfactory unilateral micro-insert location in "proximal location" (> 50% of inner coil length trailing into uterus) or "distal location" (micro-insert in fallopian tube, but proximal end of inner coil is > 30 mm from contrast filling the uterine cornua) with contralateral micro-insert in an unsatisfactory location: the patient should be counselled regarding the option to undergo incisional sterilisation. In all cases, if the micro-insert removal is deemed necessary and hysteroscopic removal is not possible, incisional surgery may be required.
6. If a patient has opted for incisional sterilisation following any of the above listed scenarios, both tubes should be occluded regardless of any remaining micro-insert that is in a satisfactory location. An attempt should be made to retrieve a micro-insert if the physician believes it can be done safely, however micro-insert retrieval may not be possible. Use of intra-operative fluoroscopy is recommended to identify the location of the micro-insert(s) prior to and during surgery. Attempted retrieval should not exceed 30 minutes.

XI. Management of cases with unsuccessful Essure micro-insert placement

In the event of unilateral or bilateral micro-insert placement failure, the patient should be informed that her permanent contraception has not been completed. If the patient chooses laparoscopic sterilisation (i.e., clip application or electrical cautery), both fallopian tubes should be clipped or cauterised even if one tube has the **Essure** micro-insert implanted in it. Clipping or cauterising of the tube or tubes should be performed distal to the **Essure** micro-insert.

If the patient does not select laparoscopic sterilisation, she may be offered an **Essure** Confirmation Test (HSG) after her next menses (pre-ovulatory: day 7-14 where day 1 represents the first day of bleeding) to determine tubal patency. If tubal patency is observed, the physician may offer the patient a second attempt at micro-insert placement. If a second attempt at micro-insert placement fails, the patient is unlikely to have success with subsequent attempts. If the patient has one micro-insert left *in vivo*, she should be counselled not to rely on the unilateral micro-insert for contraception.

If only unilateral placement was achieved and the **Essure** Confirmation Test (HSG) confirms contralateral proximal tubal occlusion (PTO), the patient should be counselled regarding the option to rely on the one micro-insert, in light of the potential for a false positive diagnosis of PTO by **Essure** Confirmation Test (HSG). Tubal occlusion is defined as the failure of dye to pass from the uterine cavity into the peritoneal cavity at the time of an **Essure** Confirmation Test (HSG). She should also be counselled about the option to have incisional sterilisation. An attempt to remove a unilaterally placed micro-insert is not recommended, unless the patient is experiencing an adverse event(s) with the micro-insert.

XII. Essure micro-insert removal

WARNING: MICRO-INSERT REMOVAL SHOULD NOT BE ATTEMPTED HYSTEROSCOPICALLY ONCE THE MICRO-INSERT HAS BEEN PLACED, UNLESS 18 OR MORE COILS OF THE ESSURE MICRO-INSERT ARE TRAILING INTO THE UTERINE CAVITY. Removal of such a micro-insert should be attempted immediately following placement. However, removal may not be possible. If removal is attempted, the following steps should be employed:

1. Introduce a grasping instrument through the hysteroscope working channel.
2. Grasp the outer coil of the **Essure** micro-insert. Try to grasp the outer and inner coil of the micro-insert together.
3. Pull back on the grasping instrument and the hysteroscope at the same time, so that the entire system is withdrawn from the uterus in concert.
4. The outer coil and/or the inner coil of the **Essure** micro-insert may stretch or elongate as micro-insert removal is being attempted.
5. As necessary, administer analgesia/anaesthesia to reduce or prevent patient discomfort.
6. If complete micro-insert removal is accomplished, an attempt should be made to place another **Essure** micro-insert.
7. If the physician is not completely satisfied that the entire **Essure** micro-insert has been removed from the fallopian tube, another micro-insert should **NOT** be placed in that tube and a post-placement x-ray should be taken to determine if a micro-insert fragment remains *in vivo*.

Other than the above described scenario, micro-insert removal should only be attempted if a patient is experiencing an adverse event(s) with the micro-insert or if she demands micro-insert removal.

Should micro-insert removal be deemed necessary, a transabdominal approach (i.e., laparotomy or laparoscopy) is required.

A cornual resection of the proximal fallopian tube will be required if the micro-insert is properly located across the utero-tubal junction (UTJ).


An **Essure** micro-insert that has been improperly placed or has migrated beyond the UTJ should be removed with traditional linear salpingotomy or salpingectomy accomplished via laparoscopy or laparotomy.

1. To perform a linear salpingotomy, a small incision (approximately 2 cm in length) is made along the antimesenteric border of the fallopian tube, directly overlying the micro-insert.
2. Total or partial salpingectomy can be performed to retrieve the micro-insert along with, or independent of, the performance of a traditional tubal sterilisation procedure.

XIII. Patient identification card

Each patient who has had **Essure** micro-insert(s) implanted should be given a laminated, wallet-sized card stating that she has **Essure** micro-insert(s) in place. **The card is enclosed in this package.** The card will additionally state that there may be risks associated with the participant undergoing future intrauterine procedures or surgery of the reproductive organs.

XIV. Legend of Symbols

	Sterilised using ethylene oxide
	Batch code
	Do not reuse
	Catalogue number
	Attention, see Instructions for Use
	Use by
	Keep away from heat
	Do not use if package is open or damaged

MR Conditional

Authorised European Representative

Device complies with European Directive 93/42/EC

Keep dry

Content

Manufactured by:
Bayer Healthcare LLC
 1011 McCarthy Boulevard
 Milpitas, CA 95035
 USA
 essure.com

European Authorised Representative:

Bayer Pharma AG
 Müllerstr. 178
 D-13353 Berlin
 Germany

For product ordering or reporting of adverse events, please contact local Bayer Healthcare LLC representative.

© 2002 Bayer Healthcare LLC
 All Rights Reserved
 Printed in USA

Essure® and **DryFlow**® are trademarks or registered trademarks of Bayer Essure Inc. in the United States and/or other countries. All other trademarks are the property of their respective owners.

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PN-84731145, ART Rev. B

10.7 Appendix 7: Patient Information Brochure

**Your complete Guide to the Essure[®]
Procedure**

When you're done
having children,
consider Essure[®]
permanent birth control

YOUR COMPLETE GUIDE TO THE ESSURE[®] PROCEDURE

essure[®]

FOR MORE INFORMATION

Visit essure.com
or call 1-888-842-2937

TABLE OF CONTENTS

Glossary..... 3

What Is Essure®?..... 4

The Benefits of Essure..... 5

Essure Procedure Overview 6

Warnings, Precautions and Other Potential Risks 7

What to Expect With Essure 9

Patient Identification Card..... 10

Frequently Asked Questions 10

Questions to Ask Your Doctor 11

How Essure Performed in Clinical Studies..... 12

How Essure Has Performed Since FDA Approval in 2002..... 13

Comparison of Permanent Birth Control Methods..... 15

Comparison of Temporary Birth Control Methods..... 17

GLOSSARY

Anesthesia: Medically-induced partial or complete loss of sensation in all or part of the body. Loss of sensation may occur with or without loss of consciousness.

Cervix: The passageway that connects the vagina with the uterus.

Contraceptive: Any process, device, or method that reduces the likelihood of pregnancy.

Ectopic Pregnancy: The development of a fertilized egg outside the uterus, such as in a fallopian tube. Ectopic pregnancies can be dangerous and possibly life-threatening.

Essure Insert: The small, soft, flexible device that is placed in your fallopian tubes for permanent pregnancy prevention.

Fallopian Tubes: The tubes that carry the eggs from the ovaries to the uterus.

General Anesthesia: Medication that induces total loss of consciousness and sensation.

Hysterosalpingogram (HSG): An X-ray of the uterus and fallopian tubes after they have been filled with contrast dye.

Hysteroscope: A telescopic instrument that is used to view the inside of the uterus.

In Vitro Fertilization (IVF): Fertilization of an egg outside the body. Once fertilized, the egg is placed into the uterus.

Intrauterine Device (IUD)/Intrauterine Contraceptive (IUC): A medical device that is placed in the uterus for temporary prevention of pregnancy.

Local Anesthetic: Medication that is applied or injected to numb a certain part of the body.

NovaSure® Endometrial Ablation: A procedure that removes the lining of the uterus to lighten or stop your periods.

Occlusion: An obstruction or a closure of a passageway or a vessel.

Tubal Ligation: A form of permanent birth control by means of cutting, tying, burning or clipping the fallopian tubes so that they are blocked.

Uterus: The womb, where a developing fetus grows.

Vasectomy: Permanent birth control for men that involves cutting or blocking a segment of the vas deferens (the tubes that carry the sperm).

WHAT IS ESSURE?

Essure is a permanent birth control procedure that works with your body to create a natural barrier against pregnancy. The Essure procedure involves placing soft, flexible inserts into your fallopian tubes. Over a period of about three months, tissue forms around the inserts. The build-up of tissue creates a barrier that keeps sperm from reaching the eggs and prevents conception.

Essure may be right for you if:

- You are certain you do not want any more children.
- You desire a permanent form of birth control.
- You would like to stop worrying about getting pregnant.
- You prefer a method or procedure that
 - Is simple and does not take a lot of time.
 - Can be done in your doctor's office.
 - Does not require surgery or exposure to its potential risks.
 - Does not contain any hormones.

Essure is NOT right for you if:

- You suspect you are pregnant.
- You have only one fallopian tube.
- You have one or both fallopian tubes closed or obstructed.
- You have had your “tubes tied” (tubal ligation).
- You are allergic to contrast dye used during x-ray exams.
- You are unwilling to undergo the Essure Confirmation Test.
- You are uncertain about ending your fertility.

You should delay having the Essure procedure if:

- You are or have been pregnant within the past 6 weeks.
- You have had a recent pelvic infection.
- You are in the second half (weeks 3 and 4) of your menstrual cycle. During that time, there is an increased risk of being pregnant prior to having the Essure procedure.

You should speak to your doctor if:

- You are taking or receiving therapy that suppresses your immune system. Examples include chemotherapy or corticosteroids, such as prednisone. Therapy that suppresses the immune system may make the Essure procedure less effective for birth control.
- You have, or think you may have, a nickel allergy.

Talk to your Doctor about the Essure procedure and whether it is right for you.

IMPORTANT: Essure inserts do not protect against HIV or other sexually transmitted diseases.

THE BENEFITS OF ESSURE

Highly Effective

The Essure procedure is 99.83% effective based on five-year clinical study data.

Available for Over 10 Years

Since FDA approval in 2002, over half a million women have relied on Essure.

Non-Surgical Procedure

Essure is a simple procedure that can be done in 10 minutes in your doctor's office.

Non-Hormonal

Essure inserts do not contain or release hormones.

No General Anesthesia Required

You can remain fully conscious during the procedure.

Quick Recovery

Most women return to normal activity within one to two days.

Benefit of Confirmation

An Essure Confirmation Test will verify that the inserts are placed correctly so that you can rely on Essure for birth control.

ESSURE PROCEDURE OVERVIEW

Step 1

Placing the Essure inserts

The Essure procedure is usually performed in your doctor's office. During the procedure, the doctor will place a tiny insert into each of your fallopian tubes. The inserts are soft and flexible, and are delivered with a tube through your vagina and cervix, and into your fallopian tubes. No incisions are needed.

Step 2

Waiting for the natural barrier to form

Over the next 3 months, your body will form tissue around the Essure inserts. The tissue forms a natural barrier within the fallopian tubes. The barrier prevents sperm from reaching the eggs that are produced every month. During the 3-month period, you **must** continue using another form of birth control to prevent pregnancy.

Since Essure does not contain hormones that interfere with your body's menstrual cycle, your ovaries will continue to release eggs. Since the eggs cannot be fertilized, they are simply absorbed back into your body.

Step 3

Essure Confirmation Test

After 3 months, a doctor will administer the Essure Confirmation Test. The test will verify that the inserts are in their correct location and your fallopian tubes are blocked. The doctor will use contrast dye and a special type of x-ray during the confirmation test.

IMPORTANT: FOR SOME WOMEN, IT MAY TAKE LONGER THAN 3 MONTHS FOR ESSURE TO COMPLETELY BLOCK THE FALLOPIAN TUBES, REQUIRING A REPEAT CONFIRMATION TEST AT 6 MONTHS. YOU MUST SEE YOUR DOCTOR FOR THE ESSURE CONFIRMATION TEST BEFORE YOU CAN RELY ON ESSURE FOR BIRTH CONTROL. UNTIL YOU RECEIVE CONFIRMATION FROM YOUR DOCTOR, YOU MUST CONTINUE TO USE ANOTHER FORM OF BIRTH CONTROL TO PREVENT PREGNANCY.

WARNINGS, PRECAUTIONS AND OTHER POTENTIAL RISKS

WARNING: Be sure you are done having children before you undergo the Essure procedure. Essure is a permanent method of birth control.

- The younger a woman is when she chooses to end her fertility, the more likely she is to regret her choice later.

WARNING: You must continue to use another form of birth control until you have your Essure Confirmation Test and your doctor tells you that you can rely on Essure for birth control.

- You can rely on Essure for birth control only after your doctor has reviewed your Essure Confirmation Test results and told you that you may rely. If you rely on Essure for birth control before having your Essure Confirmation Test, you are at risk of getting pregnant.
- Talk to your doctor about which method of birth control you should use for the 3 months after the procedure. Some women can remain on their current birth control. Other women, such as those using an intrauterine device or contraceptive (IUD or IUC), will need to switch to another method.
- It can take longer than three months for the Essure procedure to be effective. In rare cases, it has taken up to 6 months. Make sure to continue using an alternate form of birth control up until your doctor has reviewed your Essure Confirmation Test results and confirmed that you can rely on Essure for birth control.

Risks: During the Essure procedure

- You may experience mild to moderate pain.
- Your doctor may be unable to place one or both Essure inserts correctly.
- In rare cases, part of an Essure insert may break off. Your doctor may remove the piece or let it leave your body during your period.
- In rare cases, part of an Essure insert may puncture the fallopian tube. Surgery may be necessary to repair the puncture.
- Your body may absorb a large amount of the salt water solution used during the procedure.
- Your doctor may recommend a local anesthesia, which numbs the cervix. Ask your doctor about the risks associated with this type of anesthesia.

Risks: Immediately following the procedure

- You may experience mild to moderate pain and/or cramping, vaginal bleeding, and pelvic or back discomfort for a few days after the procedure. Some women experience nausea and/or vomiting or fainting. You should arrange to have someone available to take you home after the procedure.
- In rare instances, an Essure insert may be expelled from the body. This is usually detected during the Essure Confirmation Test.

Risks: During the Essure Confirmation Test

- Because the Essure Confirmation Test requires an x-ray, you will be exposed to very low levels of radiation. This is standard with most x-rays.
- In rare instances, women may experience spotting and/or infection.

Risks: Long-term

- **There are reports of chronic pelvic pain in women, possibly related to Essure.**
- **There are reports of the Essure insert migrating into the lower abdomen and pelvis. If this happens, it may be necessary to surgically remove the migrated device.**
- No birth control method is 100% effective. There is a chance that you can become pregnant after completing the Essure procedure. In the original premarketing studies for Essure, no pregnancies were reported for women who had the Essure inserts for up to 5 years. Although successful pregnancies have been reported with Essure devices in place, if you do become pregnant, the risks to you, the fetus, the pregnancy and childbirth are unknown.
- Women who have the Essure procedure are more likely to have an ectopic pregnancy if they get pregnant. Ectopic pregnancy is when the pregnancy occurs outside of the uterus. The pregnancy usually happens in one of the fallopian tubes. Ectopic pregnancies can be very serious or life-threatening.
- If you have the NovaSure® procedure, a procedure that removes the lining of the uterus to lighten or stop menstrual bleeding, after the Essure procedure, it is unknown if this will affect the blockage in your tubes, and your risk of pregnancy may increase.
- The Essure insert is made of materials that include a nickel-titanium alloy. Once placed inside the body, small amounts of nickel are released from the inserts. Patients who are allergic to nickel may have an allergic reaction to the inserts. Symptoms include rash, itching and hives.

Unknown Risks:

- The safety and effectiveness of Essure has not been established in women under 21 or over 45 years old.
- The safety and effectiveness of reversing the Essure procedure are not known.
- The safety and effectiveness of in vitro fertilization (IVF) after the Essure procedure are not known.
- The risks to you and your fetus if you get pregnant after the Essure procedure are not known.

WHAT TO EXPECT WITH ESSURE**Preparing for your procedure**

Your doctor will schedule your Essure procedure for a time soon after your menstrual period ends. This will make it easier for your doctor to see the openings of your fallopian tubes and place the inserts.

The day of your procedure

You will take a pregnancy test before or on the day of your procedure. This will ensure you are not pregnant. Your doctor may also give you medication to take 1 to 2 hours before your procedure to help you relax and open (dilate) your fallopian tubes. This will also help to reduce cramping. Talk to your doctor about the types of medications that are right for you.

During your procedure

Your doctor will first insert an instrument called a speculum inside your vagina. The speculum helps the doctor widen the opening and see inside. Then your doctor will insert a narrow, telescope-like instrument (hysteroscope) through your cervix and into your uterus. A camera attached to the hysteroscope sends video images to a video monitor that lets your doctor see inside your uterus. A salt water solution is used to expand the uterus. This makes it easier for your doctor to find the openings of your fallopian tubes.

The Essure insert is attached to the end of a small, flexible tube that passes through the hysteroscope and into your fallopian tube. Once the insert is placed, the tube is removed.

The procedure is then repeated to place an insert into your other fallopian tube. The entire process usually takes less than ten minutes.

Your doctor will schedule you for an Essure Confirmation Test for 3 months after the procedure.

IMPORTANT: Not all women will achieve successful placement of both Essure inserts. Fewer than 1 out of 12 women are not able to have one or both of the inserts placed. If this occurs, talk to your doctor about a second Essure procedure.

After your procedure

Most women are able to leave the doctor's office about 45 minutes after the procedure is completed. Most return to normal activities within one to two days. Call your doctor if you experience pain, bleeding, fever, vaginal discharge, or other symptoms following the procedure.

It takes about 3 months (sometimes longer) for your body to produce tissue around the inserts and form a barrier to prevent pregnancy. During that time you can still get pregnant. You must rely on another type of birth control to prevent pregnancy during this time period.

IMPORTANT: YOU MUST SEE YOUR DOCTOR FOR THE ESSURE CONFIRMATION TEST BEFORE YOU CAN RELY ON ESSURE FOR BIRTH CONTROL. UNTIL YOU RECEIVE CONFIRMATION FROM YOUR DOCTOR, YOU MUST CONTINUE TO USE ANOTHER FORM OF BIRTH CONTROL TO PREVENT PREGNANCY.

The Essure Confirmation Test

The Essure Confirmation Test verifies that the inserts are in the correct location and that your fallopian tubes are completely blocked.

The Essure Confirmation Test is a type of x-ray exam involving little to no pressure or discomfort. Although the test is similar to the HSG (hysterosalpingogram) test used to diagnose and treat infertility, it is considerably more comfortable.

During the test, a radiologist injects a special contrast dye into your uterus. The dye is visible on x-rays. This lets the radiologist look at your fallopian tubes to confirm that the inserts are properly placed and that your tubes are blocked.

ESSURE PATIENT ID CARD AND PROCEDURES AFTER ESSURE

After your Essure procedure, you will be given an Essure ID card. The ID card tells doctors and others that you have Essure inserts. Show the card when undergoing any procedure involving your uterus or fallopian tube. These include an MRI, D&C, hysteroscopy, endometrial biopsy, or endometrial ablation. Body areas near the inserts may be obscured when they are seen on x-rays, MRIs and other imaging.

IMPORTANT: Women with Essure inserts who undergo MRI procedures should tell their doctor they have been implanted with Essure inserts.

FREQUENTLY ASKED QUESTIONS**Can I trust Essure to prevent pregnancy?**

Yes, the Essure procedure is 99.83% effective, based on the original premarketing studies for Essure.

Pregnancies have been reported among the hundreds of thousands of women who have completed the procedure since it became commercially available. Many of those pregnancies were a result of not having completed the procedure (for example not undergoing the Essure Confirmation Test) or not following instructions. In some cases, the Essure Confirmation Test results were misinterpreted by the person reading the test.

Is Essure painful?

There may be some pain associated with placing Essure. Some women report mild discomfort, pain, and cramping during or after the placement procedure. Symptoms may be similar to what they might experience in their normal monthly cycle. There are reports of chronic pelvic pain in women, possibly related to Essure.

Is Essure reversible?

No, the Essure procedure is not reversible. Like having your tubes tied or a vasectomy for men, Essure is permanent birth control. You need to be sure you are done having children before you decide to have the Essure procedure.

Will I still get my period after the Essure procedure?

Yes, you will still have a period. Some women find that their period may become slightly lighter or heavier after the procedure. These changes are often temporary. They may also be due to you stopping your previous hormonal birth control, rather than the Essure procedure.

What are the Essure inserts made of?

The inserts are made from polyester fibers, nickel-titanium and stainless steel. These same materials have been used for many years in cardiac stents and other medical devices placed in other parts of the body.

Is Essure covered by my insurance?

Some or all of the costs of the Essure procedure are covered by most insurance providers. When the Essure procedure is done in a doctor's office, your cost may be as low as the usual co-pay amount for an office visit/procedure. This depends on your insurance plan.

Review your insurance coverage with your doctor and provider before having the procedure. When you speak with your insurance provider, ask for your plan's specific coverage and reimbursement guidelines for hysteroscopic sterilization (code 58565). This is the code that most plans use to categorize the Essure procedure. Specify if the procedure will be done in a doctor's office, clinic, hospital, etc. The location may affect the amount you will need to cover out of your own pocket.

QUESTIONS TO ASK YOUR DOCTOR

If you are considering having the Essure procedure, here are some questions you might ask your doctor:

- Is Essure right for me compared to other birth control methods?
- Where will my Essure procedure be performed?
- What type of medications will be used before and/or during my procedure?
- How should I prepare?
- What are my options if both inserts cannot be placed on the first attempt?
- How do I schedule my Essure Confirmation Test?
- Can I continue to use my current method of birth control until I have the results of my Essure Confirmation Test?

HOW ESSURE PERFORMED IN CLINICAL STUDIES

The effectiveness and safety of Essure was measured in clinical studies in the United States, Australia, and Europe. Over 700 women between the ages of 21 and 45 were studied.

Was Essure effective in preventing pregnancy?

The primary goal of any birth control method is to prevent pregnancy. Every birth control method has a measured effectiveness rate.

In the original Essure clinical studies, zero (0) pregnancies were reported in women who had the Essure inserts for up to 5 years. The effectiveness rate with Essure was estimated to be 99.83%.

Was the Essure placement procedure successful?

In only a few instances (fewer than 1 out of 12 women), the doctor was unable to place one or both Essure inserts in the fallopian tubes.

The majority of women (96.5%) were able to achieve complete blockage of their fallopian tubes within the expected 3-month period. A small percentage of women (3.5%) took up to 6 months.

Was the Essure procedure safe?

In the original premarketing study, some women reported mild to moderate adverse events during and after the procedure. During the procedure, the most common problem reported was mild to moderate pain (9.3% of women). Some of the women in the study reported moderate pain (12.9% of women) and/or cramping (29.6% of women) on the day of the procedure. A smaller percentage of women reported nausea/vomiting (10.8%) and vaginal bleeding (6.8%). Eighty-eight percent (88%) of women rated tolerance of the placement procedure as good, very good, or excellent.

HOW ESSURE HAS PERFORMED OUTSIDE OF CLINICAL TRIALS BETWEEN 2002 AND 2010

Pregnancies have been reported with “real world” use of Essure. The table below shows the number of **reported** pregnancies and the most likely reason the pregnancy occurred. The actual number of pregnancies in women with Essure may be higher.

Pregnancies Reported With Essure, 2002-2010*

Potential Contributing Factor	Within The U.S. Number (%)***	Outside the U.S. Number (%)****	Total Number (%)
Patient did not comply with instructions	213 (32%)	16 (18%)	229 (31%)
Essure insert perforated fallopian tube**	91 (14%)	4 (5%)	95 (13%)
Essure inserts were not placed correctly**	32 (5%)	13 (15%)	45 (6%)
Doctor did not comply with instructions	22 (3%)	13 (15%)	35 (5%)
Pregnant at time of placement	26 (4%)	6 (7%)	32 (4%)
Essure Confirmation Test misinterpreted or incomplete**	28 (4%)	0 (0%)	28 (4%)
Essure insert was expelled from the fallopian tube**	20 (3%)	4 (5%)	24 (3%)
Fallopian tube was open/unobstructed**	19 (3%)	1 (1%)	20 (3%)
Unable to determine cause	209 (32%)	31 (35%)	240 (32%)
Total	660 (100%)	88 (100%)	748 (100%)

* Number of pregnancies reported worldwide since commercial launch in 2001 through the end of 2010. Over that time period, 497,306 Essure kits were sold. This works out to be a “reported” pregnancy rate of 0.15%. The “true” overall pregnancy rate is unknown, because the actual number of pregnancies and the number of devices actually implanted is not known.

** Most of these pregnancies are due to Essure Confirmation Tests that were misinterpreted. In most cases, the Essure inserts were not properly placed but the x-ray indicated that fallopian tubes were blocked.

*** Only HSG was used for these procedures.

**** TransVaginal Ultrasound, HSG, and/or x-ray were acceptable confirmation tests.

Many pregnancies reported by women who had the Essure procedure were preventable. For example, the patient or physician did not follow the Instructions for Use. So it is important that you follow the instructions provided to you by your doctor.

Make sure to:

- Continue to use an alternative form of birth control for at least 3 months following the placement of the Essure inserts.
- See your doctor to have your Essure Confirmation Test 3 months after your Essure inserts are placed.
- Have your doctor confirm that you can rely on Essure inserts for your birth control.

COMPARISON OF PERMANENT BIRTH CONTROL METHODS

The following table provides information about a variety of permanent birth control methods. It includes data on the percentage of women likely to become pregnant within a year and five years while utilizing that method. For a complete list, visit the FDA website at www.fda.gov and search for the Birth Control Guide.

ESSURE

Soft, flexible inserts are delivered through the vagina and uterus and placed in each fallopian tube. A natural barrier forms around the inserts and prevents sperm from reaching the eggs.

No incision is necessary to deliver or place the inserts. General anesthesia is not necessary during the procedure.

Failure rate*	Recovery Time	Pain/discomfort	Proof that method was effective
<ul style="list-style-type: none"> • 0.3 out of 1,000 women at 1 year • 1.7 out of 1,000 women at 5 years 	1-2 days or sooner	<ul style="list-style-type: none"> • Cramping • Discharge 	Yes. Three months after the procedure, the Essure Confirmation Test confirms correct insert location and blockage of the fallopian tubes.

TUBAL LIGATION⁵

The fallopian tubes are blocked so that sperm is unable to reach the eggs. One of three methods is used to block the tubes:

- Clamping with metal clips or plastic rings that remain in the body.
- Cutting away a section of the tube.
- Burning a portion of the tube.

The procedure requires an incision and is performed under general anesthesia. Gas is used to expand the abdomen. Stitches or staples are then used to close the incision.

Failure rate*	Recovery time	Pain/discomfort	Proof that method was effective
<ul style="list-style-type: none"> • 5.5 out of 1,000 women at 1 year • 13.1 out of 1,000 women at 5 years • 18.5 out of 1,000 women at 10 years 	4-6 days	<ul style="list-style-type: none"> • Cramping • Discharge • Pain at the incision area • Bruising near the incision area • Bloating abdomen and/or sharp pains in the neck or shoulder (due to gas used) • Tired and achy feeling 	None.

*Expected number of pregnancies with this method over time.

VASECTOMY (MEN)

The two vas deferens tubes that propel sperm through the urethra are tied in two places with permanent sutures. Between the ties, the tubes are severed using one of three methods:

- Burning a portion of the tube.
- Cutting the tube.
- Blocking the tube with clips or clamps that remain in the body.

The scrotal area is shaved and cleaned with an antiseptic solution. An incision or puncture is made into the scrotum (the sac containing the testicles). Stitches or staples are used to close the cuts.

Failure rate*	Recovery time	Pain/discomfort	Proof that method was effective
<ul style="list-style-type: none"> • 7.4 out of 1,000 women at 1 year • 11.3 out of 1,000 women at 5 years 	2-3 days	<ul style="list-style-type: none"> • Bruising • Pain and swelling in the testicles 	Yes. A follow-up sperm count test is performed 3 months after the vasectomy to confirm no sperm are evident.

COMPARISON OF TEMPORARY BIRTH CONTROL METHODS

The following table provides information on a variety of temporary birth control methods. It includes data on the percentage of women likely to become pregnant within a year while utilizing that method. For a complete list, visit the FDA website at www.fda.gov and search for the Birth Control Guide.

Not all temporary methods of birth control listed below can be used during the 3-month waiting period before the Essure Confirmation Test. Talk to your doctor about which form of temporary birth control is right for you.

ORAL CONTRACEPTIVES (COMBINATION ESTROGEN/PROGESTIN PILL)

An estrogen/progestin-based pill that suppresses ovulation.

Failure rate*	Risks	Routine
80 out of 1,000 women	Dizziness, nausea, changes in menstruation, mood, and weight gain. Rare events include cardiovascular disease, including high blood pressure, blood clots, heart attack, and stroke.	Must be taken daily.

ORAL CONTRACEPTIVES (PROGESTIN-ONLY PILL)

A progestin-based pill that inhibits fertilization.

Failure rate*	Risks	Routine
80 out of 1,000 women	Irregular bleeding, weight gain, breast tenderness, and less protection against ectopic pregnancy.	Must be taken daily.

INJECTION (DEPO PROVERA®)

A progestin-containing injection that inhibits ovulation and fertilization.

Failure rate*	Risks	Routine
30 out of 1,000 women	Irregular bleeding, weight gain, breast tenderness, and headaches.	One injection every 1-3 months.

VAGINAL CONTRACEPTIVE RING (NUVARING®)

A flexible ring inserted in the vagina that releases progestin and estrogen to prevent ovulation and fertilization.

Failure rate*	Risks	Routine
80 out of 1,000 women	Vaginal discharge, vaginitis, irritation, and other risks similar to those posed by oral contraceptives.	Inserted by the woman and kept in place for 3-week intervals. If expelled for more than 3 hours during the 3-week interval, another method of birth control must be used.

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A patch worn on the body that releases progestin and estrogen to prevent ovulation and fertilization.

Failure rate* 80 out of 1,000 women	Risks Similar to the oral estrogen-progestin pill	Routine A new patch must be applied every week other than the week of the menstrual period.
---	---	---

IUD/IUC (MIRENA®)

A device placed in the uterus (by a doctor) that emits hormones, preventing ovulation.

Failure rate* 2 out of 1,000 women	Risks Ovarian cysts, pelvic inflammatory disease, perforation of the uterus, embedding into the uterus, cramps, bleeding, miscarriage, premature birth, breast cancer, nausea, mood swings, headaches, nervousness, inflammation/pain of vagina/uterus, back pain, weight gain, acne, hypertension, and changes in menstrual cycle.	Routine Remains in place for 1 to 5 years.
--	---	--

IUD/IUC (PARAGUARD®)

A device placed in the uterus (by a doctor) that releases copper, preventing ovulation and fertilization.

Failure rate* 8 out of 1,000 women	Risks Pelvic inflammatory disease, perforation of the uterus, embedding into the uterus, cramps, bleeding, vaginal discharge, allergic reaction, expulsion, anemia, ectopic pregnancy, life-threatening infection, miscarriage, premature birth, Wilson’s disease, vaginal infection, inflammation/pain of vagina/uterus, back pain, pain during sex, fainting, and changes in menstrual cycle.	Routine Remains in place for 1 to 10 years.
--	---	---

MALE CONDOM

A sheath placed over the penis that prevents passage of sperm.

Failure rate* 150 out of 1,000 women	Risks Irritation, allergic reactions, and reduced effectiveness if used with oil-based lubricants.	Routine Applied immediately before intercourse and used only once.
--	--	--

FEMALE CONDOM

A lubricated sheath placed in the vagina to prevent sperm from entering the uterus.

Failure rate* 210 out of 1,000 women	Risks Irritation and allergic reactions.	Routine Applied immediately before intercourse and used only once.
--	--	--

DIAPHRAGM WITH SPERMICIDE

A dome-shaped rubber disk with a flexible rim that covers the cervix. The disk prevents sperm from reaching the uterus. A spermicide must be applied to the dome of the diaphragm before insertion.

Failure rate* 160 out of 1,000 women	Risks Irritation, allergic reactions, urinary tract infection, and risk of toxic shock syndrome.	Routine Inserted before intercourse and left in place for 6 to 24 hours afterward. For repeated intercourse, spermicide must be added without removing the diaphragm.
--	--	---

SPERMICIDE

A foam, cream, jelly, film, suppository or tablet containing a sperm-killing chemical (nonoxynol-9).

Failure rate* 290 out of 1,000 women	Risks Irritation, allergic reactions, and urinary tract infections.	Routine Instructions vary. Inserted 5 to 90 minutes before intercourse and usually left in place for at least 6 to 8 hours afterward.
--	---	---

PERIODIC ABSTINENCE/RHYTHM METHOD

Deliberately refraining from having sexual intercourse during times when pregnancy is more likely.

Failure rate* 250 out of 1,000 women	Risks None	Routine Requires continuous monitoring of ovulation cycle and body temperature.
--	----------------------	---

Data adapted from Trussell J. Contraceptive efficacy. In: Hatcher RA, Trussell J, Nelson AL, Cates W, Stewart FH, Kowal D. *Contraceptive Technology: Nineteenth Revised Edition*. New York, NY: Ardent Media; 2007.

essure®

REFERENCES

1. Cooper JM, Carignan CS, Cher D, Kerin JF. Microinsert nonincisional hysteroscopic sterilization. *Obstet Gynecol.* 2003;102:59-67. **2.** Kerin JF, Cooper JM, Price T, et al. Hysteroscopic sterilization using a micro-insert device: results of a multicentre phase II study. *Hum Rep.* 2003;18:1223-1230. **3.** Syed R, Levy J, Childers ME. Pain associated with hysteroscopic sterilization. *JLS.* 2007;11:63-65. **4.** PMA: P020014/S9. **5.** Peterson HB, Xia Z, Hughes JM, et al. The risk of pregnancy after tubal sterilization: findings from the U.S. Collaborative Review of Sterilization. *Am J Obstet Gynecol.* 1996;174:1161-1170. **6.** Essure ESS305: Instructions For Use 2011: 1-10.

Essure, Bayer and the Bayer Cross are registered trademarks of Bayer.
© 2002 Bayer, All rights reserved.

August 2015